



**Passion in Health.  
Since 1988.**

Product range

## Meditrade UK

Meditrade UK is a member of the Meditrade Group located in central London. Meditrade Group has always had a large number of UK customers, however Brexit meant that the quality of our overseas relationships were going to decline, naturally. Because of this Meditrade made the decision to open an office on the ground in the UK, which would allow Meditrade UK to offer favourable business solutions to the existing UK customer base. In addition, we can offer Irish customers services from Europe (Germany – HQ).

At Meditrade UK we have developed a logistic concept that offers a flexible and consistent supply of our products into the UK and Irish market. Part of our logistics concept is the multi-stage logistics hub, which is being built in Gainsborough, Lincolnshire. The strategically central geographical placement establishes an efficient distribution network which allows us to promptly assist our customer's demands at all times.





## Meditrade – 30 years of innovation

In 1988 Meditrade was founded under the maxim 'Passion in Health' - Meditrade Group combines more than 30 years of expertise in manufacturing, distributing and supplying medical consumables to healthcare facilities, manufacturers, distributors, retailers and wholesalers. Starting as a manufacturer of medical disposable gloves, Meditrade has now become one of the leading manufacturers of medical equipment in Europe with a wide range of products from wound care, protective clothing, operating theatre sets, incontinence products and disinfection products for the skin, hands, and surfaces.

As manufacturers and system suppliers of medical consumables, we see ourselves as partners and service providers for the basic care of all medical facilities. Quality, innovation and beneficial care options are our daily motivation to make a small but essential contribution to modern and future-oriented healthcare. The employees of Meditrade stand by this core essence and fulfil these demanding tasks for a higher standard of health with competence and professionalism. With factories all over the world and logistics centres in Germany and the United Kingdom, Meditrade meets the increasing demands and requirements of its customers and always assists its customer's requests wherever they are.

## Our product range

- Gloves
- Protective clothing
- Wound care
- Sets and operating theatre supplies
- Incontinence
- Disinfection and care
- Hygiene, cleaning and cleaning products
- Ward and patient accessories



## Logistics

Our approach to logistics ensures the flexible and reliable supply of our products to the entire European market.

- Main site with storage space for 20,000 pallets
- Our expanded logistics centre in Neuenkirchen offers space for 28,000 pallets
- 12.5-metre high racking with a semi-automatic picking system
- Completion of our logistics centre in Gainsborough in 2022




## About us

- Meditrade products are manufactured in our own facilities all over the world in order to ensure the high quality of our products
- We offer a comprehensive product portfolio to supply our customers with a full range of medical consumables
- Leading provider of PPE and medical consumables in Europe, with over thirty years of experience in the industry
- A global company that stays true to the values of a family-run business
- >150 employees



## Locations



## Our core competencies

- Products**  
High quality and innovative  
Quality products for  
healthcare provision
- Knowledge**  
Industry expertise backed by  
thirty years of experience  
in the healthcare sector
- Organisation**  
A major industry supplier  
with several  
Logistics centres in Europe
- Service**  
A responsive customer service,  
which lets us adapt to  
match the needs of our customers

## Gloves

<b>Latex examination gloves</b>	
Reference™	9
Gentle Skin® sensitive	9
Gentle Skin® grip	10
Gentle Skin® compact +	10
Gentle Skin® Aloecare™	11
Gentle Skin® black	12
<b>Nitrile examination gloves</b>	
Nitril® 3000	13
Nitril® NextGen®	13
Nitril® Sensory®	14
Nitril® 3000 x-long blue	15
Nitril® black	16
Nitril® Viola	17
Nitril® Magenta	18
Nitril® LIOX	20 / 21
<b>Vinyl examination gloves</b>	
Vinyl 2000 PF	22
Vinyl Blue Gen	23
<b>Cotton and hygienic gloves</b>	
Cotton glove	24
Froggy™	25
<b>Finger cots</b>	
Latex finger cots	28
<b>Household gloves</b>	
MyClean® easy touch	27
<b>Sterile examination gloves</b>	
Copolymed®	28
Gentle Skin® steril	29
Nitril® steril	30
<b>Surgical gloves</b>	
Reference™ OP	33
Gentle Skin® Premium OP	34
Gentle Skin® Superior OP	35
Gentle Skin® Smooth OP	36
Gentle Skin® Micro OP®	37
Maxitex® UG	38
Gentle Skin® Securitex® OP	39
Gentle Skin® Isopretex® OP	40
Neopretex® OP	41
Finessis® Aegis®	43
Finessis® Zero®	44
Finessis® Corium®	45

## Protective clothing

<b>Masks</b>	
Suavel® Protec Plus	48
Suavel® Protec	49
Suavel® Protec Black	49
Suavel® Comfort Plus	50
Suavel® Comfort	50
Suavel® Antifog	51
Suavel® Antifluid	51
Suavel® Sensima	52
Respima FFP2	52
<b>Caps</b>	
Suavel® Apollo	53
Suavel® Astrid	53
<b>Scrubs / gowns / aprons</b>	
BeeSana® PP gown	54
BeeSana® PP/PE gown	55
BeeSana® SMS gown	55 / 56
Suavel® Yeti Coverall	56
PE disposable apron	57
Suavel® OP PRO Standard / Expert	58 / 59

## Wound care

<b>Gauze swabs</b>	
BeeSana® Gauze cutting	62
BeeSana® Medical gauze	62
<b>Gauze Swabs &amp; Dressing pads</b>	
BeeSana® Gauze swab	63 - 65
BeeSana® Slit gauze swab	66
BeeSana® Non-woven swab	67
BeeSana® Slit non-woven swab	68
ABE® drainal	69
ABE® tracheal	69
BeeSana® Absorbent dressing pad	70
BeeSana® Eye pad	71
ABE® nose bandage	72
<b>Swabs and tamponades</b>	
BeeSana® Gauze balls	73
BeeSana® Pointed swabs	74
ABE® ribbon gauze	74
BeeSana® Cellulose swabs	75
Alcohol swabs	75
<b>Treatment care and accessories</b>	
Cotton buds	76
Wooden tongue depressor	77
<b>Wound plasters and bandages</b>	
ABE® injection plaster	78
ABE® light	78
ABE® Absorbent adhesive dressing	79
ABE® derm	80
ABE® derm plus	80
<b>Fixing plasters/non-woven materials</b>	
ABE® cannula retention dressing	81
ABE® non-woven material	81
ABE® silk	82
<b>Elastic fixing bandages, self-adhering, crepe</b>	
BeeSana® SafeCut	83
ABE® gauze bandage	84
ABE® Last	85
ABE® Last crepe	86
ABE® fix color	87
ABE® Last cohesive	88
<b>Compression &amp; Support bandages</b>	
ABE® Ideal 100	89
ABE® universal	90
ABE® lastic	91
ABE® padded bandage	92
ABE® extra	92
ABE® extra cohesive	93
ABE® Ven	93
ABE® tubular bandage	94
ABE® tubular net bandage	95
ABE® Gilchrist's bandage	96
ABE® body bandage	97
ABE® heart catheter pressure bandage	98
ABE® limb holder	98
ABE® pressure bandage roll	99
Scrotum roll	99
<b>Wound care X-ray</b>	
BeeSana® X-Ray lap sponge	100 / 101
BeeSana® X-Ray gauze swab	102 / 103
BeeSana® X-Ray gauze balls	104
BeeSana® X-Ray preparation swabs	105
BeeSana® X-Ray pointed swabs	106
BeeSana® X-Ray gauze balls	107
ABE® X-Ray ribbon gauze	107

## Set solutions

BeeSana® Ward catheter set	112
BeeSana® PEG bandage set	113
BeeSana® Bandage change set	114
BeeSana® Suture removal set	114
BeeSana® Ward protection set	115

## Incontinence products

<b>Incontinence briefs</b>	
Superslip®	116
Superslip® Air Comfort	117
BeeSana® Medi-Pants	117
<b>Fixation briefs</b>	
Light fixation briefs	118
Panty fixation briefs	119
<b>Incontinence liners/pads</b>	
Liberty Light®	120
<b>Incontinence pads</b>	
Superform®	120
Superform® Air	121
BeeSana® Absorbent cushion	121
BeeSana® Absorbent pad	122
BeeSana® Rectangular pad	122
BeeSana® Male incontinence pad	123
BeeSana® Female incontinence pad	123
<b>Incontinence undersheets</b>	
Layered undersheet	124
BeeSana® Comflock Super	125
<b>Nursing pads</b>	126

## Disinfection / care

<b>Disinfection</b>	
Ethasept®	130
Ethasept® Gel	131
Medizid® Rapid+	132
Medizid® Rapid+ / Rapid+ XXL disinfectant wipes	133
Medizid® Rapid QF Flowpack	134
Medizid® AF+	135
Medizid® Alpha+	136
Medizid® Alpha+ XXL disinfectant wipes	137
Medizid® Alpha+ Flowpack	138
Instrusol® AF+	140 / 141
<b>Care</b>	
Gentle Med® wash lotion	143
Gentle Med® skin care lotion (O/W)	144
Gentle Med® cream lotion (W/O)	144
Gentle Med® active gel	145
Gentle Med® oil bath	145
Gentle Med® skin cleansing and care foam	146
Gentle Med® moist care wipes	146
<b>Application aids</b>	
Metering pumps/attachments	147
(Foam) spray heads	147
Canister tap/canister key	147
Dosing bottle/graduated flask	147
Disinfectant dispenser	148 / 149

## Ward and patient supplies

<b>Patient and personal hygiene</b>	
BagBath® Haircare	152
BagBath®	152
BeeSana® Wash glove	154
Protective towel	154
Protective bib	155
BeeSana® Disposable briefs	155
Anti-slip sock universal	156
Surgical sock	156
PE shoe covers	157
Cellulose sponge sheets	157
Layered cellulose	158
Rollicel® / medical crepe	158
Protective sheets	158
Protective covering	159
Multi-purpose cloths	159
Protective bed coverings	159
Universal cloth	160
Wipe cloths	160
Cosmetic wipes	161
V-fold paper hand towels	161
Sick bag (BaSick Bag)	162
Kidney dish	162
Disposable comb	162
Disposable razor	163
Urine bottle	163
Urine bottle holder	163
<b>Diagnosis and accessories</b>	
Nurse's stethoscope	164
Tourniquet	164
Dosicair®	165
Urine beaker	165
Ultrasound scan gel	165
<b>Ward supplies and accessories</b>	
Hot-Cold Pack	166
Hot-Cold Pack Sleeve	166
Rescue blanket	167
Medicine Pod	167
Medication dispenser	168
Sharps Bin	168
Glove dispenser	169
Rubbish bags	170
Kitchen towel	171
Toilet paper	171





**We know gloves like the back of our hand**

When making our examination and surgical gloves, we at Meditrade understand that it is vital to provide gloves that combine mechanical strength with good touch sensitivity. Above all, our first priority is to protect the user and the patient, in accordance with regulatory requirements. We also invest a lot in product development know-how and in our production sites. Our huge range of latex glove models are complimented with extremely cost-effective gloves in vinyl and other alternative materials, including gloves from pure natural rubber and nitrile gloves made from nitrile butadiene rubber, ensuring we can always offer the right product for the user.

# Reference™

Our powdered latex examination glove

**Features**

- Latex examination glove, natural (DIN EN 455 standard)
- lightly powdered
- micro-roughened
- elastic and tear-proof
- ambidextrous



**Non-sterile**

Size	REF	Pack contents	Shipping unit	PZN
XS	1121XS	100 pcs	10 packs	06177515
S	1121S	100 pcs	10 packs	06177521
M	1121M	100 pcs	10 packs	06177538
L	1121L	100 pcs	10 packs	06177544
XL	1121XL	100 pcs	10 packs	02243899

**Usage**

In all areas where the user is exposed to a risk of contamination. Suitable for highly infectious situations in hospitals and in general practice. Suited to use in the manufacturing sector.

# Gentle Skin® sensitive

Our most popular latex examination glove

**Features**

- Latex examination glove, natural (DIN EN 455 standard)
- PPE Cat III
- powder-free, low protein
- micro-roughened
- elastic and tear-proof
- ambidextrous



**Non-sterile**

Size	REF	Pack contents	Shipping unit	PZN
XS	1221RXST	100 pcs	10 packs	02243563
S	1221RST	100 pcs	10 packs	02243586
M	1221RMT	100 pcs	10 packs	02243592
L	1221RLT	100 pcs	10 packs	02243600
XL	1221RXLT	100 pcs	10 packs	02243617

**Usage**

In all areas where the user is exposed to a risk of contamination. Suitable for highly infectious situations in hospitals and in general practice. Suitable for use in nursing, by lab personnel, in the pharmaceutical industry, in cleaning work and in manufacturing.

# Gentle Skin<sup>®</sup> grip

Our easy-grip latex glove

### Features

- Latex examination glove, natural (DIN EN 455 standard)
- PPE Cat III
- good, secure grip
- powder-free, low-protein
- micro-roughened
- elastic and tear-proof
- ambidextrous



### Non-sterile

Size	REF	Pack contents	Shipping unit	PZN
XS	1221GRIP-XS	100 pcs	10 packs	02470253
S	1221GRIP-S	100 pcs	10 packs	02470282
M	1221GRIP-M	100 pcs	10 packs	02470299
L	1221GRIP-L	100 pcs	10 packs	02470307
XL	1221GRIP-XL	100 pcs	10 packs	02470313

### Usage

Suitable for use in dental surgery: its easy-grip qualities enable the safe, non-slip handling of instruments and machinery.

# Gentle Skin<sup>®</sup> compact +

Our touch sensitive latex examination glove

### Features

- Latex examination glove, natural (DIN EN 455 standard)
- our popular disposable glove for manufacturing
- good, secure grip
- powder-free
- micro-roughened
- elastic and tear-proof
- ambidextrous



### Non-sterile

Size	REF	Pack contents	Shipping unit	PZN
XS	1221XSI	100 pcs	10 packs	05522453
S	1221SI	100 pcs	10 packs	05522447
M	1221MI	100 pcs	10 packs	05522430
L	1221LI	100 pcs	10 packs	05522424
XL	1221XLI	100 pcs	10 packs	05521198

### Usage

Usable in manufacturing and ideal for household and cleaning work.

# Gentle Skin<sup>®</sup> Aloecare<sup>™</sup>

Our latex glove with an Aloe vera coating

### Features

- Latex examination glove, green (DIN EN 455 standard)
- inner coating with Aloe vera
- powder-free, low protein
- micro-roughened
- elastic and tear-proof
- ambidextrous



### Usage

Latex glove for applications such as cosmetics. The special inner coating with Aloe vera offers benefits for users with stressed skin.

### Non-sterile

Size	REF	Pack contents	Shipping unit	PZN
XS	1231XS	100 pcs	10 packs	05030230
S	1231S	100 pcs	10 packs	02244516
M	1231M	100 pcs	10 packs	02244580
L	1231L	100 pcs	10 packs	02244597
XL	1231XL	100 pcs	10 packs	05030247



# Gentle Skin® black

Our black latex glove

**Features**

- Latex examination glove, black (DIN EN 455 standard)
- powder-free, low protein
- micro-roughened
- elastic and tear-proof
- ambidextrous



**Usage**

The latex glove for hair stylists, beauticians and tattoo/piercing artists. Usable by caterers and chefs, and for general applications in cleaning.

**Non-sterile**

Size	REF	Pack contents	Shipping unit	PZN
XS	1224XS	100 pcs	10 packs	11352676
S	1224S	100 pcs	10 packs	11352682
M	1224M	100 pcs	10 packs	11352736
L	1224L	100 pcs	10 packs	11352771
XL	1224XL	100 pcs	10 packs	11352802

# Nitril® 3000

The white nitrile examination and protective glove

**Features**

- Nitrile examination glove, white (DIN EN 455 standard)
- PPE Cat III
- powder-, protein- and latex-free
- micro-roughened
- elastic and tear-proof
- ambidextrous, with rolled edge
- suitable for use as a food contact material



**Non-sterile**

Size	REF	Pack contents	Shipping unit	PZN
XS	1280XS	100 pcs	10 packs	06966703
S	1280S	100 pcs	10 packs	06967967
M	1280M	100 pcs	10 packs	06968056
L	1280L	100 pcs	10 packs	06968300
XL	1280XL	100 pcs	10 packs	06968381

**Usage**

Suitable for handling cleaning products and disinfectants, and for other applications in laboratories, the pharmaceutical industry, for professional cleaning, in the catering business and for handling foodstuffs in general. Suitable for users with a diagnosed (type I) latex allergy.

# Nitril® NextGen®

The stretchy nitrile examination and protective glove

**Features**

- Nitrile examination glove, light blue (DIN EN 455 standard)
- PPE Cat III
- powder-, protein- and latex-free
- micro-roughened
- elastic and tear-proof
- ambidextrous, with rolled edge
- suitable for use as a food contact material



**Non-sterile**

Size	REF	Pack contents	Shipping unit	PZN
XS	1283XS	100 pcs	10 packs	09102239
S	1283S	100 pcs	10 packs	09102245
M	1283M	100 pcs	10 packs	09102251
L	1283L	100 pcs	10 packs	09102268
XL	1283XL	100 pcs	10 packs	09102274

**Usage**

Suitable for handling cleaning products and disinfectants, and for other applications in laboratories, the pharmaceutical industry, for professional cleaning, in the catering business and for handling foodstuffs in general. Suitable for users with a diagnosed (type I) latex allergy.



# Nitril<sup>®</sup> Sensory<sup>®</sup>

Our touch-sensitive nitrile examination glove

### Features

- Nitrile examination glove, blue, white and violet (DIN EN 455 standard)
- PPE Cat III
- powder-, protein- and latex-free
- micro-roughened
- elastic and tear-proof
- ambidextrous, with rolled edge
- suitable for use as a food contact material

### Usage

Suitable for handling cleaning products and disinfectants, and for other applications in laboratories, the pharmaceutical industry, for professional cleaning, in the catering business and for handling foodstuffs in general. Suitable for users with a diagnosed (type I) latex allergy.



## Non-sterile

	Size	REF	Pack contents	Shipping unit	PZN
Blue	XS	2283XS	200 pcs	10 packs	14310075
	S	2283S	200 pcs	10 packs	14272305
	M	2283M	200 pcs	10 packs	14272311
	L	2283L	200 pcs	10 packs	14272328
	XL	2283XL	200 pcs	10 packs	14272334
White	XS	2280XS	200 pcs	10 packs	14306783
	S	2280S	200 pcs	10 packs	14272268
	M	2280M	200 pcs	10 packs	14272274
	L	2280L	200 pcs	10 packs	14272280
	XL	2280XL	200 pcs	10 packs	14272297
Violet	XS	2285XS	200 pcs	10 packs	-
	S	2285S	200 pcs	10 packs	-
	M	2285M	200 pcs	10 packs	-
	L	2285L	200 pcs	10 packs	-
	XL	2285XL	200 pcs	10 packs	-

# Nitril<sup>®</sup> 3000 x-long blue

Our long cuff examination and protective glove

### Features

- Nitrile examination glove, blue (DIN EN 455 standard)
- PPE Cat III
- with long cuff (approx. 29 cm) for better protection
- powder-, protein- and latex-free
- micro-roughened
- elastic and tear-proof
- ambidextrous
- suitable for use as a food contact material



### Usage

Thanks to its long cuff, this glove offers effective (splash) protection both to the hands and to parts of the forearms when handling liquids. This makes it suited to the handling of chemicals, cytostatic agents, detergents and disinfectants. The glove can also be used in the lab, pharmaceuticals and catering. Also for latex allergy sufferers with a diagnosed type I allergy.

## Non-sterile

Size	REF	Pack contents	Shipping unit	PZN
S	1282S	100 pcs	10 packs	00420825
M	1282M	100 pcs	10 packs	00420802
L	1282L	100 pcs	10 packs	00400857
XL	1282XL	100 pcs	10 packs	00400834

# Nitril® black

Our black nitrile examination and protective glove

### Features

- Nitrile examination glove, black (DIN EN 455 standard)
- PPE Cat III
- powder-, protein- and latex-free
- micro-roughened
- elastic and tear-proof
- ambidextrous, with rolled edge
- suitable for use as a food contact material



### Usage

Suitable for handling cleaning products and disinfectants, and for other applications in laboratories, the pharmaceutical industry, for professional cleaning, in the catering business and for handling foodstuffs in general. Suitable for users with a diagnosed (type I) latex allergy.

### Non-sterile

Size	REF	Pack contents	Shipping unit	PZN
S	1284S	100 pcs	10 packs	12614730
M	1284M	100 pcs	10 packs	12614747
L	1284L	100 pcs	10 packs	12614753
XL	1284XL	100 pcs	10 packs	12614782

# Nitril® Viola

The violet nitrile examination and protective glove

### Features

- Nitrile examination glove, violet (DIN EN 455 standard)
- PPE Cat III
- powder-, protein- and latex-free
- with silk-sericin coating
- micro-roughened
- elastic and tear-proof
- ambidextrous, with rolled edge
- suitable for use as a food contact material



### Usage

Suitable for handling cleaning products and disinfectants, and for other applications in laboratories, the pharmaceutical industry, for professional cleaning, in the catering business and for handling foodstuffs in general. Suitable for users with a diagnosed (type I) latex allergy.

### Non-sterile

Size	REF	Pack contents	Shipping unit	PZN
XS	1285XS	100 pcs	10 packs	-
S	1285S	100 pcs	10 packs	-
M	1285M	100 pcs	10 packs	-
L	1285L	100 pcs	10 packs	-
XL	1285S	100 pcs	10 packs	-



# Nitril<sup>®</sup> Magenta

Our pink examination and protection glove with silk sericin coating

## Features

- Nitrile examination glove, pink (DIN EN 455 standard)
- PPE Cat III
- powder-, protein- and latex-free
- with silk-sericin coating
- micro-roughened
- elastic and tear-proof
- ambidextrous, with rolled edge
- suitable for use as a food contact material



## Usage

A favourite in the cosmetics industry. Suitable for handling cleaning products and disinfectants, and for applications in laboratories, in the catering business and for handling foodstuffs in general. The silk-sericin inner coating moisturises dry and sensitive skin.

## Non-sterile

Size	REF	Pack contents	Shipping unit	PZN
S	1287S	100 pcs	10 packs	11584731
M	1287M	100 pcs	10 packs	11584694
L	1287L	100 pcs	10 packs	11584783

# Nitril LIOX®

Our first antimicrobial examination and protective glove

## Features

- Nitrile examination glove, indigo (DIN EN 455 standard)
- PPE Cat III
- powder-, protein- and latex-free
- micro-roughened
- elastic and tear-proof
- ambidextrous, with rolled edge
- suitable for use as a food contact material



## Usage

An examination and protective glove for universal use, for applications in any areas where the user is exposed to a risk of contamination. Suitable for emergency situations, rescue services and ICU staff. Suitable for highly infectious situations in hospitals and in general practice.

## Benefits:

- Protects against a wide variety of hospital-acquired infections
- Kills up to 99.999% of selected microorganisms
- Skin compatibility dermatologically tested

## Mechanism of action

An active dye incorporated into the outer surface of the glove material converts conventional oxygen into a compound known as singlet oxygen  $^1O_2$  as soon as the glove material comes into contact with both light and oxygen. The highly reactive singlet oxygen attacks the cell membranes of gram-positive bacteria (including MRSA, VRE, etc.), which ultimately leads to a dying-off of the affected microbes. Singlet oxygen is already in widespread use within medicine.



Size	REF	Pack contents	Shipping unit	PZN
XS	2284XS	200 pcs	10 packs	15377756
S	2284S	200 pcs	10 packs	15377762
M	2284M	200 pcs	10 packs	15377779
L	2284L	200 pcs	10 packs	15377785
XL	2284XL	180 pcs	10 packs	15377816





# Vinyl 2000 PF

Our vinyl examination glove

### Features

- latex-free examination glove (DIN EN 455 standard)
- with rolled edge
- ambidextrous
- vinyl is suitable for use as a latex substitute for users and patients with a diagnosed (type I) latex allergy.



### Usage

Suitable for daily patient care in retirement homes and hospitals, in general practice, for general healthcare work, and for handling foodstuffs that do not contain fat. Also suitable for use in the manufacturing sector. Thanks to its features and cost-effectiveness, it is ideal for applications with a very high consumption of disposable gloves.

### Powder-free, non-sterile

Size	REF	Pack contents	Shipping unit	PZN
XS	1251XS	100 pcs	10 packs	02243907
S	1251S	100 pcs	10 packs	02243913
M	1251M	100 pcs	10 packs	02243936
L	1251L	100 pcs	10 packs	02243942
XL	1251XL	100 pcs	10 packs	02243994

# Vinyl Blue Gen

Our blue vinyl examination glove

### Features

- latex-free examination glove (DIN EN 455 standard)
- beaded Cuff
- ambidextrous
- vinyl is suitable for use as a latex substitute for users and patients with a diagnosed (type I) latex allergy



### Usage

Suitable for daily patient care in retirement homes and hospitals, in general practice, for general healthcare work, and for handling foodstuffs that do not contain fat. Also suitable for use in the manufacturing sector. Thanks to its features and cost-effectiveness, it is ideal for applications with a very high consumption of disposable gloves.

### Powder-free, non-sterile

Size	REF	Pack contents	Shipping unit	PZN
XS	1255XS	100 pcs	10 packs	-
S	1255S	100 pcs	10 packs	-
M	1255M	100 pcs	10 packs	-
L	1255L	100 pcs	10 packs	-
XL	1255XL	100 pcs	10 packs	-

# Cotton glove

Our soft glove from pure cotton

## Features

- pure cotton glove
- air permeable
- anatomical fit



## Usage

For use by skincare practitioners. Especially when treating a patient with a very greasy or slow-absorbing ointment, for example for atopic dermatitis, for very dry skin or eczema and to cover gauze dressings (avoid direct contact with open wounds!). Also usable for the reprocessing of sterilised surgical instruments. Not a medical device.

## Non-sterile

Size	REF	Pack contents	Shipping unit	PZN
Size 6/XS	3611	10 pairs	30 packs	04354686
Size 7/XS	3612	10 pairs	30 packs	04354692
Size 8/S	3613	10 pairs	30 packs	04354700
Size 9/S	3614	10 pairs	30 packs	04354717
Size 10/M	3619	10 pairs	30 packs	04354568
Size 11/M	3620	10 pairs	30 packs	04354574
Size 12/L	3621	10 pairs	30 packs	04354580
Size 13/L	3622	10 pairs	30 packs	04354597
Size 14/L	3623	10 pairs	30 packs	04354657
Size 15/L	3624	10 pairs	30 packs	04354663

# Froggy™

Our hygienic glove made from soft polyethylene

## Features

- transparent PE hygienic glove
- powder-, protein- and latex-free
- in a handy glove dispenser for easy removal
- welded seams
- ambidextrous
- in women's and men's sizes



## Non-sterile

Size	REF	Pack contents	Shipping unit	PZN
Women	1192S	100 pcs	100 packs	06177745
Men	1192L	100 pcs	100 packs	06177751

## Usage

Economical to use for simple, day-to-day activities. Suitable for protecting the hands from general environmental dirt and soiling, and applications needing to meet hygienic standards. Usable in chemists, supermarkets, bakeries, the food industry and service stations. **IMPORTANT:** Use only once when filling motor vehicles with diesel fuel. Short-term protection against slight, superficial diesel soiling: this PE disposable glove is not classed as personal protective equipment (PPE) in accordance with Regulation (EU) 2016/425. Not a medical device.



# Latex finger cots

Our non-sterile disposable product with rolled edge and a smooth surface

### Features

- approved as safe for use as a short-term food contact material
- white natural latex finger cots, powder-, silicone- and thiuram-free
- only suitable for a minimal level of risk
- not suitable for activities involving hazardous, sharp-edged or aggressive materials



### Usage

For covering individual fingers to insert pessaries/suppositories, apply creams or conduct simple examinations. Typically used in the healthcare, hygiene and manufacturing sectors. No medical device.

### Non-sterile

Size	REF	Pack contents	Shipping unit	PZN
S	9776	100 pcs	1 pack	01333784
M	9777	100 pcs	1 pack	01333790
L	9778	100 pcs	1 pack	01333809
XL	9779	100 pcs	1 pack	01333815

# MyClean® easy touch

The reusable glove made from latex

### Features

- reusable glove made from latex, powder-free
- fully anatomical fit
- comfortable to wear thanks to cotton flocked lining
- special, easy-grip surface with honeycomb texture on glove palms
- suitable for food contact use



### Usage

Especially suitable for cleaning work, and for rinsing, washing, garden work and animal care. As a food contact material, the glove can also be used in the hospitality, catering and food processing industries.

### Non-sterile

Size	REF	Pack contents	Shipping unit
S	MAM74152	1 pair	100 packs
M	MAM74153	1 pair	100 packs
L	MAM74154	1 pair	100 packs
XL	MAM76286	1 pair	100 packs

# Copolymed®

## Our sterile copolymer examination glove

### Features

- sterile copolymer examination glove (DIN EN 455 standard)
- powder-, protein- and latex-free
- touch-sensitive
- tear-proof
- ambidextrous
- pressed on paper



### Usage

Suitable for short-term medical interventions, for sterile activities in the operating theatre and in wards, as well as in manufacturing. Can be used either as an underglove or an overglove if additional protection is required against soiling and fluids.

#### Sterile, individually packed

Size	REF	Pack contents	Shipping unit	PZN
S	8092S	100 pcs	12 packs	03896416
M	8092M	100 pcs	12 packs	03896221
L	8092L	100 pcs	12 packs	03896215

#### Sterile, packed in pairs

S	8091S	50 pairs	12 packs	11362864
M	8091M	50 pairs	12 packs	11362870
L	8091L	50 pairs	12 packs	11362887

# Gentle Skin® steril

## Our sterile latex examination glove

### Features

- sterile latex examination glove (DIN EN 455 standard)
- powder-free, low protein
- micro-roughened
- elastic and tear-proof
- ambidextrous
- packed singly or in pairs



### Usage

Suitable for use in any surgical application. Anywhere where wearers are working in a sterile environment. For smaller surgical interventions, changing dressings, taps, etc.

#### Sterile, individually packed

Size	REF	Pack contents	Shipping unit	PZN
S	8022S	50 pcs	8 packs	06402811
M	8022M	50 pcs	8 packs	06402828
L	8022L	50 pcs	8 packs	06402834

#### Sterile, packed in pairs

S	8021RSW	50 pairs	4 packs	02244479
M	8021RMW	50 pairs	4 packs	02244485
L	8021RLW	50 pairs	4 packs	02244491

# Nitril<sup>®</sup> sterile

## Our sterile nitrile examination glove

### Features

- sterile nitrile examination glove (DIN EN 455 standard)
- PPE Cat III
- micro-roughened
- elastic and tear-proof
- ambidextrous



### Usage

Usable for surgical applications. Anywhere where wearers are working in a sterile environment. Usable for smaller surgical interventions, changing dressings, taps, etc. Suitable for latex allergy sufferers with a diagnosed type I allergy.

### Sterile, packed in pairs

Size	REF	Pack contents	Shipping unit	PZN
S	8023S	50 pairs	4 packs	09544658
M	8023M	50 pairs	4 packs	09544641
L	8023L	50 pairs	4 packs	09544635
XL	8023XL	50 pairs	4 packs	11362976



# Our surgical gloves set the standards for protection and safety

During surgery, protecting medical staff and patients is of prime importance. The Meditrade surgical glove range includes a wide selection of glove models made from a variety of materials. Common to all gloves is their outstanding mechanical durability, thickness, grip, comfort and biological compatibility.



## Reference™ OP

### The powdered latex surgical glove

#### Features

- powdered latex surgical glove, natural (DIN EN 455 standard)
- PPE Cat III
- micro-roughened on the fingers
- lightly powdered
- anatomical fit modelled on natural hand posture
- left- and right-handed, with rolled edge
- packed in pairs
- sterile (ISO 11137-2 standard)

#### Usage

Suitable for use in all surgical areas that require sterility and anywhere where the user is working in a sterile environment with an identified risk of contamination. For highly infectious situations in hospitals, general practice, in laboratories and in the pharmaceutical industry. Also suitable for double-gloving.



#### Sterile, packed in pairs

Size	REF	Pack contents	Shipping unit	PZN
6	93216	50 pairs	6 packs	02243824
6.5	932165	50 pairs	6 packs	02243830
7	93217	50 pairs	6 packs	02243847
7.5	932175	50 pairs	6 packs	02243853
8	93218	50 pairs	6 packs	02243876
8.5	932185	50 pairs	6 packs	02243882
9	93219	50 pairs	6 packs	05030313



## Gentle Skin<sup>®</sup> Premium OP

The powder-free latex surgical glove

### Features

- powder-free latex surgical glove, natural (DIN EN 455 standard)
- PPE Cat III
- with polymer coating
- micro-roughened on the fingers
- anatomical fit modelled on natural hand posture
- left- and right-handed, with rolled edge
- packed in pairs
- sterile (ISO 11137-2 standard)

### Usage

Especially suited for vascular surgery, orthopaedics, abdominal surgery and in all surgical fields with sterile requirements. Anywhere there is an identified risk of contamination for the user in a sterile environment. For highly infectious situations in hospitals, general practice, in laboratories and in the pharmaceutical industry. Also suitable for double-gloving.



Sterile, packed in pairs

Size	REF	Pack contents	Shipping unit	PZN
5.5	902155W	50 pairs	4 packs	05030299
6	90216W	50 pairs	4 packs	02244114
6.5	902165W	50 pairs	4 packs	02244278
7	90217W	50 pairs	4 packs	02244290
7.5	902175W	50 pairs	4 packs	02244315
8	90218W	50 pairs	4 packs	02244338
8.5	902185W	50 pairs	4 packs	02244427
9	90219W	50 pairs	4 packs	02244433



## Gentle Skin<sup>®</sup> Superior OP

The powder-free latex surgical glove

### Features

- powder-free latex surgical glove, natural (DIN EN 455 standard)
- PPE Cat III
- anatomical fit modelled on natural hand posture
- left- and right-handed, with rolled edge
- packed in pairs
- sterile (ISO 11137-2 standard)

### Usage

Suitable for use in all surgical areas that require sterility and anywhere where the user is working in a sterile environment with an identified risk of contamination. For highly infectious situations in hospitals, general practice, in laboratories and in the pharmaceutical industry. Also suitable for double-gloving.



Sterile, packed in pairs

Size	REF	Pack contents	Shipping unit	PZN
5.5	904155	50 pairs	4 packs	-
6	904160	50 pairs	4 packs	-
6.5	904165	50 pairs	4 packs	-
7	904170	50 pairs	4 packs	-
7.5	904175	50 pairs	4 packs	-
8	904180	50 pairs	4 packs	-
8.5	904185	50 pairs	4 packs	-
9	904190	50 pairs	4 packs	-





## Gentle Skin® Smooth OP

The latex surgical glove with a smooth texture

### Features

- powder-free latex surgical glove, natural (DIN EN 455 standard)
- PPE Cat III
- with polymer coating
- anatomical fit modelled on natural hand posture
- left- and right-handed, with rolled edge
- packed in pairs
- sterile (ISO 11137-2)

### Usage

Especially suited for vascular surgery, orthopaedics, abdominal surgery and in all surgical fields with sterile requirements. Anywhere there is an identified risk of contamination for the user in a sterile environment. For highly infectious situations in hospitals, general practice, in laboratories and in the pharmaceutical industry. Also suitable for double-gloving.



Sterile, packed in pairs

Size	REF	Pack contents	Shipping unit	PZN
6	90266	50 pairs	4 packs	11154775
6.5	902665	50 pairs	4 packs	11154781
7	90267	50 pairs	4 packs	11154798
7.5	902675	50 pairs	4 packs	11154806
8	90268	50 pairs	4 packs	11154812
8.5	902685	50 pairs	4 packs	11154829
9	90269	50 pairs	4 packs	11154835



## Gentle Skin® Micro OP

The touch-sensitive latex surgical glove

### Features

- powder-free latex surgical glove, brown (DIN EN 455 standard)
- PPE Cat III
- micro-roughened on the fingers
- anatomical fit modelled on natural hand posture
- left- and right-handed, with rolled edge
- packed in pairs
- sterile (ISO 11137-2)

### Usage

The surgical glove for handling demanding applications where particularly high touch sensitivity is needed. Generally suitable for use in all surgical areas that require sterility and anywhere where the user is working in a sterile environment with an identified risk of contamination. For highly infectious situations in hospitals, general practice, in laboratories and in the pharmaceutical industry. Also suitable for double-gloving.



Sterile, packed in pairs

Size	REF	Pack contents	Shipping unit	PZN
6	90316	50 pairs	4 packs	06972460
6.5	903165	50 pairs	4 packs	06972589
7	90317	50 pairs	4 packs	06972595
7.5	903175	50 pairs	4 packs	06972626
8	90318	50 pairs	4 packs	06972655
8.5	903185	50 pairs	4 packs	06972661





## Maxitex® UG

The surgical latex underglove

### Features

- powder-free latex underglove, green (DIN EN 455 standard)
- can be combined with any primary surgical glove
- anatomical fit modelled on natural hand posture
- left- and right-handed, with rolled edge
- packed in pairs
- sterile (ISO 11137-2 standard)

### Usage

As a surgical underglove for use in double-gloving for greater safety during operations. Double-gloving increases the level of protection from contamination and infection by reducing the risk of perforation. Suitable for invasive interventions, especially for high-risk patients, as well as for emergency service use and interventions with a much higher risk of pathogen transfer.



### Sterile, packed in pairs

Size	REF	Pack contents	Shipping unit	PZN
6	1202001	50 pairs	4 packs	–
6.5	1202002	50 pairs	4 packs	–
7	1202003	50 pairs	4 packs	–
7.5	1202004	50 pairs	4 packs	–
8	1202005	50 pairs	4 packs	–
8.5	1202006	50 pairs	4 packs	–
9	1202007	50 pairs	4 packs	–



## Gentle Skin® Securitex® OP

The latex surgical glove with indicator function

### Features

- powder-free latex surgical glove, green and brown (DIN EN 455 standard)
- micro-roughened on the fingers
- anatomical fit modelled on natural hand posture
- left- and right-handed, with rolled edge
- 2x packed in pairs
- sterile (ISO 11137-2)


### Usage

Suitable for surgery and interventions with a high risk of infection that require a high level of protection, in hospitals and general practice, in labs and the pharmaceutical industry. For surgical applications that require sterility and anywhere where the user is working in a sterile environment with an identified risk of contamination.



### Sterile set, 2x packed in pairs

Size	REF	Pack contents	Shipping unit	PZN
6	90516	25 sets	4 packs	11349970
6.5	905165	25 sets	4 packs	03851224
7	90517	25 sets	4 packs	03851230
7.5	905175	25 sets	4 packs	03851247
8	90518	25 sets	4 packs	03851253
8.5	905185	25 sets	4 packs	03851276



## Gentle Skin® Isopretex® OP

The latex-free surgical glove made from polyisoprene

### Features

- powder-free polyisoprene surgical glove, beige (DIN EN 455 standard)
- PPE Cat III
- micro-roughened on the fingers
- anatomical fit modelled on natural hand posture
- Left- and right-handed, with rolled edge
- packed in pairs
- sterile (ISO 11137-2 standard)


### Usage

Latex-free glove for all surgical areas that require sterility, and anywhere where the user is working in a sterile environment with an identified risk of contamination. For highly infectious situations in hospitals, general practice, in laboratories and in the pharmaceutical industry. Also suitable for double-gloving.



Sterile, packed in pairs

Size	REF	Pack contents	Shipping unit	PZN
5.5	952255	50 pairs	4 packs	11190854
6	95226	50 pairs	4 packs	11190860
6.5	952265	50 pairs	4 packs	11190877
7	95227	50 pairs	4 packs	11190883
7.5	952275	50 pairs	4 packs	11190908
8	95228	50 pairs	4 packs	11190914
8.5	952285	50 pairs	4 packs	11190920
9	95229	50 pairs	4 packs	11190937



## Neopretex® OP

The latex-free surgical glove made from polychloroprene

### Features

- powder-free polychloroprene surgical glove, brown (DIN EN 455 standard)
- PPE Cat III
- micro-roughened on the fingers
- anatomical fit modelled on natural hand posture
- left- and right-handed, with rolled edge
- packed in pairs
- sterile (ISO 11137-2)

### Usage

Latex-free glove for all surgical areas that require sterility, and anywhere where the user is working in a sterile environment with an identified risk of contamination. For highly infectious situations in hospitals, general practice, in laboratories and in the pharmaceutical industry. Also suitable for double-gloving.



Sterile, packed in pairs

Size	REF	Pack contents	Shipping unit	PZN
6	95216	50 pairs	4 packs	02426890
6.5	952165	50 pairs	4 packs	02426909
7	95217	50 pairs	4 packs	02426915
7.5	952175	50 pairs	4 packs	02426921
8	95218	50 pairs	4 packs	02426938
8.5	952185	50 pairs	4 packs	02426944
9	95219	50 pairs	4 packs	02426950



# Finessis®

An innovative surgical glove programme, that reinvents

**PURITY, COMFORT and REACTIVE SAFETY.**

Made from **FLEXYLON™** – a high-performance elastomer



Flexylon™ is a compact poured film that contains no accelerators or cross-linkers.

This approach eliminates virtually all stitch holes and other micro-defects, thereby achieving the very highest quality standards (AQL = 0.1)!

## Reactive protection

Reactive antimicrobial intermediate layer


Outer layer made from FLEXYLON™



Inner layer made from FLEXYLON™



**Figure 1:** A sharp object presses against the outer layer, but without breaking through.



**Figure 2:** As pressure is increased, the thin membranes tear that separate the gel drops from one another. A larger, contiguous quantity of the antimicrobial, gel-like liquid is created.



**Figure 3:** As soon as the sharp object breaks through, it is enveloped by the solution and the viruses or bacteria are effectively neutralised.

# Finessis® Aegis®

Reactive protection

Trauma and high-risk Orthopaedics

FLEXYLON™ | Without accelerators | Synthetic surgical glove

Aegis® is a revolutionary protective system. It is a surgical glove that is equipped with a reaction-triggered disinfection system. Finessis® Aegis® was designed as a glove interior featuring an antimicrobial layer between the sheath, which consists of high-performance elastomers. The disinfection system is activated and effective only IF the glove is PERFORATED. In this case, the activated gel is released in the RIGHT QUANTITY and at the EXACT LOCATION where the perforation has occurred.

### Special characteristics:

- powder-free surgical glove (DIN EN 455 standard)
- PPE glove (DIN EN 374 standard)
- reaction-initiated disinfection system
- WITHOUT accelerators
- antimicrobial intermediate layer, which releases activated gel in the event of a perforation
- best AQL in the industry (AQL = 0.10)
- manufactured from FLEXYLON™ (layer technology)
- high electrical resistance
- anatomical fit

Size	REF	Pack contents	Shipping unit
6	90286	40 pairs	4 packs
6.5	902865	40 pairs	4 packs
7	90287	40 pairs	4 packs
7.5	902875	40 pairs	4 packs
8	90288	40 pairs	4 packs
8.5	902885	40 pairs	4 packs
9	90289	40 pairs	4 packs



COLOUR

Light green/white





# Finessis® Zero®

**Performance and safety**  
General surgery



**FLEXYLON™ | Without accelerators | Synthetic surgical glove**

Finessis® Zero® is the next generation of surgical gloves that are manufactured without any accelerators to achieve a new level of PURITY. By using FLEXYLON™ layer technology, the risk of a type I skin allergy is completely eliminated and the risk of a type IV skin allergy is reduced to a previously unattainable minimum.

**Special characteristics:**

- powder-free surgical glove (DIN EN 455)
- PPE glove (DIN EN 374)
- NO accelerators
- best AQL in the industry (AQL = 0.10)
- manufactured from FLEXYLON™ (layer technology)
- high electrical resistance
- anatomical fit

COLOUR

○ White




Size	REF	Pack contents	Shipping unit
6	90276	50 pairs	4 packs
6.5	902765	50 pairs	4 packs
7	90277	50 pairs	4 packs
7.5	902775	50 pairs	4 packs
8	90278	50 pairs	4 packs
8.5	902785	50 pairs	4 packs
9	90279	50 pairs	4 packs

Finessis®

# Finessis® Corium®

**Microsurgery and underglove**



**FLEXYLON™ | Without accelerators | Synthetic surgical glove**

Finessis® Corium® is the glove with a high tactility for double-gloving and microsurgery. The glove was developed to ensure improved touch sensitivity. Corium® works as an underglove, to achieve optimal comfort without losing grip or traction. When combined with the Finessis® Zero® for the perfect fit, the Finessis® double-gloving system presents an impermeable barrier for any surgical application.

**Special characteristics:**

- powder-free surgical glove (DIN EN 455)
- PPE glove (DIN EN 374)
- NO accelerators
- best AQL in the industry (AQL = 0.10)
- manufactured from FLEXYLON™ (layer technology)
- smooth surface
- high wearer comfort with improved tactility
- high electrical resistance
- anatomical fit

COLOUR

● Green




Size	REF	Pack contents	Shipping unit
6	90246	50 pairs	4 packs
6.5	902465	50 pairs	4 packs
7	90247	50 pairs	4 packs
7.5	902475	50 pairs	4 packs
8	90248	50 pairs	4 packs
8.5	902485	50 pairs	4 packs
9	90249	50 pairs	4 packs

Finessis®

# Improving safety in operating theatres and wards

## Features of Suavel® surgical masks



## Improving safety in operating theatres and wards

Our products make an important contribution to stopping infection pathways and the transmission of germs. Personnel, patients and visitors should feel safe and comfortable as they interact with each other. Wearer comfort, breathability and an adjustable fit are important aspects for suitable clothing.

For protective clothing in particular, it is not enough to simply rely on general or universal products: instead, clothing needs to meet the high standards in the medical field. When making Suavel face masks, for example, our techniques include an ultrasonic sealing method that prevents stitches coming loose in the final garment.

\* Suavel Antifluid: 4-ply



# Suavel® Protec Plus

# Suavel® Protec Plus

# Suavel® Protec

# Suavel® Protec Black

Our premium quality surgical mask with elastic ear loops, made in Germany

Our high filtration surgical mask with elastic ear loops

Our high quality surgical mask with elastic ear loops

Our high quality surgical mask with elastic ear loops in black

**Features**

- Type II/II R surgical mask (EN 14683 standard)
- bacterial filter efficiency  $\geq 98\%$
- high-performance filter with low respiratory resistance
- 3-ply, green exterior
- elastic ear loops
- genuine surgical non-woven
- cool and hypoallergenic
- disposable
- non-sterile

**Features**

- Type II surgical mask (EN 14683 standard)
- bacterial filter efficiency  $\geq 98\%$
- 3-ply, green exterior
- elastic ear loops
- genuine surgical non-woven
- hypoallergenic
- disposable
- non-sterile

**Features**

- Type II/II R surgical mask (EN 14683 standard)
- bacterial filter efficiency  $\geq 98\%$
- 3-ply, green or blue exterior
- elastic ear loops
- disposable
- non-sterile

**Features**

- Type II surgical mask (EN 14683 standard)
- bacterial filter efficiency  $\geq 98\%$
- 3-ply, black exterior
- elastic ear loops
- disposable
- non-sterile



**Usage**

For all surgical applications and for ward use, for the protection of both staff and patients.

**Usage**

For all surgical applications and for ward use, for the protection of both staff and patients.

**Usage**

For all surgical applications and for ward use, for the protection of both staff and patients.

**Usage**

For all surgical applications and forward use, for the protection of both staff and patients.

**Non-sterile**

REF	Type	Colour	Pack contents	Shipping unit	PZN
80-956	II	Green	50 pcs	20 packs	16760256
80-957	II R	White	50 pcs	20 packs	-

**Non-sterile**

REF	Colour	Pack contents	Shipping unit	PZN
80-955MP	Green	50 pcs	18 packs	02244605

**Non-sterile**

REF	Type	Colour	Pack contents	Shipping unit	PZN
80-900	II	Green	50 pcs	18 packs	00068771
80-901	II	Blue	50 pcs	20 packs	16731071
80-902	II R	Blue	50 pcs	40 packs	16731088

**Non-sterile**

REF	Colour	Pack contents	Shipping unit	PZN
80-904	Black	50 pcs	20 packs	-



# Suavel® Comfort Plus Suavel® Comfort

## Our tie-on surgical mask with superior filtering

### Features

- Type II surgical mask (EN 14683 standard)
- bacterial filter efficiency ≥98%
- 3-ply, blue exterior
- tie-on mask
- long ties
- genuine surgical non-woven
- hypoallergenic
- disposable
- non-sterile



### Usage

For all surgical applications and for ward use, for the protection of both staff and patients.

### Non-sterile

REF	Colour	Pack contents	Shipping unit	PZN
80-410M	Blue	50 pcs	12 packs	02244611

## Our high quality tie-on surgical mask

### Features

- Type II/II R surgical mask (EN 14683 standard)
- bacterial filter efficiency ≥98%
- 3-ply, blue exterior
- tie-on mask
- long ties
- disposable
- non-sterile



### Usage

For all surgical applications and for ward use, for the protection of both staff and patients.

### Non-sterile

REF	Type	Colour	Pack contents	Shipping unit	PZN
80-400	II	Blue	50 pcs	12 packs	00068788
80-402	II R	Blue	50 pcs	20 packs	–

# Suavel® Antifog

## Our surgical mask for spectacle wearers

### Features

- with effective anti-fogging insulation
- Type II surgical mask (EN 14683 standard)
- bacterial filter efficiency ≥98%
- 3-ply, blue exterior
- tie-on mask
- long ties
- disposable
- non-sterile



### Usage

Specially designed for spectacle wearers. For all surgical applications and for ward use, for the protection of both staff and patients.

### Non-sterile


REF	Colour	Pack contents	Shipping unit	PZN
80-470M	Blue	50 pcs	12 packs	00559457

# Suavel® Antifluid

## Our surgical mask with face shield

### Features

- High degree of protection against infection with integrated face shield for mouth, eyes and nose
- liquid stop layer
- Type II R surgical mask (EN 14683 standard)
- Bacterial filter efficiency ≥98%
- 4-ply, blue exterior
- tie-on mask
- long ties
- disposable
- non-sterile



### Usage

For highly infectious situations, protects user against splashes of blood/wound secretions, also for rinsing operations. For all surgical applications and for ward use.

### Non-sterile

REF	Colour	Pack contents	Shipping unit	PZN
80-455M	Blue	50 pcs	12 packs	00559517

# Suavel<sup>®</sup> Sensima

# Respima FFP2

## Our hypoallergenic mask for allergy sufferers

### Features

- hypoallergenic
- Type II surgical mask (EN 14683 standard) for users with sensitive skin
- bacterial filter efficiency ≥98%
- 3-ply, white exterior
- tie-on mask
- long ties
- disposable



Hypoallergenic laminated material



### Usage

For allergy sufferers and users with sensitive skin. For all surgical applications and for ward use, for the protection of both staff and patients.

### Non-sterile

REF	Colour	Pack contents	Shipping unit	PZN
80-440M	White	50 pcs	12 packs	11362700

## Disposable Protective Respirator Mask

### Features

- Personal protective equipment category III (EN 149:2001+A1:2009)
- 5 ply, white
- elastic ear loops
- disposable
- non-sterile
- without valve



### Usage

The filtering half mask is intended to protect against solid and oily particulates and liquids such as coal dust or cement dust.

### Non-sterile

REF	Type	Characteristics	Packaging	PZN
80-804	FFP2 NR, KSC001	packed individually	20 pcs. per package; 40 packages per carton	-

# Suavel<sup>®</sup> Apollo

## Astronaut-shaped surgical cap

### Features

- surgical cap in astronaut helmet style, made from a tear-resistant polypropylene non-woven
- with long ties and neck protection
- breathable
- with integrated sweatband



### Usage

For operations, for covering the hair and preventing excessive perspiration. Suitable for beard wearers.

### Non-sterile

REF	Colour	Pack contents	Shipping unit	PZN
70-131G	Green	100 pcs	10 packs	11362597

# Suavel<sup>®</sup> Astrid

## Our universal bouffant cap

### Features

- bouffant cap made from soft, light-weight non-woven polypropylene
- with double-sewn rubber drawstring
- breathable



### Usage

For use in hospitals and in the healthcare sector, for both nursing staff and patients. Usable anywhere where the hair has to be covered, also including applications in the food and catering sector.

### Suavel<sup>®</sup> Astrid, non-sterile

REF	Size	Colour	Pack contents	Shipping unit	PZN
70-121B	51 cm	Blue	100 pcs	10 packs	11362611
70-121G	51 cm	Green	100 pcs	10 packs	11362628
70-121W	51 cm	White	100 pcs	10 packs	11362605

### Suavel<sup>®</sup> Astrid XL, non-sterile

70-124G	60 cm	Green	100 pcs	10 packs	11362634
70-124W	60 cm	White	100 pcs	10 packs	11362640



## BeeSana® PP gown 35 g

**The high-performance**  
scrubs for personnel and visitors

**Features**

- scrubs made from 35 g/m<sup>2</sup> PP non-woven
- fastening ties at the neck and waist
- cuffs at the end of the sleeves
- liquid-repellent
- non-sterile
- disposable

**Usage**

Ward and visitor gowns to prevent the spread of germs in hospitals and in geriatric nursing; for personnel, patients and visitors.



**Non-sterile**

REF	Size	Colour	Pack contents	Shipping unit	PZN
3700	140 x 140 cm	Green	10 pcs	10 packs	09747462

## BeeSana® PP gown 23 g

**The industry-standard**  
protective gown for personnel and visitors

**Features**

- protective gown made from 23 g/m<sup>2</sup> PP non-woven
- fastening ties at the neck and waist
- sleeves come with either sewn-in, elastic rubber ties or with cuffs
- lightweight material
- moisture-repellent
- non-sterile
- disposable

**Usage**

Ward and visitor gowns to prevent the spread of germs in hospitals and in geriatric nursing; for personnel, patients and visitors.



**Non-sterile, with elastic rubber ties**

REF	Size	Colour	Pack contents	Shipping unit	PZN
3684	110 x 140 cm	White	10 pcs	10 packs	11509391
3685	110 x 140 cm	Green	10 pcs	10 packs	11509416
3699	120 x 140 cm	Yellow	10 pcs	10 packs	09747479
3729	130 x 150 cm	White	10 pcs	10 packs	11527609

**Non-sterile, with elastic cuffs**

4685	110 x 140 cm	Green	10 pcs	10 packs	11363036
4699	120 x 140 cm	Yellow	10 pcs	10 packs	11363042

## BeeSana® PP/PE gown 26 g

**The waterproof**  
semi-laminated protective gown for  
personnel and visitors

**Features**

- protective gown made from 26 g/m<sup>2</sup> PP non-woven with PE coating covering front/arm areas
- fastening ties at the neck and waist
- cuffs at the end of the sleeves
- waterproof across the coated areas
- non-sterile
- disposable

**Usage**

Ward and visitor gowns for applications requiring a high level of protection to prevent the spread of germs in hospitals and in geriatric nursing; for personnel, patients and visitors.



**Non-sterile**

REF	Size	Colour	Pack contents	Shipping unit	PZN
4698	136 x 136 cm	Yellow	10 pcs	10 packs	11642842
4698-XL	150 x 170 cm	Yellow	10 pcs	10 packs	17376146

## BeeSana® SMS gown 40 g

**High protection**  
Protective gown offering a high level of  
protection against liquids and particulates

**Features**

- protective gown made from soft, 3-ply SMS composite material
- combined layers provide a highly effective barrier against fluids, germs and light splashes of chemicals
- fastening ties at the neck and waist and cuffs at the end of the sleeves
- non-sterile
- disposable

**Usage**

Suitable for highly infectious situations. To prevent the spread of germs on emergency and isolation wards, in care homes, and in hospitals, for personnel, patients and visitors.



**Non-sterile**

REF	Size	Colour	Pack contents	Shipping unit	PZN
3597	115 x 150 cm	Blue	10 pcs	5 packs	11652154



# BeeSana<sup>®</sup> SMS gown 15 g

# Suavel<sup>®</sup> Yeti Coverall

**The lightweight**  
protective gown  
for personnel and visitors

**Features**

- protective gown made from soft, 3-ply 15 g/m<sup>2</sup> SMS composite material
- fastening ties at the neck and waist and cuffs at the end of the sleeves
- non-sterile
- disposable

**Usage**

For highly infectious situations. To prevent the spread of germs on emergency and isolation wards, in care homes, and in hospitals, for personnel, patients and visitors.



Protective suit

**Features**

- protective suit made from polypropylene with an air-permeable coating, offering protection against harmful dusts (type 5) and splashes of liquid (type 6)
- elastic rubber ties at wrists and ankles
- concealed zip closure at front
- non-sterile
- disposable

**Usage**

Category III protective suit for use in hospitals, e.g. in the ongoing pandemic situation. Offers protection against harmful dusts and splashes of liquid.



**Non-sterile**

REF	Size	Height	Chest	Packing unit
4798	M	168 - 176 cm	92 - 100 cm	1 pc./box 50 packs
4799	L	174 - 182 cm	100 - 108 cm	1 pc./box 50 packs
4800	XL	182 - 188	108 - 116 cm	1 pc./box 50 packs
4801	XXL	188 - 194	116 - 124 cm	1 pc./box 50 packs

**Non-sterile**

REF	Size	Colour	Pack contents	Shipping unit	PZN
3598	140 x 140 cm	Blue	10 pcs	10 packs	11518680

# PE disposable apron

Our PE apron in various models

**Features**

- disposable LDPE apron
- tested to the EN 14126 standard
- with closing ties at the waist to ensure a good fit
- lightweight material
- hard-wearing and moisture-repellent
- in a handy pouch or on hanger block

**Usage**

Protection for clothing in nursing/care activities, and in all areas with stringent hygiene requirements, including large caterers and hospitals.



**Non-sterile, in pouch, for individual removal**

REF	Size	Thickness	Colour	Pack contents	Shipping unit	PZN
3591	140 x 75 cm	30 µm	Blue	50 pcs	10 packs	11857643
3690	125 x 81 cm	16 µm	White	100 pcs	10 packs	04353379
3691	140 x 81 cm	16 µm	White	100 pcs	10 packs	04353422
3692	160 x 81 cm	16 µm	White	100 pcs	10 packs	04353474
3696	117 x 69 cm	16 µm	White	100 pcs	10 packs	-

**Non-sterile, on hanger block, for tearing-off individually**

3687	125 x 81 cm	16 µm	White	50 pcs	10 packs	04353385
3688	140 x 81 cm	16 µm	White	50 pcs	10 packs	04353468
3694	160 x 81 cm	16 µm	White	50 pcs	10 packs	04353480
3697	117 x 69 cm	16 µm	White	200 pcs	5 rolls	-

## Suavel® OP PRO Standard SMMS surgical gown, sterile

### Features

- surgical gowns made from 35 g/m<sup>2</sup> SMMS composite material (DIN EN 13795 standard)
- combined layers provide a highly effective barrier against fluids, germs and light chemical spills
- fastening ties with sterile card at the waist and cuffs at the end of the sleeves; hook-and-pile fastener at neck
- SMS wrapping cloth plus two paper towels
- counting cards on the peel pack
- sterile and disposable

### Usage

During surgical interventions, wearing the Suavel® OP PRO Standard prevents the transfer of particulates and infectious material from the surgical team to the surgical wound, and vice versa.

### Sterile

REF	Size	Colour	Pack contents	Shipping unit	PZN
5541	M 120 x 140 cm	Blue	1 pc	30 packs	15295662
5542	L 130 x 150 cm	Blue	1 pc	30 packs	15295679
5543	XL 150 x 160 cm	Blue	1 pc	30 packs	15295685
5544	XXL 170 x 160 cm	Blue	1 pc	30 packs	15295691



## Suavel® OP PRO Expert SMMS surgical gown, sterile, reinforced

### Features

- surgical gowns made from 43 g/m<sup>2</sup> SMMS composite material (DIN EN 13795 standard)
- combined layers provide a highly effective barrier against fluids, germs and light chemical spills
- fastening ties with sterile card at the waist and cuffs at the end of the sleeves; hook-and-pile fastener at neck
- reinforced chest, abdominal area and sleeves
- SMS wrapping cloth plus two paper towels
- counting cards on the peel pack
- sterile and disposable

### Usage

During surgical interventions, wearing the Suavel® OP PRO Expert prevents the transfer of particulates and infectious material from the surgical team to the surgical wound, and vice versa.

### Sterile

REF	Size	Colour	Pack contents	Shipping unit	PZN
5551	M 120 x 140 cm	Blue	1 pc	30 packs	15295716
5552	L 130 x 150 cm	Blue	1 pc	30 packs	15295722
5553	XL 150 x 160 cm	Blue	1 pc	30 packs	15295739
5554	XXL 170 x 160 cm	Blue	1 pc	30 packs	15295745





# Wound treatment and dressings – the pillars that shoulder the weight of medical care

Alongside our wound dressings, we aim to offer the ideal range of products to treat every type of wound at each stage of recovery. After all, a wound dressing alone cannot hope to satisfy every aspect or goal for treatment.

Our product range includes all non-sterile and sterile gauze/non-woven swabs as well as absorbent dressing pads. In the field of therapeutic dressings, we aim to offer a comprehensive set of products for an array of applications that include dressing, supporting, and relieving the wound.

Patient safety is at the forefront of everything we do, and all our sterile products are fully compliant with the relevant safety standards. Features here include woven-in x-ray contrast threads and the counting check card.



# BeeSana® Gauze cutting

For individual wound coverings

## Features

- pre-cut gauze (EN 14079 standard)
- 20-thread
- absorbent



## Usage

For individual wound covering applications. For the external wound care of exuding and/or contaminated wounds. For blood and exudate absorption.

### Non-sterile

REF	Size (width x length)	Pack contents	Shipping unit	PZN
1330	20 x 20 cm (approx.)	500 pcs	10 packs	–
1335	30 x 40 cm (approx.)	500 pcs	6 packs	–

# BeeSana® Medical gauze

For individual wound coverings

## Features

- medical gauze (EN 14079 standard)
- 20-thread

## Usage

For preparing individual bandages and wound coverings.



### Non-sterile

80 cm wide, 1-ply and 4-ply

REF	Size (width x length)	Pack contents	Shipping unit	PZN
1301	80 cm x 40 m (approx.)	5 pcs	2 packs	04852674

### Double sterile, with self-adhesive counting cards

80 cm wide, 5-ply, folded to 16 x 20 cm, 3 sterile pcs per pack

REF	Length	Pack contents	Shipping unit	PZN
4179	100	8 x 3 pcs	2 packs	13581004

# BeeSana® Gauze swab

Our standard gauze swab with selvedge cut edges

## Features

- gauze swab, medical gauze (EN 14079 standard)
- selvedge cut edges, 100% cotton
- folded, 17-thread
- absorbent



## Usage

For general-purpose wound care of exuding wounds and for absorbing bodily fluids. For applying ointments. For smaller-scale outpatient surgical interventions, in wards and for first aid applications.

### Non-sterile

	REF	Size (width x length)	Pack contents	Shipping unit	PZN	
<b>8-ply</b>	1000	5 x 5 cm (approx.)	100 pcs	10 packs	04852705	
	1006	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	04852711	
	1012	10 x 10 cm (approx.)	100 pcs	10 packs	04852728	
	1053	10 x 20 cm (approx.)	100 pcs	5 packs	01330633	
	10 pcs, pre-counted For set packers	2712	10 x 10 cm (approx.)	200 x 10 pcs	1 box	–
<b>12-ply</b>	1002	5 x 5 cm (approx.)	100 pcs	10 packs	04852734	
	1008	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	04852740	
	1014	10 x 10 cm (approx.)	100 pcs	10 packs	04852757	
	1038	10 x 13 cm (approx.)	100 pcs	10 packs	04852763	
	1023	10 x 20 cm (approx.)	100 pcs	10 packs	04852786	
		10 pcs, pre-counted For set packers	2708	7.5 x 7.5 cm (approx.)	200 x 10 pcs	1 box
		2714	10 x 10 cm (approx.)	200 x 10 pcs	1 box	–
		2723	10 x 20 cm (approx.)	100 x 10 pcs	1 box	–
<b>16-ply</b>	1004	5 x 5 cm (approx.)	100 pcs	10 packs	04852792	
	1010	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	04852800	
	1016	10 x 10 cm (approx.)	100 pcs	10 packs	04852817	
	1052	10 x 13 cm (approx.)	100 pcs	10 packs	04852823	
	1025	10 x 20 cm (approx.)	100 pcs	5 packs	11849589	
<b>32-ply</b>	1018	10 x 10 cm (approx.)	50 pcs	10 packs	–	
<b>64-ply</b>	20-thread	1235	10 x 13 cm (approx.)	25 pcs	8 packs	11849626

# BeeSana® Gauze swab

Our standard gauze swab with selvedge cut edges

### Features

- gauze swab, medical gauze (EN 14079 standard)
- selvedge cut edges, 100% cotton
- folded, 17-thread
- absorbent



### Single sterile

		REF	Size (width x length)	Pack contents	Shipping unit	PZN
<b>8-ply</b>	Pack of 2 pcs	1043	5 x 5 cm (approx.)	25 x 2 pcs	27 packs	04852846
		1044	7.5 x 7.5 cm (approx.)	25 x 2 pcs	16 packs	04852852
	Pack of 5 pcs	1045	10 x 10 cm (approx.)	25 x 2 pcs	16 packs	04852869
		1046	10 x 20 cm (approx.)	25 x 2 pcs	8 packs	04852875
		4600	5 x 5 cm (approx.)	80 x 5 pcs	8 packs	13458367
		4601	7.5 x 7.5 cm (approx.)	80 x 5 pcs	8 packs	01330662
Pack of 10 pcs	1062	10 x 10 cm (approx.)	35 x 5 pcs	6 packs	01330679	
	4606	10 x 10 cm (approx.)	50 x 10 pcs	4 packs	12348949	
	<b>12-ply</b>	Pack of 2 pcs	4130	5 x 5 cm (approx.)	25 x 2 pcs	24 packs
4131			7.5 x 7.5 cm (approx.)	25 x 2 pcs	18 packs	07216335
4132			10 x 10 cm (approx.)	25 x 2 pcs	12 packs	07216341
4133			10 x 20 cm (approx.)	25 x 2 pcs	6 packs	07216358
Pack of 5 pcs		1064	5 x 5 cm (approx.)	30 x 5 pcs	12 packs	01330722
		1065	7.5 x 7.5 cm (approx.)	30 x 5 pcs	6 packs	01330739
	1063	10 x 10 cm (approx.)	30 x 5 pcs	6 packs	01330745	
Pack of 10 pcs	1066	10 x 20 cm (approx.)	15 x 5 pcs	6 packs	01330751	
	4616	5 x 5 cm (approx.)	50 x 10 pcs	8 packs	12348955	
	4617	7.5 x 7.5 cm (approx.)	30 x 10 pcs	8 packs	12349044	
	1026	10 x 10 cm (approx.)	40 x 10 pcs	4 packs	12348694	
Pack of 20 pcs	4619	10 x 20 cm (approx.)	20 x 10 pcs	4 packs	12349015	
	4622	10 x 10 cm (approx.)	20 x 20 pcs	4 packs	12349021	
Pack of 30 pcs	4640	7.5 x 7.5 cm (approx.)	25 x 30 pcs	4 packs	12349185	
<b>16-ply</b>	Pack of 4 pcs	1092	7.5 x 7.5 cm (approx.)	64 x 4 pcs	8 packs	12348783
<b>32-ply</b>	Pack of 2 pcs	4628	10 x 10 cm (approx.)	30 x 2 pcs	8 packs	12349038



### Single sterile, rolled gauze dressing pad

		REF	Size (width x length)	Pack contents	Shipping unit	PZN
<b>16-ply</b>	Pack of 2 pcs	1146	20 x 20 cm (approx.)	20 x 2 pcs	4 packs	13582587



### Double sterile, with self-adhesive counting cards

<b>16-ply</b>	Pack of 10 pcs	1058	10 x 10 cm (approx.)	25 x 10 pcs	4 packs	13582601
	Pack of 60 pcs	4644	10 x 10 cm (approx.)	8 x 60 pcs	2 packs	12349216
<b>32-ply</b>	Pack of 5 pcs	4636	12.5 x 20 cm (approx.)	12 x 5 pcs	4 packs	12349239

# BeeSana® Slit gauze swab

Our gauze swabs with Y-slit

## Features

- gauze swab, medical gauze (EN 14079 standard)
- with Y-slit
- selvedge cut edges
- folded, 17-thread



## Usage

For wound care for drainage tubes, cannulas, tracheotomies and tractions.

### Non-sterile

	REF	Size (width x length)	Pack contents	Shipping unit	PZN
<b>8-ply</b>	4291	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	12491084
<b>12-ply</b>	4192	5 x 5 cm (approx.)	100 pcs	10 packs	12491090
	4122	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	12491109
	4160	10 x 10 cm (approx.)	100 pcs	10 packs	12491115
<b>16-ply</b>	4159	5 x 5 cm (approx.)	100 pcs	10 packs	12491121



### Single sterile, pack of 2 pcs

<b>8-ply</b>	4158	10 x 10 cm (approx.)	100 x 2 pcs	8 packs	03870954
<b>12-ply</b>	4115	5 x 5 cm (approx.)	25 x 2 pcs	16 packs	01916669
	4116	7.5 x 7.5 cm (approx.)	25 x 2 pcs	16 packs	01916652
	4117	10 x 10 cm (approx.)	25 x 2 pcs	16 packs	01916646

# BeeSana® Non-woven swab

Our non-woven swab for external wound care

## Features

- viscose and polyester non-woven swab
- absorbent
- virtually no fraying



## Usage

For general-purpose wound care of exuding wounds and for absorbing bodily fluids.

For applying ointment in smaller-scale outpatient surgical interventions, in wards and for first aid applications.

### Non-sterile

	REF	Size (width x length)	Pack contents	Shipping unit	PZN	
<b>4-ply</b> 30 g	1789	5 x 5 cm (approx.)	100 pcs	10 packs	01450254	
	1790	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	01450225	
	1791	10 x 10 cm (approx.)	100 pcs	10 packs	01450231	
	1792	10 x 20 cm (approx.)	100 pcs	5 packs	01445773	
10 pcs, pre-counted For set packers	2791	10 x 10 cm (approx.)	200 x 10 pcs	1 box	–	
	40 g	1124	5 x 5 cm (approx.)	100 pcs	10 packs	04852881
		1125	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	04852898
		1126	10 x 10 cm (approx.)	100 pcs	10 packs	04852906
	Large box	1127	10 x 10 cm (approx.)	100 pcs	76 packs	–
		11261	10 x 20 cm (approx.)	100 pcs	5 packs	04852912
	<b>6-ply</b> 30 g	1796	5 x 5 cm (approx.)	100 pcs	10 packs	01330886
1797		7.5 x 7.5 cm (approx.)	100 pcs	10 packs	01330900	
1798		10 x 10 cm (approx.)	100 pcs	10 packs	01330917	
1799		10 x 20 cm (approx.)	100 pcs	5 packs	01330923	
10 pcs, pre-counted For set packers		2798	10 x 10 cm (approx.)	200 x 10 pcs	1 box	–
	2799	10 x 20 cm (approx.)	100 x 10 pcs	1 box	–	

### Single sterile

<b>4-ply</b> 30 g	Pack of 2 pcs	1137	5 x 5 cm (approx.)	25 x 2 pcs	27 packs	03036843
		1138	7.5 x 7.5 cm (approx.)	25 x 2 pcs	16 packs	03036895
	Pack of 5 pcs	1139	10 x 10 cm (approx.)	25 x 2 pcs	16 packs	03036903
		1140	10 x 20 cm (approx.)	25 x 2 pcs	8 packs	03036949
Pack of 10 pcs	4376	7.5 x 7.5 cm (approx.)	100 x 5 pcs	8 packs	12350768	
	4377	10 x 10 cm (approx.)	60 x 5 pcs	8 packs	12350774	
	1219	10 x 10 cm (approx.)	25 x 10 pcs	8 packs	12350691	
<b>6-ply</b> 30 g	Pack of 2 pcs	1800	5 x 5 cm (approx.)	25 x 2 pcs	27 packs	01330834
		1801	7.5 x 7.5 cm (approx.)	25 x 2 pcs	16 packs	01330840
		1802	10 x 10 cm (approx.)	25 x 2 pcs	16 packs	01330857
		1803	10 x 20 cm (approx.)	25 x 2 pcs	8 packs	03036955
	Pack of 5 pcs	4225	5 x 5 cm (approx.)	120 x 5 pcs	8 packs	12350662
		4226	7.5 x 7.5 cm (approx.)	70 x 5 pcs	8 packs	12350679
		4227	10 x 10 cm (approx.)	40 x 5 pcs	8 packs	12350685



# BeeSana® Slit non-woven swab

Our non-woven swab with Y-slit

**Features**

- viscose and polyester non-woven slit swab
- with Y-slit
- absorbent
- virtually no fraying



**Usage**

For wound care for drainage tubes, cannulas, tracheotomies and tractions.

**Non-sterile**

	REF	Size (width x length)	Pack contents	Shipping unit	PZN
<b>4-ply, 30 g</b>	1793	5 x 5 cm (approx.)	100 pcs	10 packs	01330768
	1794	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	01330774
	1795	10 x 10 cm (approx.)	100 pcs	10 packs	01330780



**Single sterile, pack of 2 pcs**

<b>4-ply, 30 g</b>	1232	5 x 5 cm (approx.)	25 x 2 pcs	27 packs	01364684
	1234	7.5 x 7.5 cm (approx.)	25 x 2 pcs	16 packs	01364678
	1233	10 x 10 cm (approx.)	25 x 2 pcs	16 packs	01364661
<b>6-ply, 30 g</b>	4105	10 x 10 cm (approx.)	25 x 2 pcs	28 packs	13507290

# ABE® drainal

Our slit/hole dressing pad with metallic coating

**Features**

- non-woven slit dressing pad with hole
- metallic coating lessens the risk of the material adhering to the wound
- hole diameter: 5.5–11 mm (approx.)



**Usage**

For external wound care with drainage tubes and infusion catheters. Place the metallic coating against the skin.

**Single sterile, pack of 1 pc.**

REF	Size (width x length)	Pack contents	Shipping unit	PZN
4118	6 x 7 cm (approx.)	50 x 1 pc.	6 packs	01333867

# ABE® tracheal

Our slit/hole dressing pads with metal coating

**Features**

- non-woven slit dressing pad with hole
- metallic coating lessens the risk of the material adhering to the wound
- sterile packed
- hole diameter: 13–30 mm (approx.)



**Usage**

Slit/hole dressing pad for would care with catheters, stomata and drainage tubes. The hole in the pad also has a slit, so the ABE® Tracheal can create a perfect seal around a drainage tube or stoma, and the exudate can be quickly absorbed by the the non-woven dressing pad. Place the metallic coating against the skin.

**Single sterile, pack of 1 pc.**

REF	Size (width x length)	Pack contents	Shipping unit	PZN
4119	8 x 10 cm (approx.)	50 x 1 pc.	6 packs	—

# BeeSana® Absorbent dressing pad

Our absorbent pad for superior fluid absorption

## Features

- wound-friendly absorbent dressing pad
- absorbent core made from 100% cellulose
- fully enveloped in non-woven material
- blue top surface for identifying the reverse side and as a hydrophobic barrier layer to avoid soaking through, and protect underwear and clothing from contamination

## Usage

For wound care of exuding wounds in all areas of health-care. Absorbent dressing pads are also ideal for padding wounds.



### Non-sterile, in polybag

REF	Size (width x length)	Pack contents	Shipping unit	PZN
1227	10 x 10 cm (approx.)	25 pcs	16 packs	01363816
1228	10 x 20 cm (approx.)	25 pcs	8 packs	07607202
1224	10 x 30 cm (approx.)	25 pcs	15 packs	12146619
1225	15 x 20 cm (approx.)	25 pcs	12 packs	12146625
1231	15 x 25 cm (approx.)	25 pcs	7 packs	00137727
1229	20 x 20 cm (approx.)	25 pcs	12 packs	08696337
1244	20 x 30 cm (approx.)	25 pcs	5 packs	12146631
1230	20 x 40 cm (approx.)	25 pcs	6 packs	08442499

### Single sterile, shrink-wrapped pack of 1 pc.

1354	10 x 10 cm (approx.)	60 x 1 pc.	10 packs	08416059
1355	10 x 20 cm (approx.)	30 x 1 pc.	12 packs	06193514
1388	15 x 25 cm (approx.)	30 x 1 pc.	8 packs	12146602
1356	20 x 20 cm (approx.)	30 x 1 pc.	6 packs	08416042
1357	20 x 40 cm (approx.)	25 x 1 pc.	4 packs	07737904

# BeeSana® Eye pad

Our eye dressing pad with cushioning effect

## Features

- oval dressing pad with all-round seal
- absorbent inner filling



## Usage

For absorbing fluids from the eye. For padding and for eye bandages. Important: For injuries to the external eye and for post-operative applications, a sterile eye dressing pad/wound covering must be used. The oval design is ideal for fitting the eye socket.

### Non-sterile

REF	Size (width x length)	Pack contents	Shipping unit	PZN
4195	5.6 x 7.2 cm (approx.)	50 pcs	10 packs	02214099
4155	5.6 x 7.2 cm (approx.)	500 pcs	1 pack	12348872
4194	7 x 8.5 cm (approx.)	50 pcs	10 packs	02435297

### Single sterile, shrink-wrapped pack of 1 pc.

4196	5 x 7 cm (approx.)	25 x 1 pc.	9 packs	02214107
4156	7 x 8.5 cm (approx.)	25 x 1 pc.	6 packs	02214136

# ABE® nose bandage

Our flexible nasal sling

## Features

- bi-elastic nose bandage
- with elastic ear loops or an elastic strap
- with integrated dressing pad in centre



## Usage

For post-operative ENT applications with exudate.

Worn with ear loops or by being tied at the back of the head.

### Non-sterile, with elastic strap

REF	Pack contents	Shipping unit	PZN
4039	20 pcs	10 packs	01333749

### Non-sterile, with ear loops

REF	Pack contents	Shipping unit	PZN
9814	20 pcs	10 packs	11885562
9818	500 pcs	1 pack	13579533

# BeeSana® Gauze balls

Our super-absorbent gauze balls

## Features

- gauze balls, medical gauze (EN 14079 standard)
- innovative wrapping system from just one piece of pre-cut gauze
- inner-facing cut edges
- 20-thread
- absorbent



## Usage

For cleaning wounds and staunching bleeding. Available in various swab sizes to suit each wound. Non-sterile,

### Non-sterile, in polybag

	REF	Size (width x length)	Shape	Pack contents	Shipping unit	PZN
	1102	10 x 10 cm (approx.)	Round	500 pcs	8 packs	01364017
	4460	15 x 16 cm (approx.)	Round	500 pcs	2 packs	00822216
	4461	20 x 20 cm (approx.)	Round	250 pcs	16 packs	00822239
	1106	20 x 20 cm (approx.)	Thimble	250 pcs	16 packs	00822222
	4462	25 x 26 cm (approx.)	Round	250 pcs	10 packs	00822251
	1109	25 x 26 cm (approx.)	Thimble	250 pcs	10 packs	00822245
	1348	30 x 30 cm (approx.)	Round	100 pcs	20 packs	11862578
	1113	30 x 40 cm (approx.)	Round	100 pcs	15 packs	01363940
	1116	40 x 50 cm	Round	100 pcs	10 packs	00822268

### Single sterile, in peel pack

Pack of 2 pcs	Plum-sized	4356	20 x 20 cm (approx.)	Round	100 x 2 pcs	8 packs	13506801
Pack of 3 pcs	Plum-sized	4357	20 x 20 cm (approx.)	Round	80 x 3 pcs	8 packs	00137696
	Egg-sized	4308	25 x 26 cm (approx.)	Round	60 x 3 pcs	8 packs	13506818
	Fist-sized	4085	45 x 45 cm (approx.)	Round	35 x 3 pcs	4 packs	13506824
2+2 pcs per pack Packed in soft peel pack with 2 separate compartments (tear-off perforation)	Plum-sized	1084	20 x 20 cm (approx.)	Round	50x (2+2) pcs	8 packs	01364129
	Plum-sized	1164	20 x 20 cm (approx.)	Thimble	50x (2+2) pcs	8 packs	01364106
2+3 pcs per pack Packed in soft peel pack with 2 separate compartments (tear-off perforation)	Walnut-sized	1074	15 x 16 cm (approx.)	Round	45x (2+3) pcs	8 packs	13506830
	Plum-sized	1131	20 x 20 cm (approx.)	Round	40x (2+3) pcs	8 packs	01364135
	Plum-sized	1189	20 x 20 cm (approx.)	Thimble	40x (2+3) pcs	8 packs	01364075
	Egg-sized	1141	25 x 26 cm (approx.)	Round	35x (2+3) pcs	8 packs	–
	XX-sized	4216	30 x 40 cm (approx.)	Round	35x (2+3) pcs	4 packs	13506847
Pack of 5 pcs	Plum-sized	1184	20 x 20 cm (approx.)	Round	50 x 5 pcs	8 packs	13506853
	Egg-sized	4209	25 x 26 cm (approx.)	Round	40 x 5 pcs	8 packs	13506876
	Egg-sized	1165	25 x 26 cm (approx.)	Thimble	40 x 5 pcs	8 packs	01364069
	XX-sized	4211	30 x 40 cm (approx.)	Round	40 x 5 pcs	4 packs	13506882
Pack of 10 pcs	Plum-sized	4215	20 x 20 cm (approx.)	Round	30 x 10 pcs	8 packs	07393706





## BeeSana® Pointed swabs

Our pointed swabs made from pure cotton

### Features

- gauze swab, 100% cotton (EN 14079 standard)
- 20-thread
- tapered, formed from just one piece of gauze



### Usage

For blood and exudate absorption, especially in ENT applications. For staunching bleeding or as a tamponade.

Double sterile, with counting cards

	REF	Size (unrolled) (width x length)	Length	Pack contents	Shipping unit	PZN
Pack of 10 pcs	4170	Large, 15 x 16 cm (approx.)	13.5 cm (approx.)	40 x 10 pcs	8 packs	–

## ABE® ribbon gauze

Our tamponade made from pure cotton

### Features

- medical gauze (EN 14079 standard)
- 100% cotton
- woven edges
- white
- rolled



### Usage

As a tamponade in natural bodily orifices.

Single sterile, pack of 1 pc.

	REF	Size (width x length)	Pack contents	Shipping unit	PZN
24-thread	4300	1 cm x 5 m (approx.)	100 x 1 pc.	8 packs	13574814
	4301	2 cm x 5 m (approx.)	70 x 1 pc.	8 packs	13574820

## BeeSana® Cellulose swabs

Our cellulose swab on a roll

### Features

- swab made from highly bleached cellulose (German pharmacopoeia (DAB) standard)
- ready-to-use
- pre-cut



### Usage

For local disinfection before injections and infusions.

Non-sterile

	REF	Size (width x length)	Pack contents	Shipping unit	PZN
	1190	4 x 5 cm (approx.)	2 rolls x 500 pcs	16 packs	07210048
	1278	8 x 10 cm (approx.)	2 rolls x 125 pcs	10 packs	01364253
Reduced bioburden	1177	4 x 5 cm (approx.)	1 roll x 500 pcs	20 packs	01331070
Cellulose swab box	9871		1 pc.	1 pack	01363756

## Alcohol swabs

Our ready-to-use swabs for skin disinfection

### Features

- non-woven alcohol swabs
- impregnated with 70% isopropyl alcohol
- ready-to-use

### Usage

For disinfecting skin before subcutaneous injections and capillary blood sampling.



Non-sterile

	REF	Size (width x length)	Pack contents	Shipping unit	PZN
	4465	Approx. 3 x 3 cm (folded) Approx. 6 x 3 cm (unfolded)	100 pcs	60 packs	09391893

# Cotton buds

## Features

- various designs and variations for different types of applications
- sterile and non-sterile formats
- various head sizes
- wood or flexible plastic variants



## Usage

For a wide range of uses in medical applications, such as cleaning, applying liquid substances and ointments as well as for laboratory work.

### Non-sterile, in ziploc bag

REF	Length	Head	Handle	Pack contents	Shipping unit	PZN
4046	7.5 cm	Small head	Plastic	10 x 200 pcs	10 packs	11857620
1260	15 cm	Small head	Wood	10 x 100 pcs	20 packs	06177225
1261	15 cm	Medium head	Wood	10 x 100 pcs	10 packs	06177231
1263	15 cm	Large head	Wood	10 x 50 pcs	10 packs	01363762
1268	15 cm	Large head	Plastic	10 x 50 pcs	10 packs	11857494

### Single sterile, 2 pcs per pack

1262	15 cm	Small head	Wood	100 x 2 pcs	6 packs	06177254
1264	15 cm	Large head	Wood	100 x 2 pcs	3 packs	01363779
1270	15 cm	Large head	Plastic	100 x 2 pcs	3 packs	13711168

# Wooden tongue depressor

Our tongue depressor for examinations and for applying ointment



### Non-sterile

REF	Size (width x length)	Pack contents	Shipping unit	PZN
9802	150 x 18 mm (approx.)	100 pcs	20 packs	00822334

### Single sterile, 1 pc. per pack

9665	150 x 18 mm (approx.)	100 x 1 pc.	10 packs	01333821
------	-----------------------	-------------	----------	----------

## ABE® injection plaster

Our fabric injection plaster

### Features

- ready-to-use injection plaster
- with wound pad
- polyester net reduces risk of dressing pad becoming stuck to the wound
- with polyacrylate adhesive
- woven fabric structure



### Usage

For minor injuries such as scrapes/grazes and shallow cuts, and for wound care after taking blood.

### Non-sterile

REF	Size (width x length)	Pack contents	Shipping unit	PZN
9525	2 x 4 cm (approx.)	500 pcs	20 packs	01664860

## ABE® light

Our elastic wound bandage with wound dressing pad

### Features

- white, slightly elastic wound bandage
- pad placed in the middle
- polyester net reduces risk of dressing pad becoming stuck to the wound
- with polyacrylate adhesive



### Usage

For first aid/follow-up treatment of minor injuries, also on sensitive parts of the body. Good fit for wounds to joints.

### Non-sterile

REF	Size (width x length)	Pack contents	Shipping unit	PZN
9547	4 cm x 5 m (approx.)	1 pc.	14 packs	06790233
9548	6 cm x 5 m (approx.)	1 pc.	10 packs	06790256
9549	8 cm x 5 m (approx.)	1 pc.	8 packs	06790262

## ABE® Absorbent adhesive dressing

Our absorbent and hypoallergenic instant wound dressing

### Features

- self-adhesive wound dressing made from white non-woven material
- with polyacrylate adhesive
- with central wound pad
- reduces risk of sticking to the wound



### Usage

For post-operative wound care or in the ward as a sterile instant wound bandage for absorption of blood and exudate, and to protect wounds from (cross-)infection.

### Single sterile, shrink-wrapped pack of 1 pc.

REF	Size (width x length)	Pack contents	Shipping unit	PZN
4750	5 x 7 cm (approx.)	100 x 1 pc.	12 packs	01452997
4752	6 x 10 cm (approx.)	50 x 1 pc.	6 packs	01126022
4754	8 x 10 cm (approx.)	50 x 1 pc.	6 packs	01446011
4760	15 x 10 cm (approx.)	50 x 1 pc.	6 packs	01445980
4762	20 x 10 cm (approx.)	50 x 1 pc.	6 packs	01445968
4764	25 x 10 cm (approx.)	25 x 1 pc.	6 packs	01445951
4768	35 x 10 cm (approx.)	25 x 1 pc.	6 packs	01446057



## ABE<sup>®</sup> derm

Transparent film without absorbant pad

### Features

- elastic, semi-permeable polyurethane film
- coated with polyacrylate adhesive
- wound environment visible through the transparent film



### Usage

Suitable for postoperative wound care and for monitoring the wound environment, wound care and catheter fixation. Can also be placed directly on primary wounds for wound monitoring.

Single sterile, shrink-wrapped pack of 1 pc.

REF	Size (width x length)	Pack contents	Shipping unit	PZN
9571	6 x 7 cm (approx.)	100 x 1 pc.	40 packs	06790204
9572	10 x 12 cm (approx.)	50 x 1 pc.	12 packs	06790210
9573	15 x 20 cm (approx.)	50 x 1 pc.	15 packs	06790227

## ABE<sup>®</sup> derm plus

Transparent film with absorbant pad

### Features

- elastic, semi-permeable polyurethane film
- coated with polyacrylate adhesive
- integrated wound pad
- wound environment visible through the transparent film



### Usage

Suitable for postoperative wound care and for monitoring the wound environment, wound care and catheter fixation. Can also be placed directly on primary wounds for wound monitoring.

Single sterile, shrink-wrapped pack of 1 pc.

REF	Size (width x length)	Pack contents	Shipping unit	PZN
9580	5 x 7 cm (approx.)	50 x 1 pc.	24 packs	04589917
9581	6 x 10 cm (approx.)	50 x 1 pc.	12 packs	04589923
9584	15 x 10 cm (approx.)	50 x 1 pc.	12 packs	04656720
9585	20 x 10 cm (approx.)	50 x 1 pc.	24 packs	04656737

## ABE<sup>®</sup> cannula retention dressing

Our cannula fixing plaster

### Features

- self-adhesive fixing dressing
- with central slit
- less risk of the integrated wound pad sticking to the wound
- coated with polyacrylate adhesive



### Usage

For sterile wound dressing of the puncture site and the fixing of venous access ports.

Single sterile, shrink-wrapped pack of 1 pc.

REF	Size (width x length)	Pack contents	Shipping unit	PZN
4457	8 x 6 cm (approx.)	50 x 1 pc.	10 packs	12589854

## ABE<sup>®</sup> non-woven material

Our self-adhesive non-woven fixing material

### Features

- self-adhesive non-woven fixing material, white
- stretchable widthways (min. 100%)
- less risk of the integrated wound pad sticking to the wound
- coated with polyacrylate adhesive



### Usage

For covering and fixing large bandaged areas in medical applications. For securing cannulas, drainage tubes, probes, catheters, measuring instruments, etc. For securing bandages to joints, and other round/tapered parts of the body.

Non-sterile

REF	Size (width x length)	Pack contents	Shipping unit	PZN
4470	5 cm x 10 m (approx.)	1 pc.	50 packs	06168663
4471	10 cm x 10 m (approx.)	1 pc.	24 packs	06168686
4472	15 cm x 10 m (approx.)	1 pc.	24 packs	06168692
4473	20 cm x 10 m (approx.)	1 pc.	52 packs	06168700

# ABE<sup>®</sup> silk

Our synthetic silk plaster on a roll

### Features

- fixing plaster made from white cellulose acetate
- coated with polyacrylate adhesive
- zigzag edge, easily torn off widthways or lengthways



### Usage

For fixing wound coverings, bandages, drainage tubes, catheters and cannulas in place.

### Non-sterile, with protective ring

REF	Size (width x length)	Pack contents	Shipping unit	PZN
9531	2.5 cm x 5 m (approx.)	24 pcs	10 packs	11732887

### Non-sterile, without protective ring

9535	1.25 cm x 5 m (approx.)	36 pcs	10 packs	11734544
9536	2.5 cm x 5 m (approx.)	24 pcs	10 packs	11745536
9540	1.25 cm x 9.14 m (approx.)	36 pcs	10 packs	03916834
9541	2.5 cm x 9.14 m (approx.)	24 pcs	10 packs	03916840

# BeeSana<sup>®</sup> SafeCut

Our self-adhering finger bandage

### Features

- self-adhering, elastic fixation bandage
- fine latex coating prevents slippage
- stretches by approx. 100%



### Usage

For personal first aid for finger and joint injuries, and for torn fingernails. For fixing wound coverings, and for use as a light support and compression bandage. For immobilising injured parts of the body and limbs.

### Non-sterile

REF	Size (width x length, stretched)	Colour	Pack contents	Shipping unit	PZN
1844	2.5 cm x 4.5 m (approx.)	Blue	1 pc.	252 packs	10544882

### Fold the support pad (fixation possible with or without pad)



1. Fold the end of the bandage over itself a few times until the support pad created is adequate.
2. Fix the support pad onto the intended location with the rest of the bandage until the bandage is stable enough.
3. Apply the bandage firmly and then tear off.

# ABE® gauze bandage

Our traditional fixing bandage

## Features

- rigid gauze bandage (DIN 61631 standard)
- smooth structure
- woven edges
- 65% viscose/35% cotton



## Usage

Broad field of use for medical applications, e.g. as an emergency bandage for first aid or as a carrier material for zinc paste bandages.

## Non-sterile

REF	Size (width x length)	Pack contents	Shipping unit	PZN
1400	4 cm x 4 m (approx.)	20 pcs	16 packs	01331029
1401	6 cm x 4 m (approx.)	20 pcs	16 packs	01331035
1402	8 cm x 4 m (approx.)	20 pcs	12 packs	01331041
1403	10 cm x 4 m (approx.)	20 pcs	8 packs	01331058

# ABE® Last

Our elastic conforming bandage with a smooth structure

## Features

- elastic conforming bandage (DIN 61634 standard)
- fine, smooth structure
- fully stretched length 4 m (approx.)
- 54% polyester/46% viscose



## Usage

For fixing wound coverings, for fixing support splints, and for use as a light support and compression bandage. Suitable for use on joints, tapered and rounded body parts. For immobilising injured parts of the body and limbs.

## Non-sterile, pack of 20 pcs.

REF	Size (width x length, stretched)	Pack contents	Shipping unit	PZN
1409	4 cm x 4 m (approx.)	20 pcs	28 packs	04853248
1410	6 cm x 4 m (approx.)	20 pcs	24 packs	01364282
1411	8 cm x 4 m (approx.)	20 pcs	20 packs	04853260
1412	10 cm x 4 m (approx.)	20 pcs	15 packs	04853277
1413	12 cm x 4 m (approx.)	20 pcs	6 packs	04853283



## ABE® Last crepe

Our elastic conforming bandage with a crepe structure

### Features

- elastic conforming bandage (DIN 61634 standard)
- crepe structure
- adhesion of individual layers
- fully-stretched length 4 m (approx.)
- 70% viscose/30% polyamide



### Usage

For fixing wound coverings, for fixing support splints, and for use as a light support and compression bandage. Suitable for use on joints, tapered and rounded body parts.

### Non-sterile

REF	Size (width x length, stretched)	Pack contents	Shipping unit	PZN
1749	4 cm x 4 m (approx.)	40 pcs	12 packs	14439863
1750	6 cm x 4 m (approx.)	20 pcs	16 packs	07198324
1751	8 cm x 4 m (approx.)	20 pcs	12 packs	07198330
1752	10 cm x 4 m (approx.)	20 pcs	10 packs	07198347
1753	12 cm x 4 m	20 pcs	8 packs	07198353

### Double sterile, pack of 1 pc. with 2 counting cards

4741	6 cm x 4 m (approx.)	30 x 1 pc.	4 packs	13655731
4742	8 cm x 4 m (approx.)	30 x 1 pc.	4 packs	13655642
4588	12 cm x 4 m (approx.)	30 x 1 pc.	4 packs	13655725

## ABE® fix color

Our elastic and cohesive coloured conforming bandage, smooth

### Features

- cohesive and elastic conforming bandage
- smooth structure
- distensibility (DIN 61632 standard) min. 40%
- 61% polyamide/39% cotton



### Usage

For fixing wound coverings and bandages on body parts that are moved frequently. Such as for isolating wrist or ankle joints, for example.

### Non-sterile

REF	Size (width x length, stretched)	Colour	Pack contents	Shipping unit	PZN
19461	6 cm x 20 m (approx.)	Blue	1 pc.	15 packs	08520615
19463	6 cm x 20 m (approx.)	Green	1 pc.	36 packs	08520638
19465	6 cm x 20 m (approx.)	Red	1 pc.	15 packs	08520650

# ABE® Last cohesive

Our elastic and cohesive conforming bandage

## Features

- cohesive, elastic conforming bandage
- crepe structure
- no additional fixation required
- microdot latex impregnation
- large-pored
- stretches by approx. 60%
- 70% polyester/30% polyamide



## Usage

For fixing wound coverings, for fixing support splints, and for use as a light support and compression bandage. Suitable for use on joints, tapered and rounded body parts. For immobilising injured parts of the body and limbs.

### Non-sterile

REF	Size (width x length, stretched)	Pack contents	Shipping unit	PZN
1289	4 cm x 4 m (approx.)	1 pc.	30 packs	04019746
1391	6 cm x 4 m (approx.)	1 pc.	36 packs	07198436
1392	8 cm x 4 m (approx.)	1 pc.	36 packs	07198442
1393	10 cm x 4 m (approx.)	1 pc.	36 packs	07198459
1394	12 cm x 4 m (approx.)	1 pc.	36 packs	07198465
1288	4 cm x 20 m (approx.)	1 pc.	22 packs	02435498
1383	6 cm x 20 m (approx.)	1 pc.	16 packs	04853308
1384	8 cm x 20 m (approx.)	1 pc.	12 packs	04853314
1385	10 cm x 20 m (approx.)	1 pc.	10 packs	04853320
1386	12 cm x 20 m (approx.)	1 pc.	10 packs	04853337

### Double sterile, pack of 1 pc. with counting card

4793	12 cm x 4 m	15 x 1 pc.	4 packs	–
------	-------------	------------	---------	---

# ABE® Ideal 100

Our elastic short-stretch fabric compression bandage

## Features

- elastic fabric bandage (DIN 61632 standard)
- natural white
- 100% cotton
- fully stretched length 5 m (approx.)
- selvedge
- washable



## Usage

For relief bandages and follow-up treatment of dislocations or fractures.

### Non-sterile

REF	Size (width x length, stretched)	Pack contents	Shipping unit	PZN
1589	4 cm x 5 m (approx.)	10 pcs	10 packs	01364299
1423	6 cm x 5 m (approx.)	10 pcs	10 packs	01364477
1424	8 cm x 5 m (approx.)	10 pcs	10 packs	01364448
1425	10 cm x 5 m (approx.)	10 pcs	5 packs	04853484
1426	12 cm x 5 m (approx.)	10 pcs	5 packs	04853490
1576	15 cm x 5 m (approx.)	10 pcs	5 packs	01364460
1577	20 cm x 5 m (approx.)	10 pcs	5 packs	01364454

### Double sterile, pack of 1 pc. with counting card

4771	6 cm x 5 m (approx.)	25 x 1 pc.	4 packs	–
4772	8 cm x 5 m (approx.)	25 x 1 pc.	4 packs	–
4773	10 cm x 5 m (approx.)	20 x 1 pc.	4 packs	–
4774	12 cm x 5 m (approx.)	15 x 1 pc.	4 packs	–

# ABE® universal

Our permanently-elastic, medium-stretch compression bandage

## Features

- permanently elastic bandage for medium compression
- open, air-permeable structure
- fully stretched length: approx. 5 m
- latex-free
- 67% cotton/30% polyester/3% elastane



## Usage

For bandages where medium compression is required, and especially for compressive wound bandages. For fixing wound coverings. For immobilising injured parts of the body and limbs, and for supporting and relieving pressure on limbs.

### Non-sterile

REF	Size (width x length, stretched)	Pack contents	Shipping unit	PZN
1414	6 cm x 5 m (approx.)	10 pcs	16 packs	00761561
1415	8 cm x 5 m (approx.)	10 pcs	12 packs	00761590
1416	10 cm x 5 m (approx.)	10 pcs	10 packs	00761526
1432	12 cm x 5 m (approx.)	10 pcs	8 packs	00761532
1433	15 cm x 5 m (approx.)	10 pcs	6 packs	00761549
1434	20 cm x 5 m (approx.)	10 pcs	5 packs	00761555

### Double sterile, pack of 1 pc. with counting card

4722	6 cm x 5 m (approx.)	25 x 1 pc.	4 packs	13655636
4723	8 cm x 5 m (approx.)	25 x 1 pc.	4 packs	13655694
4572	10 cm x 5 m (approx.)	20 x 1 pc.	4 packs	13655748
4560	12 cm x 5 m (approx.)	15 x 1 pc.	4 packs	13655719

# ABE® lastic

Our permanently elastic, light short-stretch compression bandage

## Features

- permanently elastic bandage for light compression
- open, air-permeable structure
- fully stretched length: approx. 5 m
- latex-free
- 60% cotton/38% polyester/2% elastane



## Usage

For bandages where light compression is required, and especially for compressive wound support and pressure relief bandages. For immobilising injured parts of the body and limbs. For fixing wound coverings.

### Non-sterile

REF	Size (width x length, stretched)	Pack contents	Shipping unit	PZN
1780	4 cm x 5 m (approx.)	10 pcs	10 packs	02423265
1781	6 cm x 5 m (approx.)	10 pcs	10 packs	02423259
1782	8 cm x 5 m (approx.)	10 pcs	10 packs	02423242
1783	10 cm x 5 m (approx.)	10 pcs	10 packs	02423236
1784	12 cm x 5 m (approx.)	10 pcs	5 packs	02423213
1785	15 cm x 5 m (approx.)	10 pcs	5 packs	02423207
1786	20 cm x 5 m (approx.)	10 pcs	5 packs	02423182

### Double sterile, pack of 1 pc. with counting card

4779	6 cm x 5 m (approx.)	30 x 1 pc.	4 packs	13655671
4780	8 cm x 5 m (approx.)	25 x 1 pc.	4 packs	13655659
4781	10 cm x 5 m (approx.)	20 x 1 pc.	4 packs	13655702
4782	12 cm x 5 m (approx.)	15 x 1 pc.	4 packs	13655665



# ABE<sup>®</sup> padded bandage

Our synthetic wadding for padding

### Features

- synthetic wadding bandage made from polyester fibre
- latex-free
- can be sterilised

### Usage

As a padding material under plaster casts, and support/compression bandages.



### Non-sterile

REF	Size (width x length)	Pack contents	Shipping unit	PZN
1701	6 cm x 3 m (approx.)	1 pc.	108 packs	06627958
1806	8 cm x 3 m (approx.)	1 pc.	72 packs	04351788
1702	10 cm x 3 m (approx.)	1 pc.	72 packs	06627964
1703	15 cm x 3 m (approx.)	1 pc.	36 packs	06627970

### Single sterile, pack of 1 pc.

1705	10 cm x 3 m (approx.)	15 x 1 pc.	4 packs	13715901
------	-----------------------	------------	---------	----------

# ABE<sup>®</sup> extra

Our short-stretch support and compression bandages

### Features

- powerful, fabric elastic compression bandage
- short-stretch
- open, air-permeable structure
- fully stretched length 5 m (approx.)
- 100% cotton
- with bandage clips



### Usage

For strong compression. Suitable for opposing bandaging techniques, according to Pütter. For stabilising, relief and fixation during the treatment of injuries to the musculoskeletal system.

### Non-sterile

REF	Size (width x length, stretched)	Pack contents	Shipping unit	PZN
1421	8 cm x 5 m (approx.)	1 pc.	72 packs	01364402
1422	10 cm x 5 m (approx.)	1 pc.	36 packs	01364394
1735	6 cm x 5 m (approx.)	10 pcs	10 packs	08490759
1736	8 cm x 5 m (approx.)	10 pcs	10 packs	08490765
1737	10 cm x 5 m (approx.)	10 pcs	5 packs	08490771
1738	12 cm x 5 m (approx.)	10 pcs	5 packs	08490788

# ABE<sup>®</sup> extra cohesive

Our cohesive, elastic short-stretch compression bandage

### Features

- strong, permanently elastic compression bandage
- short-stretch
- self-adhering
- selvedge
- can be sterilised with gas and steam at 121 °C, x-ray and gamma
- fully-stretched length 5 m (approx.)
- 88% cotton/11% polyamide/1% polyurethane



### Usage

For stabilising, relief and fixation during the treatment of injuries to the musculoskeletal system.

### Non-sterile

REF	Size (width x length, stretched)	Pack contents	Shipping unit	PZN
1777	8 cm x 5 m (approx.)	1 pc.	64 packs	01340933
1778	10 cm x 5 m (approx.)	1 pc.	32 packs	01340956
1779	12 cm x 5 m (approx.)	1 pc.	32 packs	01340979

# ABE<sup>®</sup> Ven

Our elastic, long-stretch compression bandage

### Features

- permanently elastic
- with latex-free end fixation
- 93% cotton/7% polyurethane
- fully stretched length 7 m (approx.)
- washable
- with bandage clips



### Usage

For stabilising, relief and fixation during the treatment of injuries to the musculoskeletal system.

### Non-sterile

REF	Size (width x length, stretched)	Pack contents	Shipping unit	PZN
1563	8 cm x 7 m (approx.)	1 pc.	64 packs	04854147
1564	10 cm x 7 m (approx.)	1 pc.	32 packs	04854153
1565	12 cm x 7 m (approx.)	1 pc.	32 packs	04854176

# ABE<sup>®</sup> tubular bandage

Our tubular support bandage for universal use

## Features

- tube bandage made from 100% cotton
- stretches lengthways and widthways
- easy to use
- can be cut anywhere to fit
- seamless to avoid creating pressure points



## Usage

As a pressure-free bandage and an underlay or overlay to protect the skin for zinc paste/plaster bandages. For immobilising and traction bandages. As an overlay for splints.

### Non-sterile

REF	Size (width x length)	Usage	Pack contents	Shipping unit	PZN
5181	Size 1 = 1.5 cm x 20 m (approx.)	Finger/toe	1 pc.	24 packs	01331354
5182	Size 2 = 2.5 cm x 20 m (approx.)	Several fingers/toes	1 pc.	16 packs	01331360
5183	Size 3 = 4 cm x 20 m (approx.)	Hand/forearm/child's arm	1 pc.	12 packs	01331377
5184	Size 4 = 6 cm x 20 m (approx.)	Upper arm/lower leg	1 pc.	28 packs	01331383
5185	Size 5 = 8 cm x 20 m (approx.)	Thigh	1 pc.	22 packs	01331420
5186	Size 6 = 10 cm x 20 m (approx.)	Head/child's torso	1 pc.	18 packs	01331437
5187	Size 7 = 12 cm x 20 m (approx.)	Torso	1 pc.	16 packs	01331443

# ABE<sup>®</sup> tubular net bandage

Our elastic tubular net bandage for universal use

## Features

- elastic tube bandage
- simple to use
- can be cut anywhere to fit
- binding prevents bandage unravelling further
- seamless to avoid creating pressure points
- follows anatomical contours to ensure a good fit



## Usage

For fixing bandages and padding materials. Scissors can be used to cut a hole in the mesh material, so as to permit patients to exclude parts of the body (such as the thumb with a hand bandage) from being wrapped in the bandage. Suitable for use on joints, tapered and rounded body parts.

### Non-sterile

REF	Size (width x length)	Pack contents	Shipping unit	PZN
5164	Size 3 = forearm (approx., 2.2 cm x 25 m)	1 pc.	10 packs	03065804
5166	Size 4 = upper arm (approx., 2.5 cm x 25 m)	1 pc.	8 packs	03065810
5170	Size 5 = leg (approx., 4 cm x 25 m)	1 pc.	5 packs	03065827
5172	Size 6 = head (approx., 5.5 cm x 25 m)	1 pc.	4 packs	03065833
5176	Size 7 = torso (approx., 6.2 cm x 25 m)	1 pc.	4 packs	03065856
5180	Size 8 = torso XXL (approx., 8 cm x 25 m)	1 pc.	3 packs	03065862

## ABE® Gilchrist's bandage

Our padded shoulder immobiliser sling

### Features

- universal shoulder bandage



### Usage

For injuries to the shoulder and upper arm area, for immobilising or fixing the shoulder joint.

### Non-sterile

REF	Size	Pack contents	Shipping unit	PZN
1914	Universal	1 pc.	12 packs	01363437

## ABE® body bandage

Our post-operative support bandage

### Features

- support and relief bandage
- wide application padding
- broad, adjustable hook-and-pile fastener
- compression can be adjusted to therapeutic requirements



### Usage

For post-operative support and treatment of fractured ribs.

### Non-sterile

REF	Size	(circumference)	Pack contents	Shipping unit	PZN
1890	Size 1	65–75 cm (approx.)	1 pc.	15 packs	01363727
1891	Size 2	75–85 cm (approx.)	1 pc.	15 packs	01363710
1892	Size 3	85–95 cm (approx.)	1 pc.	15 packs	01363704
1893	Size 4	95–105 cm (approx.)	1 pc.	15 packs	01363696
1894	Size 5	105–115 cm (approx.)	1 pc.	15 packs	01363609
1895	Size 6	115–125 cm (approx.)	1 pc.	15 packs	01363555
1904	Size 7	125–135 cm (approx.)	1 pc.	15 packs	02435357



## ABE® heart catheter pressure bandage

Our cardiac catheter pressure bandage

### Features

- adjustable hook-and-pile fastener
- compression can be adjusted to therapeutic requirements
- colour-coded



### Usage

For compression after femoral artery/vein punctures for invasive diagnostic work/treatment.

### Non-sterile

REF	Size (universal)	Pack contents	Shipping unit	PZN
2907	Left (colour code: red)	1 pc.	12 packs	04351742
2908	Right (colour code: blue)	1 pc.	12 packs	04351759

## ABE® limb holder

Our limb support cuff

### Features

- cuff from special, soft material
- hypoallergenic, thanks to cotton contact area
- hook-and-pile fastener for simple fixation
- no additional fixing strap required, since already integrated into the cuff itself



### Usage

For short-term immobilisation.

### Non-sterile

REF	Size	Pack contents	Shipping unit	PZN
1879	For adults	10 pcs	12 packs	12553046

## ABE® pressure bandage roll

Special bandage

### Features

- pressure bandage roll made from polyurethane, in knitted fabric tube
- diameter 7 cm (approx.)
- latex-, PVC- and DEHP-free
- disposable



### Usage

As a pressure padding, for example after heart catheterisation.

### Non-sterile

REF	Size	Pack contents	Shipping unit	PZN
4037	7 cm Ø x 10 cm (approx.)	12 pcs	9 packs	06620850

## Scrotum roll

### Features

- scrotum roll from foamed material
- with a knitted fabric tube covering

### Usage

For elevating the scrotum after operations or with oedemas.



### Non-sterile

REF	Size	Pack contents	Shipping unit	PZN
4035	Ø 8 cm x 10 cm (approx.)	6 pcs	8 packs	01333755

# BeeSana® X-Ray lap sponge

## Our abdominal swabs with x-ray chip

### Features

- pre-washed abdominal swabs from medical gauze (EN 14079 standard)
- with X-ray contrast chip
- 20-thread
- absorbent
- specified formats: Dimensions PRIOR to washing in production



### Usage

For absorbing and staunching blood and bodily fluids during surgical procedures.

The abdominal swabs come in various sizes, layer counts and colours to suit a wide range of purposes.

### Non-sterile, pack of 5 pcs, bundled with wrapper

		REF	Size (width x length)	Loop	Pack contents	Shipping unit	PZN
4-ply	White	1613	20 x 30 cm (approx.)	Yes	4 x 5 pcs	36 packs	13506907
		1547	40 x 40 cm (approx.)	Yes	4 x 5 pcs	18 packs	13506913
		1616	45 x 45 cm (approx.)	Yes	4 x 5 pcs	12 packs	13506936
		1648	45 x 45 cm (approx.)	No	4 x 5 pcs	12 packs	13506942
		1670	50 x 60 cm (approx.)	Yes	4 x 5 pcs	8 packs	13506959
	Unwashed	1602	50 x 60 cm (approx.)	No	4 x 5 pcs	10 packs	13506965
		1465	8 x 90 cm (approx.)	Yes	4 x 5 pcs	30 packs	13506994
	Green	1614	20 x 30 cm (approx.)	Yes	4 x 5 pcs	36 packs	13507002
		1617	45 x 45 cm (approx.)	Yes	4 x 5 pcs	12 packs	13507019
		1607	45 x 45 cm (approx.)	No	4 x 5 pcs	12 packs	13507025
1696		8 x 90 cm (approx.)	Yes	4 x 5 pcs	30 packs	13507031	
Blue	4593	45 x 45 cm (approx.)	Yes	4 x 5 pcs	12 packs	13507048	
6-ply	White	1612	20 x 30 cm (approx.)	Yes	4 x 5 pcs	30 packs	13507054
		1502	40 x 40 cm (approx.)	Yes	4 x 5 pcs	15 packs	13507083
		1618	45 x 45 cm (approx.)	Yes	4 x 5 pcs	10 packs	13507108
	Unwashed	1520	45 x 45 cm (approx.)	Yes	4 x 5 pcs	10 packs	13507120
		1659	45 x 45 cm (approx.)	No	4 x 5 pcs	10 packs	13507137
		1671	50 x 60 cm (approx.)	Yes	4 x 5 pcs	6 packs	13507143
		1449	8 x 90 cm (approx.)	Yes	4 x 5 pcs	27 packs	13507166
	Green	1615	20 x 30 cm (approx.)	Yes	4 x 5 pcs	30 packs	13507172
		1496	20 x 30 cm (approx.)	Yes	4 x 5 pcs	10 packs	13507189
		1608	45 x 45 cm (approx.)	Yes	4 x 5 pcs	10 packs	13507203
	24-thread	1651	45 x 45 cm (approx.)	No	4 x 5 pcs	10 packs	13507226
		1583	50 x 60 cm (approx.)	Yes	2 x 5 pcs	10 packs	13507232
		1603	50 x 60 cm (approx.)	No	4 x 5 pcs	6 packs	13507249
	Blue	4058	8 x 90 cm (approx.)	Yes	4 x 5 pcs	27 packs	13507255
4592		45 x 45 cm (approx.)	Yes	4 x 5 pcs	10 packs	13507261	
4594		45 x 45 cm (approx.)	No	4 x 5 pcs	10 packs	13507278	
8-ply	White	1506	45 x 45 cm (approx.)	Yes	2 x 5 pcs	10 packs	13507284

### Double sterile, with self-adhesive counting cards

				REF	Size (width x length)	Loop	Pack contents	Shipping unit	PZN
4-ply	Pack of 2 pcs	White	Double-peel	1664	45 x 45 cm (approx.)	Yes	15 x 2 pcs	4 packs	13166417
			Double-peel	4515	50 x 60 cm (approx.)	Yes	24 x 2 pcs	2 packs	13166423
			Double-peel	1682	8 x 90 cm (approx.)	Yes	30 x 2 pcs	4 packs	13166446
	Green	Double-peel	1665	45 x 45 cm (approx.)	Yes	15 x 2 pcs	4 packs	13166452	
		Double-peel	4504	50 x 60 cm (approx.)	Yes	24 x 2 pcs	2 packs	13166469	
		Double-peel	1684	8 x 90 cm (approx.)	Yes	30 x 2 pcs	4 packs	13166475	
	Pack of 5 pcs	White	Double-peel	1571	20 x 30 cm (approx.)	Yes	15 x 5 pcs	4 packs	13614399
			Double-peel	1593	45 x 45 cm (approx.)	Yes	15 x 5 pcs	2 packs	13166498
	Unwashed, 24-thread	Green	Double-peel	1597	30 x 30 cm (approx.)	No	18 x 5 pcs	2 packs	13166297
Double-peel			1572	20 x 30 cm (approx.)	Yes	15 x 5 pcs	4 packs	13166280	
Double-peel			1594	45 x 45 cm (approx.)	Yes	15 x 5 pcs	2 packs	13166334	
6-ply	Pack of 2 pcs	White	Double-peel	1666	45 x 45 cm (approx.)	Yes	15 x 2 pcs	4 packs	13166357
			Double-peel	1677	45 x 45 cm (approx.)	No	15 x 2 pcs	4 packs	13166363
			Double-peel	1685	8 x 90 cm (approx.)	Yes	30 x 2 pcs	4 packs	14218665
	Pack of 2 pcs	Green	Double-peel	1667	45 x 45 cm (approx.)	Yes	15 x 2 pcs	4 packs	13166386
			Double-peel	4501	50 x 60 cm (approx.)	Yes	20 x 2 pcs	2 packs	13166185
			Double-peel	1680	8 x 90 cm (approx.)	Yes	30 x 2 pcs	4 packs	13166191
Pack of 5 pcs	White	Double-peel	1595	45 x 45 cm (approx.)	Yes	12 x 5 pcs	2 packs	13166222	
		Double-peel	1591	45 x 45 cm (approx.)	No	12 x 5 pcs	2 packs	13166239	
		Double-peel	1669	8 x 90 cm (approx.)	Yes	24 x 5 pcs	2 packs	13166245	
Pack of 5 pcs	Green	Double-peel	1574	20 x 30 cm (approx.)	Yes	9 x 5 pcs	4 packs	13166251	
		Double-peel	1587	40 x 40 cm (approx.)	Yes	12 x 5 pcs	2 packs	13166268	
		Double-peel	1596	45 x 45 cm (approx.)	Yes	12 x 5 pcs	2 packs	13166400	

# BeeSana® X-Ray gauze swab

Our classic gauze swab with X-ray contrast thread

## Features

- gauze swabs (EN 14079 standard)
- woven-in X-ray contrast threads (BaSO<sub>4</sub>)
- selvedge cut edges
- folded, 17-thread
- absorbent
- non-sterile in sterilisable paper packaging
- double sterile packaged with self-adhesive counting cards for secure documentation during operations

## Usage

For (invasive) surgical interventions in the operating theatre.  
For staunching bleeding and wound care.



## Non-sterile

	REF	Size (width x length)	Pack contents	Shipping unit	PZN
<b>8-ply</b>	1013	10 x 10 cm (approx.)	100 pcs	10 packs	12348599
	1037	10 x 20 cm (approx.)	100 pcs	5 packs	12348719
<b>12-ply</b>	1003	5 x 5 cm (approx.)	100 pcs	10 packs	12348582
	1009	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	06613703
	1015	10 x 10 cm (approx.)	100 pcs	10 packs	06613726
	1039	10 x 13 cm (approx.)	100 pcs	10 packs	12348731
	1024	10 x 20 cm (approx.)	100 pcs	5 packs	06613732
	10 pcs, pre-counted For set packers	2715	10 x 10 cm (approx.)	200 x 10 pcs	1 box
	2724	10 x 20 cm (approx.)	100 x 10 pcs	1 box	–
<b>16-ply</b>	1005	5 x 5 cm (approx.)	100 pcs	10 packs	06613749
	1011	7.5 x 7.5 cm (approx.)	100 pcs	10 packs	06613755
	1017	10 x 10 cm (approx.)	100 pcs	10 packs	06613761
	1054	10 x 13 cm (approx.)	100 pcs	10 packs	12348665
	1034	10 x 20 cm (approx.)	100 pcs	5 packs	12348702
	10 pcs, pre-counted For set packers	2717	10 x 10 cm (approx.)	200 x 10 pcs	1 box
	2727	20 x 20 cm (approx.)	60 x 10 pcs	1 box	–
<b>24-ply</b>	1090	10 x 20 cm (approx.)	50 pcs	5 packs	12348642
<b>32-ply</b>	1019	10 x 10 cm (approx.)	50 pcs	10 packs	12348607
	20-thread	1021	10 x 10 cm (approx.)	50 pcs	10 packs

## Double sterile, with self-adhesive counting cards

		REF	Size (width x length)	Pack contents	Shipping unit	PZN
<b>12-ply</b>	Pack of 5 pcs	1149	10 x 10 cm (approx.)	25 x 5 pcs	8 packs	12348754
		1150	10 x 20 cm (approx.)	30 x 5 pcs	4 packs	12348760
	Pack of 10 pcs	4392	5 x 5 cm (approx.)	30 x 10 pcs	8 packs	12348932
		1174	7.5 x 7.5 cm (approx.)	40 x 10 pcs	4 packs	12348843
		1049	10 x 10 cm (approx.)	30 x 10 pcs	4 packs	11703472
		1168	10 x 20 cm (approx.)	15 x 10 pcs	4 packs	12348777
		1176	7.5 x 7.5 cm (approx.)	20 x 20 pcs	4 packs	12348866
	Pack of 20 pcs	4314	10 x 10 cm (approx.)	15 x 20 pcs	4 packs	12348820
		1041	10 x 20 cm (approx.)	20 x 20 pcs	2 packs	12348659
		1142	10 x 20 cm (approx.)	15 x 10 pcs	4 packs	12348808
<b>16-ply</b>	Pack of 10 pcs	4368	5 x 5 cm (approx.)	30 x 10 pcs	8 packs	12348990
		1071	7.5 x 7.5 cm (approx.)	30 x 10 pcs	4 packs	12348613
		1072	10 x 10 cm (approx.)	25 x 10 pcs	4 packs	12348636
		1142	10 x 20 cm (approx.)	15 x 10 pcs	4 packs	12348808
	Pack of 20 pcs	4351	10 x 10 cm (approx.)	12 x 20 pcs	4 packs	12348984
<b>32-ply</b>	Pack of 10 pcs	4318	10 x 10 cm (approx.)	15 x 10 pcs	4 packs	12348837



# BeeSana® X-Ray gauze balls

## Features

- gauze ball, medical gauze (EN 14079 standard)
- woven-in X-ray contrast thread (BaSO<sub>4</sub>)
- innovative wrapping system from just one piece of pre-cut gauze
- inner-facing cut edges
- 20-thread
- absorbent
- for documentation during the operation with counting cards (for sterile gauze swabs with X-ray contrast only)



### Non-sterile, in polybag

	REF	Size (width x length)	Shape	Pack contents	Shipping unit	PZN
Plum-sized	1154	20 x 20 cm (approx.)	Round	250 pcs	16 packs	00822297
Egg-sized	1155	25 x 26 cm (approx.)	Round	250 pcs	10 packs	00822311
Egg-sized	1110	25 x 26 cm (approx.)	Thimble	250 pcs	10 packs	00822305
X-sized	1349	30 x 30 cm (approx.)	Round	100 pcs	20 packs	13508243
XX-sized	1114	30 x 40 cm (approx.)	Round	100 pcs	15 packs	01363868
XXX-sized	1115	40 x 40 cm (approx.)	Round	100 pcs	10 packs	13508266
Fist-sized	1117	40 x 50 cm (approx.)	Round	100 pcs	10 packs	00822328
Abdominal swab	1119	60 x 80 cm (approx.)	Round	50 pcs	8 packs	13508295



### Double sterile, with self-adhesive counting cards

Pack of 5 pcs	Egg-sized	4324	25 x 26 cm (approx.)	Round	25 x 5 pcs	8 packs	13165145
	XX-sized	4326	30 x 40 cm (approx.)	Round	30 x 5 pcs	4 packs	13165151
	Fist-sized	4325	40 x 50 cm (approx.)	Round	20 x 5 pcs	4 packs	13165168
Pack of 10 pcs	Walnut-sized	1180	15 x 16 cm (approx.)	Round	30 x 10 pcs	8 packs	13165174
	Plum-sized	1181	20 x 20 cm (approx.)	Round	18 x 10 pcs	8 packs	13165180
	Egg-sized	1075	25 x 26 cm (approx.)	Round	30 x 10 pcs	4 packs	11703503
	Egg-sized	1186	25 x 26 cm (approx.)	Thimble	30 x 10 pcs	4 packs	13165197
	X-sized	1182	30 x 30 cm (approx.)	Round	20 x 10 pcs	4 packs	13165211
	XX-sized	4319	30 x 40 cm (approx.)	Round	18 x 10 pcs	4 packs	13165292
Pack of 20 pcs	Fist-sized	4333	40 x 50 cm (approx.)	Round	12 x 10 pcs	4 packs	13165228
	Plum-sized	4329	20 x 20 cm (approx.)	Round	25 x 20 pcs	4 packs	13165234
	Egg-sized	4337	25 x 26 cm (approx.)	Round	18 x 20 pcs	4 packs	13165240
	Egg-sized	4012	25 x 26 cm (approx.)	Thimble	18 x 20 pcs	4 packs	13165257
	XX-sized	4320	30 x 40 cm (approx.)	Round	20 x 20 pcs	2 packs	13165263

# BeeSana® X-Ray preparation swabs

## Features

- gauze swab, 100% cotton (EN 14079 standard)
- woven-in X-ray contrast thread (BaSO<sub>4</sub>)
- 24-thread
- innovative wrapping system from just one piece of pre-cut gauze
- highly absorbent



## Usage

For free preparation of tissue layers and anatomical structures.

### Non-sterile

	REF	Size (width x length)	Pack contents	Shipping unit	PZN
Packed in polybag	1169	4 x 4 cm (approx.)	100 pcs	50 packs	13583440
	1170	6 x 6 cm (approx.)	100 pcs	40 packs	01363839
	1171	8 x 8 cm (approx.)	100 pcs	30 packs	01363822
	1173	10 x 10 cm (approx.)	100 pcs	20 packs	12386275



### Double sterile, with self-adhesive counting cards

Pack of 10 pcs In paper bag	4089	6 x 6 cm (approx.)	60 x 10 pcs	2 packs	13165731
	4311	8 x 8 cm (approx.)	60 x 10 pcs	2 packs	13165754
	4339	12 x 12 cm (approx.)	50 x 10 pcs	2 packs	13165777
Pack of 10 pcs In counting box	4343	4 x 4 cm (approx.)	20 x 10 pcs	9 packs	13585373
	4340	6 x 6 cm (approx.)	20 x 10 pcs	9 packs	12731105
	4341	8 x 8 cm (approx.)	20 x 10 pcs	9 packs	12731111
	4342	12 x 12 cm (approx.)	10 x 10 pcs	4 packs	12731128

## BeeSana® X-Ray pointed swabs

Our pointed swabs made from pure cotton

### Features

- gauze swab, 100% cotton (EN 14079 standard)
- 20-thread
- tapered, formed from just one piece of gauze
- woven-in X-ray contrast thread (BaSO<sub>4</sub>)

### Usage

For blood and exudate absorption, especially in ENT applications.  
For staunching bleeding or as a tamponade.



### Non-sterile

REF	Size (width x length, unrolled)	Length	Pack contents	Shipping unit	PZN
1238	Large, 15 x 16 cm (approx.)	13.5 cm (approx.)	500 pcs	20 packs	13509403



### Double sterile, with self-adhesive counting cards

	REF	Size (width x length, unrolled)	Length	Pack contents	Shipping unit	PZN
Pack of 10 pcs	4004	Small, 6 x 8 cm (approx.)	8 cm (approx.)	60 x 10 pcs	8 packs	13165783
	4161	Large, 15 x 16 cm (approx.)	13.5 cm (approx.)	40 x 10 pcs	8 packs	13165808

## BeeSana® X-Ray gauze balls

Our gauze balls made from pure cotton

### Features

- gauze swab, medical gauze (EN 14079 standard)
- with green retrieval thread
- innovative wrapping system from just one piece of pre-cut gauze
- inner-facing cut edges
- 20-thread
- absorbent

### Usage

For blood and exudate absorption.  
The swab can be safely retrieved using the thread.



### Non-sterile

	REF	Size (width x length)	Shape	Ball	Pack contents	Shipping unit	PZN
Plum-sized	4043	20 x 20 cm (approx.)	Round	Double	200 pcs	5 packs	12390911

## ABE® X-Ray ribbon gauze

Our tamponade made from pure cotton

### Features

- medical gauze (EN 14079 standard)
- 100% cotton
- woven edges
- 20-thread
- 4-ply
- rolled
- woven-in X-ray contrast thread (BaSO<sub>4</sub>)

### Usage

For tamponading wound cavities or fistula tracts. Absorption of bodily fluids in hard-to-reach wounds and body cavities.



### Non-sterile

	REF	Size (width x length)	Pack contents	Shipping unit	PZN
20-thread	1446	6 cm x 2 m (approx.)	40 pcs	10 packs	12386298

### Double sterile, pack of 1 pc. with counting card

	REF	Size (width x length)	Pack contents	Shipping unit	PZN
20-thread	4312	6 cm x 2 m (approx.)	25 x 1 pc.	4 packs	02214231



## BeeSana® Sets – extremely cost-effective yet offering maximum customisation

Our set solutions offer a wide range of combinations of components in a single pack enabling our clients to use our products safely and efficiently. Material use, preparation time, waste volume and logistics costs are all reduced with our sets that come in a wide range of individual components ensuring high quality patient care. All components are ready for use in the amount precisely allocated to the patient and procedure concerned.

Of course, all our sets for the operating theatre are fitted with double, self-adhesive counting cards to ensure reliable recording of the details in the patient files. The use of set solutions leads to an optimisation of process and great costeffectiveness for the client. Set solutions are available both as prefabricated standard sets and as bespoke sets put together according to the client's wishes.

### Bespoke set solutions

The contents of each set are put together according to the specific requirements of our clients. Our trained and competent sales staff work with our clients to understand the requirements for special procedures and come up with tailored solutions to meet the highest demands. Clients can draw from our entire product range, comprising bandages, coverings, protective clothing and other materials.

### Standardised set solutions

Given our years of experience in creating bespoke set solutions for our clients, we have been able to develop a range of standard sets for various procedures both in the operating theatre and in the ward. Combinations that have been specially adapted to the requirements of various procedures make it easier for our clients to prepare and carry out standard procedures.



## BeeSana® Ward catheter set

### Features

The ward catheter set is used for hygienic intermittent bladder catheterisation and is sterile, to enable immediate use on the ward.



### Usage

For bladder catheterisation.

### Sterile

REF	Contents	Pack contents	Shipping unit	PZN
4265	1x rectangular dish, 650 ml 1x swab dish, 125 ml 2x gauze dressing pads, 7.5 x 7.5 cm, 8-ply 6x gauze swabs, 20 x 20 cm, plum size 1x pair of plastic tweezers 1x slit/hole cloth, 60 x 60 cm Wrapped in: 1 cover cloth, 60 x 60 cm	21 sets	2 packs	04136073
4267	1x rectangular dish, 650 ml 1x swab dish, 125 ml 2x gauze dressing pads, 7.5 x 7.5 cm, 8-ply 6x gauze swabs, 20 x 20 cm, plum size 2x pairs of examination gloves, size L, latex 1x pair of plastic tweezers 1x slit/hole cloth, 60 x 60 cm Wrapped in: 1 cover cloth, 60 x 60 cm	21 sets	2 packs	04136067



## BeeSana® PEG bandage set

### Features

The PEG bandage set is used to hygienically position a bandage during percutaneous endoscopic gastrostomy and is sterile for immediate use on the ward.



### Usage

For external wound care.

### Sterile

REF	Contents	Pack contents	Shipping unit	PZN
4653	2x non-woven slit dressing pads, 7.5 x 7.5 cm, 4-ply, 30 g 2x non-woven dressing pads, 7.5 x 7.5 cm, 6-ply, 30 g 1 x ABE® non-woven fixing material 12 x 15 cm	10 sets	30 packs	14399023
4577	2x slit non-woven dressing pads, 7.5 x 7.5 cm, 30 g, 4-ply 2x non-woven dressing pads, 7.5 x 7.5 cm, 30 g, 6-ply 1 x ABE® non-woven fixing material 12 x 15 cm	40 sets	8 packs	01333962
4589	2x slit non-woven dressing pads, 7.5 x 7.5 cm, 30 g, 4-ply 2x non-woven dressing pads, 7.5 x 7.5 cm, 30 g, 6-ply 1 x ABE® non-woven fixing material 12 x 15 cm	30 sets	8 packs	04152190
4836	2x gauze dressing pads, 7.5 x 7.5 cm, 8-ply 1x drainage tube dressing pad (aluminium/non-woven) 6 x 7 cm, with slit/hole 1 x ABE® non-woven fixing material 12 x 15 cm 1 x fixation holder, size 2, small, with hook-and-pile fastener	40 sets	8 packs	15995170

## BeeSana® Bandage change set

### Features

The bandage change set is used for hygienically changing bandages and is sterile for immediate use on the ward.



### Sterile

REF	Contents	Pack contents	Shipping unit	PZN
4820	2x gauze dressing pads, 10 x 10 cm, 8-ply 4x gauze swabs, 20 x 20 cm, plum size	25 sets	4 packs	00549660
5413	2x gauze dressing pads, 5 x 5 cm, 16-ply 4x gauze dressing pads, 10 x 10 cm, 12-ply	20 sets	4 packs	–

### Usage

For external wound care.

## BeeSana® Suture removal set

### Features

The suture removal set is used for the hygienic removal of sutures and is sterile, so as to enable its immediate use on the ward.



### Sterile

REF	Contents	Pack contents	Shipping unit	PZN
5450	4x gauze dressing pads, 10 x 10 cm, 12-ply 4x gauze swabs, 25 x 26 cm, egg size 1x pair of suture removal scissors, 11.5 cm 2x disposable tweezers, 130 mm, blue 1x protective covering, 45 x 75 cm, 2-ply, green	15 sets	4 packs	–

### Usage

For external wound care.

## BeeSana® Ward protection set

### Features

The ward protection set is designed to prevent patient-to-patient contamination and the transfer of an infection to medical and/or nursing staff.



### Non-sterile

REF	Contents	Pack contents	Shipping unit	PZN
4074	1x PP gown 1x nurse's cap 1x surgical mask 2x CPE shoe covers 2x examination gloves, latex, size M	80 sets	1 pack	–

### Usage

To protect against germs.



## Our incontinence products are designed to ensure quality of life for patients

A variety of conservative and surgical treatments are now available to cure or at least alleviate urinary and faecal incontinence. In cases where patients can no longer be treated, modern incontinence products provide a solution. When developing and producing Meditrade incontinence products, we strive to achieve well-designed, ergonomic solutions that provide discretion for the user.

Our comprehensive range of products enables cost-effective use with a wide range of sizes, absorbency levels and qualities, always following the principle of achieving targeted and acceptable patient care.



# Superslip®

Anatomically shaped incontinence briefs with flexible waist

### Features

Incontinence briefs with high liquid absorption, elastic back closure and double hook-and-pile fastener. Hydrophobic side seams and a dampness indicator to show when it's time to change.

### Usage

For urinary and faecal incontinence, and for people requiring all-round care with reduced mobility.



## Superslip® Classic

REF	Size	Colour code	Hip measurement	Absorbency (ISO 11948-1 standard)	Pack contents	Shipping unit	PZN	Assistive technology no.
2431	Medium	White/blue	80–115 cm	1,850 ml (approx.)	20 pcs	4 packs	07757522	15.25.03.1118
2432	Large	White/green	115–150 cm	2,100 ml (approx.)	20 pcs	4 packs	07757551	15.25.03.2112

## Superslip® Comfort

2406	Small	White	50–80 cm	1,800 ml (approx.)	15 pcs	6 packs	07757485	15.25.03.0066
2407	Medium	White/blue	80–115 cm	2,450 ml (approx.)	15 pcs	8 packs	07757580	15.25.03.1119
2408	Large	White/green	115–150 cm	2,600 ml (approx.)	15 pcs	6 packs	07757605	15.25.03.2113
2414	X-Large	White/violet	>150 cm	2,700 ml (approx.)	15 pcs	4 packs	07757634	15.25.03.2312

## Superslip® Extra

2436	Large	White/green	115–150 cm	3,100 ml (approx.)	15 pcs	4 packs	09747019	15.25.03.2314
------	-------	-------------	------------	--------------------	--------	---------	----------	---------------

## Superslip® Super

2403	Medium	White/blue	80–115 cm	3,000 ml (approx.)	15 pcs	4 packs	09754143	15.25.03.1278
2404	Large	White/green	115–150 cm	3,250 ml (approx.)	15 pcs	4 packs	00351917	15.25.03.2313

# Superslip® Air Comfort

Anatomically shaped incontinence briefs with flexible waist

### Features

Incontinence briefs with high liquid absorption, elastic back closure and double hook-and-pile fastener. Hydrophobic side seams and a dampness indicator to show when it's time to change.

### Usage

For urinary and faecal incontinence, and for people requiring all-round care with reduced mobility.



REF	Size	Colour code	Hip measurement	Absorbency (ISO 11948-1 standard)	Pack contents	Shipping unit	PZN	Assistive technology no.
2437	Medium	White/blue	80–115 cm	1,750 ml (approx.)	15 pcs	8 packs	04917043	15.25.03.1276
2438	Large	White/green	100–150 cm	2,000 ml (approx.)	15 pcs	6 packs	04918189	15.25.03.2311

# BeeSana® Medi-Pants

Discreet disposable incontinence briefs

### Features

Incontinence briefs for those with a weak bladder, easy to put on and take off like normal underwear. Provides unrestricted freedom of movement and a normal level of independence. For women and men, disposable.

### Usage

For patients suffering from a weak bladder. For mobile or dementia patients without any additional restrictions to daily life. Suitable for use with toilet training. To be worn like normal underwear.



REF	Size	Hip measurement	Pack contents	Shipping unit	PZN
2411	Medium	70–110 cm (approx.)	14 pcs	6 packs	11141726
2412	Large	100–140 cm (approx.)	14 pcs	6 packs	11141732
2413	X-Large	120–165 cm (approx.)	14 pcs	6 packs	11141749

## Light fixation briefs

For keeping incontinence pads in place

### Features

Latex-free fixation briefs with underwear-like fabric structure. Size indicated by coloured strips. Washable at 60 °C and can be tumble-dried on the delicates setting.

### Usage

For keeping incontinence pads and liners in place.



REF	Size	Colour	Hip measurement	Pack contents	Shipping unit	PZN
2221	Medium	Green	80–110 cm (approx.)	50 pcs	6 packs	11513151
2222	Large	White	90–130 cm (approx.)	50 pcs	6 packs	11513211
2223	X-Large	Blue	110–170 cm (approx.)	50 pcs	6 packs	11513228
2224	XX-Large	Orange	120–180 cm (approx.)	50 pcs	6 packs	11513234

## Panty fixation briefs

With legs for keeping incontinence pads in place

### Features

Latex-free fixation briefs, size indicated by coloured strips. Washable at 60 °C and can be tumble-dried on the delicates setting.

### Usage

For keeping incontinence pads and liners in place.



REF	Size	Colour	Hip measurement	Pack contents	Shipping unit	PZN
0711500301	Medium	Blue	70–100 cm (approx.)	100 pcs	6 packs	09081879
0711500501	Large	Brown	85–120 cm (approx.)	100 pcs	6 packs	09081885
0711500701	X-Large	Green	105–130 cm (approx.)	100 pcs	6 packs	09081891
0711500901	XX-Large	Green/blue	115–140 cm (approx.)	100 pcs	6 packs	09081916

# Liberty Light

Anatomically formed incontinence pad

### Features

Anatomical fit with hydrophobic seams in the crotch area. The secure adhesive strip is used to fix this pad into underwear.

### Usage

For patients with a weak bladder, for protection after vaginal procedures or as a pad for post-natal women.



REF	Size	Colour code	Dimensions (approx.)	Absorbency (ISO 11948-1 standard)	Pack contents	Shipping unit	PZN	Assistive technology no.
2052	Normal	Rose	28 x 10.5 cm	min. 350 ml	22 pcs	10 packs	07305197	15.25.01.5175
2053	Extra	Blue	31.5 x 12 cm	min. 400 ml	18 pcs	10 packs	07305180	15.25.01.5177
2054	Super	Green	43 x 22 cm	min. 850 ml	30 pcs	5 packs	07305174	15.25.01.5176

# Superform<sup>®</sup>

Anatomically formed incontinence pad with double cushion

### Features

Incontinence pad with anatomical fit. Hydrophobic side seams as a leak barrier and a dampness indicator to show when it's time to change.

### Usage

Suitable for both bed-bound and mobile incontinence patients.

**Note on product combination:** To fix the pad in place, we recommend using our washable fixation briefs, PANTY or LIGHT.



REF	Size	Colour code	Dimensions	Absorbency (ISO 11948-1 standard)	Pack contents	Shipping unit	PZN	Assistive technology no.
2047	Classic	Blue	63 x 30 cm	1,200 ml (approx.)	30 pcs	5 packs	07520263	15.25.01.0076
2048	Plus	Blue	63 x 30 cm	1,350 ml (approx.)	30 pcs	5 packs	07414302	15.25.01.1096
2049	Extra	Green	70 x 37 cm	1,450 ml (approx.)	30 pcs	4 packs	07414288	15.25.01.1105
2050	Super	Green	70 x 37 cm	1,850 ml (approx.)	20 pcs	4 packs	07414271	15.25.01.1065
2051	Maxi	White	70 x 37 cm	2,450 ml (approx.)	20 pcs	4 packs	07414070	15.25.01.2101

# Superform<sup>®</sup> Air

Extra light anatomically formed incontinence pad with double cushion

### Features

Incontinence pad with anatomical fit. Hydrophobic side seams as a leak barrier and a dampness indicator to show when it's time to change.

### Usage

Suitable for both bed-bound and mobile incontinence patients.

**Note on product combination:** To fix the pad in place, we recommend using our washable fixation briefs, PANTY or LIGHT.



REF	Size	Colour code	Absorbency (ISO 11948-1 standard)	Pack contents	Shipping unit	PZN	Assistive technology no.
20481	Plus	White/blue	1,250 ml (approx.)	30 pcs	5 packs	08801840	15.25.01.0098
20491	Extra	White/yellow	1,400 ml (approx.)	30 pcs	4 packs	08801857	15.25.01.1125
20501	Super	White/green	1,800 ml (approx.)	20 pcs	4 packs	08801863	15.25.01.2141
20511	Maxi	White/violet	2,200 ml (approx.)	20 pcs	4 packs	08834851	15.25.01.2142

# BeeSana<sup>®</sup> Absorbent cushion

Pad for incontinence and nursing care

### Features

Latex-free absorbent cushion from cellulose flock with non-woven sleeve and outer film. This super-absorbent material protects underwear and clothing from soiling. Simple to apply in combination with Light or Panty fixation briefs.

### Usage

For incontinence/loose stools, for catheter wearers, bed-ridden and mobile patients. For use as a universal absorbent cushion in nursing care.



	REF	Size	Pack	Pack contents	PZN
With adhesive strip	2022	12 x 38 cm	50 pcs	8 packs	02214142
Without adhesive strip	2024	18 x 40 cm	60 pcs	3 packs	02214165



# BeeSana® Absorbent pad

Our large pad for incontinence care

### Features

Latex-free absorbent cushion from cellulose flock with non-woven sleeve. Simple to apply in combination with Light or Panty fixation briefs.

### Usage

For incontinence for bed-ridden people.



REF	Size	Pack contents	Shipping unit	PZN
2119	16 x 55 cm	20 pcs	5 packs	02214219

# BeeSana® Rectangular pad

Our absorbent pad for incontinence

### Features

This latex-free absorbent cushion from cellulose flock with a non-woven sleeve is simple to apply in combination with Light or Panty fixation briefs.

### Usage

For incontinence for bed-ridden people.

For post-natal women in the period after giving birth.

Suitable for protection after vaginal procedures.



REF	Size	Pack contents	Shipping unit	PZN
2101	35 x 11 cm	25 pcs	8 packs	12147808

# BeeSana® Male incontinence pad

Anatomically formed pad for men

### Features

Latex-, PVC- and DEHP-free. Ergonomic cut to fit the male anatomy. Not noticeable when worn. With adhesive strips for fixing to underwear.

### Usage

For men with urinary incontinence, such as after surgery in the genital area.



REF	Size	Pack contents	Shipping unit	PZN
2020	8.5 x 40 cm	21 pcs	7 packs	05880164

# BeeSana® Female incontinence pad

Our hygienic female incontinence pad

### Features

Latex-free absorbent cushion from absorbent cellulose flock with non-woven sleeve.

### Usage

As protection after vaginal procedures or as a pad for post-natal women after childbirth.



REF	Size	Pack contents	Shipping unit	PZN
2118	7 x 23 cm (approx.)	20 pcs	18 packs	01091275

# Layered undersheet

An industry-standard undersheet for meeting patient needs

## Features

Undersheet with absorbent core made from environmentally friendly, 100% non-bleached recycled cellulose for fluid absorption. The dry non-woven material with outer film makes use easy in any situation.

## Usage

For live-in nursing, undersheets help to make nursing easier by simplifying day-to-day care work and, last but not least, reducing the 'laundry mountain' by protecting the bed covers. We can supply the sheets as layered or flock undersheets in a range of different qualities.



REF	Size	Number of layers	Shipping unit	PZN
104606M	40 x 60 cm	6-ply	200 pcs	02426826
104608M	40 x 60 cm	8-ply	200 pcs	00559612
104610M	40 x 60 cm	10-ply	200 pcs	00559629
104612M	40 x 60 cm	12-ply	100 pcs	06178087
104620M	40 x 60 cm	20-ply	100 pcs	00559693
106606M	60 x 60 cm	6-ply	150 pcs	00560213
106610M	60 x 60 cm	10-ply	100 pcs	00560236
106906M	60 x 90 cm	6-ply	100 pcs	00560259
106908M	60 x 90 cm	8-ply	100 pcs	00560331
106910M	60 x 90 cm	10-ply	100 pcs	00560609
106912M50	60 x 90 cm	12-ply	50 pcs	00561046
106920M	60 x 90 cm	20-ply	50 pcs	00561052

# BeeSana® Comflock Super

Undersheet made from cellulose flock

## Features

Undersheet with absorbent core made from cellulose flock for fluid absorption. The dry non-woven material with outer film makes use easy in any situation.

## Usage

For live-in nursing, undersheets help to make nursing easier by simplifying day-to-day care work and, last but not least, reducing the 'laundry mountain' by protecting the bed covers. We can supply the sheets as layered or flock undersheets in a range of different qualities.



REF	Size	Fluid absorbency	colour	Pack contents	Shipping unit	PZN
2064	40 x 60 cm	750 ml ±10%	white	50 pcs	6 packs	04136185
2065	60 x 60 cm	900 ml ±10%	white	25 pcs	6 packs	04136104
2066	60 x 90 cm	1,590 ml ±10%	blue	25 pcs	6 packs	04136156

# Nursing pad

Our pad to protect underwear during breastfeeding

## Features

Anatomically formed nursing pad. The cellulose core absorbs leaks of milk while offering a discreet fit.



## Usage

For absorbing fluids and to protect underwear during post-natal breastfeeding.

## Non-sterile

REF	Size	Pack contents	Shipping unit	PZN
2068	10.5 cm	30 pcs	12 packs	02097937





## Disinfection and hygiene products with a broad spectrum of action

The core elements of our broad range of products reduce the risk of infection, standardise prevention methods and provide lasting skin care. Please note: Our product recommendations are based on scientific testing and are offered according to the best of our knowledge. More detailed recommendations – such as in terms of material compatibility – can only be provided separately.

Our recommendations do not exempt our customers from testing our products independently to ensure they are suitable for the intended purposes and procedures. We cannot, therefore, assume any liability.

# Ethasept®

Our alcohol-based hand sanitiser

### Product type

Ethasept® is a ready-to-use, alcohol-based, hypoallergenic disinfectant without perfumes or colourings, and is suitable for disinfection applications in human hygiene.

### Applications

For hygienic and surgical hand disinfection in all areas of human hygiene.

### Usage

Clean heavily soiled hands before disinfection and let wet hands dry before starting disinfection. Ready-to-use. Keeps hands/forearms moist during application. Observe standard procedures for hygienic (EN 1500) and surgical (EN 12791) hand disinfection.

### Product properties

- fast-acting and broad-spectrum, works after just 15 s
- an effective bactericide, yeasticide, tuberculocide and virucide
- effective against MRSA
- dermatologically tested: "Excellent"
- listed by the German Association for Applied Hygiene (VAH)

### Shelf life

3 years (see label). Shelf life after breaking seal: 12 months. Keep containers tightly closed. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).

### Composition

Active substance in 100 g: 73.5 g ethanol.

### Regulatory requirements

Ethasept® is not a medicinal product and is subject to biocide legislation. Use biocides safely. Always read the label and product information before use.

Type 1 biocide product ("biocide products for human hygiene"), biocide reg. no. N-44236.

### Contact time

Activities	15 sec	30 sec	1 min / 60 sec	1.5 min / 90 sec	2 min / 120 sec
Hygienic hand disinfection: (EN 1500)	+				
Surgical hand disinfection: (EN 12791)				+	
Virucide					+
Virucide (EN 14476)					
Polio (high load)					+
Adeno (high load)			+		

REF	Pack contents	Shipping unit	PZN
01001UK	100 ml	42 bottles	-
01004UK	500 ml	18 bottles	-
01005UK	1 l	10 bottles	-

# Ethasept® gel

Hand disinfection

### Composition

73.5g Ethanol Ethasept Gel is ready-to-use alcohol-based hand sanitiser gel for hygienic hand disinfection.

### Product Type

Ethasept® gel is a ready-to-use, alcohol-based, hypoallergenic disinfectant without perfumes or colourings, and is suitable for disinfection applications in human hygiene.

### Application

Ready-to-use. Thoroughly rub hands with at least 2 x 3 ml onto the hands during 60 seconds, undiluted and without adding water. Observe standard methods for hygienic hand disinfection (EN 1500).

### Product properties

- Works quickly and comprehensively after just 60 seconds.
- Bactericidal, yeasticidal and limited virucidal activity
- Effective against MRSA
- Dermatologically tested and rated "Excellent"
- Listed by the German Association for Applied Hygiene (VAH)


### Shelf life

3 years (see label). Shelf life after breaking seal: 12 months. Keep containers tightly closed. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).

### Regulatory requirements

Ethasept® gel is not a medicinal product and is subject to biocide legislation. Use biocides safely. Always read the label and product information before use.

REF	Box contents	Shipping unit	PZN
02001UK	100 ml	42 bottles	-
02002UK	500 ml	18 bottles	-
02003UK	1 l	10 bottles	-





## Alcohol-based surface disinfection

# Medizid® Rapid+

High-performance, ready-to-use rapid disinfectant for alcohol-resistant surfaces and objects

### Product type

Medizid® Rapid+ is a ready-to-use, alcohol-based disinfectant for application using a cloth, with broad-spectrum effectiveness for the disinfection of alcohol-resistant surfaces in healthcare facilities. Ideal for rapid disinfection between treatments and procedures, especially when practical situations mean contact times need to be kept short.

### Applications

For the rapid disinfection of alcohol-resistant surfaces in the patient's immediate environment in hospitals, general practice, outpatient units, nursing homes and manufacturing.

### Usage

Remove visible soiling before performing disinfection. Wet the surfaces completely with the solution and use a hygienic, clean disposable cloth to distribute the disinfectant solution as required. Remember to observe the contact times as necessary. For professional use only. Not suitable for the disinfection of invasive medical devices. Test material compatibility on an unobtrusive portion if necessary.

### Product properties

- aldehyde-free
- effectiveness:
  - effective bactericide (incl. MRSA), yeasticide
  - effective fungicide
  - effective mycobactericide
  - effective virucide
- listed by the German Association for Applied Hygiene (VAH)

### Shelf life

3 years (see label). Shelf life after breaking seal: 28 days. Keep containers tightly closed. Use moist wipes only. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).

### Composition

Active substances in 100 g:  
63,5 g Ethanol


### Regulatory requirements

Class IIa medical device

### Contact time

Activities	fully wetted
<b>As per VAH guidelines:</b>	
Bactericide, yeasticide incl. 4-field test:	1 min. (low load)
Bactericide, yeasticide without mechanical system:	2 min. (high load)
Mycobactericide without mechanical system:	5 min. (high load)
As per EN 16615 (4-field test):	2 min. (high load)

REF	Pack contents	Shipping unit	PZN
01002+	250 ml	12 bottles	16015216
00975+	500 ml	18 bottles	16015222
00976+	1 l	12 bottles	16015239
00977+	5 l	2 canisters	16015245
<b>Accessories (see page 147)</b>			
	Canister tap		
	Canister key		
	Spray nozzle for 1 l bottle		



# Medizid® Rapid+ disinfectant wipes

Disinfectant wipes for alcohol-resistant surfaces on medical devices

### Product type

Medizid® Rapid+ is a ready-to-use, alcohol-based disinfectant for application using a cloth, with broad-spectrum effectiveness for the disinfection of alcohol-resistant surfaces in healthcare facilities. Ideal for rapid disinfection between treatments and procedures, especially when practical situations mean contact times need to be kept short.

### Applications

For the rapid disinfection of alcohol-resistant surfaces in the patient's immediate environment in hospitals, general practice, outpatient units, nursing homes, manufacturing and public organisations.

### Usage

Remove visible soiling before performing disinfection. Then wipe down the surfaces carefully with the impregnated cloth. Ensure complete wetting of the surface and observe the contact times as necessary. For professional use only. Not suitable for the disinfection of invasive medical devices. Test material compatibility on an unobtrusive portion if necessary.

### Product properties

- aldehyde-free
- effective bactericide, yeasticide, mycobactericide, fungicide, virucide
- effective against MRSA
- listed by the German Association for Applied Hygiene (VAH)


### Shelf life

3 years (see label). Shelf life after breaking seal: 28 days. Keep containers tightly closed. Use moist wipes only. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).

### Contact time

Activities	fully wetted	<b>As per EN 16615 (4-field test) - practical test (for wipes):</b>	
<b>As per VAH guidelines:</b>		Bio wipes	1 min. (high load)
Bactericide, yeasticide incl. 4-field test:	1 min. (low load)	Standard wipe	2 min. (high load)
Bactericide, yeasticide without mechanical system:	2 min. (high load)	Long-term bio wipe test (28, 60 and 90 days)	1 min. (high load)
Mycobactericide without mechanical system:	5 min. (high load)		

REF	Article	Wipe size	Pack contents	Shipping unit	PZN
00941+	Wipes in a PE dispenser tin	15 x 18 cm (approx.)	150 wipes	12 tins	16015268
00942+	Refill tin	15 x 18 cm (approx.)	150 wipes	20 bags	16015251
00946+	Wipes in a container	29 x 30 cm (approx.)	70 wipes	1 container	16015274
00947+	Refill container	29 x 30 cm (approx.)	70 wipes	6 rolls	16015280



### Composition

Active substances in 100 g:  
0.05 g didecylidimethylammonium chloride, 63.5 g ethanol

### Regulatory requirements

Class IIa medical device.



# Medizid® Rapid QF Flowpack

Alcohol-based hygiene wipes for the disinfection of medical device surfaces and objects

## Product type

Medizid® Rapid QF Flowpack is impregnated with the very fast-acting fluid from the alcohol-based product Medizid® Rapid. The hygiene wipes are therefore a good choice for disinfecting medical device surfaces and objects used in medical applications, and in the patient's immediate environment between treatments and interventions.

## Applications

For all alcohol-resistant medical devices, surfaces and objects, especially for device surfaces/objects in the patient's immediate environment in hospitals, general practice, outpatient units and nursing units.

## Usage

Medizid® Rapid QF disinfectant wipes (flowpack) are ready-to-use. The surfaces and objects to be disinfected are wiped down with the Medizid® Rapid wipe, ensuring that the surface/object is wetted thoroughly. Exercise care with alcohol-susceptible materials. Do not use for the disinfection of invasive devices.

## Product properties

- aldehyde-free
- effective bactericide, yeasticide, mycobactericide and fungicide, "Limited virucide PLUS" rating
- effective against MRSA
- listed by the German Association for Applied Hygiene (VAH)

## Shelf life

3 years (see label). Shelf life after breaking seal: 28 days. Keep containers tightly closed. Use moist wipes only. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).


## Composition

Active substances in 100 g: 50 g ethanol

## Regulatory requirements

Medizid® Rapid QF Flowpack is not a medicinal product and is subject to biocide legislation. Use biocides safely. Always read the label and product information before use. Type 1 biocide product ("biocide products for human hygiene"), biocide reg. no. N-75574.

REF	Wipe size	Pack contents	Shipping unit	PZN
00976FP	22 x 20 cm (approx.)	80 pcs	12 packs	–



## Contact time

Activities	30 sec	1 min	5 min
Activity as per EN 16615 practical test (for wipes), bactericide, yeasticide incl. 4-field test (high load): Bio wipes:		+	
<b>Activities for immersion solution:</b>			
<b>Contact times for hygiene applications</b>			
Yeasticide (EN 1650, low load)		+	
Bactericide, yeasticide (EN 13727, EN 13624, high load)		+	
Mycobactericide as per VAH method 9 (high load)			+
Tuberculocide (EN 14368, low load)		+	
<b>Contact times for food applications</b>			
Bactericide (EN 1276, low load)	+		
Yeasticide (EN 1650, low load)	+		
Bactericide/yeasticide (EN 13697, low load)			+
<b>Virucidal activities (EN 14476)</b>			
Limited virucide PLUS		+	
Vaccinia		+	
Adenovirus		+	
Norovirus		+	

## Alcohol-free surface disinfection

# Medizid® AF+

Concentrate for surface disinfection

## Product type

Medizid® AF+ is a liquid, aldehyde- and alcohol-free concentrate for the disinfection of surfaces not resistant to alcohol on medical devices, and especially device surfaces and objects in the patient's immediate environment.

## Applications

Ideal for areas with heavy soiling or odours, such as operating theatre tables, bed frames and surfaces in daily use in medical facilities, manufacturing and public organisations. Also suitable for surfaces not resistant to alcohol. Test material compatibility on an unobtrusive portion if necessary.

## Usage


Remove visible soiling before disinfection; wet surfaces completely and keep surfaces moist during the entire contact time. Use dosing aids and follow instructions; example of use: 10 l of a 1% solution = 9.9 l water + 100 ml Medizid® AF+.

## Product properties

- alcohol-free concentrate
- aldehyde-free
- effective bactericide, yeasticide, tuberculocide, "Limited virucide PLUS" rating
- high material compatibility
- listed by the German Association for Applied Hygiene (VAH)

## Contact time

Activity as per EN standards and VAH guidelines:	1 min	5 min	30 min	60 min
Bactericide, yeasticide (EN 13626, EN 16615; low load)		0.5%		
Bactericide, yeasticide (EN 16615; high load)	0.5%			
Tuberculocide (EN 14348; high load)				5%
Mycobactericide (EN 14348; incl. M. terrae and M. avium; high load)			2%	1%
<b>Virucidal activities (EN 14476)</b>				
<b>Limited virucide PLUS</b>				4%
Vaccinia Ankara (incl. HIV, HBV, HCV, influenza viruses, high load)		3%	1%	
Adeno-/noroviruses (high load)				4%
Rotaviruses (no load)			1%	
SV40 (high load)				2%



## Shelf life

3 years (see label). Shelf life after breaking seal: 12 months. Keep container tightly closed. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).

## Composition

Active substances in 100 g: 2 g 2-phenoxyethanol, 7 g alkyl (C12-16) dimethyl benzyl ammonium chloride, 5 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine

## Regulatory requirements

Class IIa medical device

REF	Pack contents	Shipping unit	PZN
00972+	2 l	6 bottles	16015185
00971+	5 l	2 canisters	16015191
<b>Accessories (see page 147)</b>			
Refillable dosing bottle			
Canister pump for dosing			
Canister tap			
Canister key			

# Medizid® Alpha+

Alcohol-free, ready-to-use disinfection for wipe-resistant surfaces on medical devices

## Product type

Medizid® Alpha+ is a ready-to-use, aldehyde- and alcohol-free spray and wipe disinfectant solution with broad-spectrum effectiveness for the disinfection of surfaces in healthcare facilities.

## Applications

For healthcare facilities, sanitary/bathroom use and nursing homes.

## Usage

Remove visible soiling before performing disinfection. Wet the surfaces completely with the solution and use a hygienic, clean disposable cloth to distribute the disinfectant solution as required. Keep the surface moist during the entire contact time. Wipe away residues of the disinfectant after the expiry of the contact time, either with a cloth soaked in water (at least drinking water quality) or a hygienically clean cloth. For professional use only. Always use undiluted.

## Product properties

- Comprehensive bactericide and yeasticide
- "Limited virucide PLUS" rating
- Material compatibility for surfaces sensitive to alcohol
- Listed by the German Association for Applied Hygiene (VAH)

## Shelf life

3 years (see label). Shelf life after breaking seal: 12 months. Keep container tightly closed. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).

## Contact time

Activities	
<b>As per EN/methods from VAH/DGHH</b>	
Bactericide, yeasticide (EN 16615; high load)	1 min
Bactericide, yeasticide (EN 13624, EN 13727, high load)	1 min
<b>Virucidal activities (EN 14476)</b>	
Limited virucide PLUS	60 min
Vaccinia Ankara (incl. HIV, HBV, HCV, high load)	1 min
Noroviruses (high load)	30 min
Adenoviruses (high load)	60 min

REF	Pack contents	Shipping unit	PZN barcode
00964+	1 l	12 bottles	16015156

Accessories (see page 147)

Spray head



## Composition

Active substances in 100 g:  
0.1 g 2-phenoxyethanol, 0.35 g alkyl (C12-16)dimethyl benzyl ammonium chloride, 0.25 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine

## Regulatory requirements

Class IIa medical device

# Medizid® Alpha+ XXL disinfectant wipes

Alcohol-free disinfectant wipes for wipe-resistant surfaces on medical devices

## Product type

Medizid® Alpha+ is a ready-to-use, aldehyde- and alcohol-free wipe disinfectant solution with broad-spectrum effectiveness for the disinfection of surfaces in healthcare facilities.

## Applications

For healthcare facilities, sanitary/bathroom use and nursing homes.

## Usage

Remove visible soiling before performing disinfection. Wet the surfaces completely with the solution and use a hygienic, clean disposable cloth to distribute the disinfectant solution as required. Keep the surface moist during the entire contact time. Wipe away residues of the disinfectant after the expiry of the contact time, either with a cloth soaked in water (at least drinking water quality) or a hygienically clean cloth. For professional use only. Always use undiluted.

## Product properties

- ready-to-use
- bactericide and yeasticide, "Limited virucide PLUS" rating
- material compatibility for surfaces sensitive to alcohol
- listed by the German Association for Applied Hygiene (VAH)

## Shelf life

3 years (see label). Shelf life after breaking seal: 12 months. Keep container tightly closed. Use moist wipes only. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).

## Contact time

Activities	
<b>As per EN/methods from VAH/DGHH</b>	
Bactericide, yeasticide (EN 16615; high load)	1 min
Bactericide, yeasticide (EN 13624, EN 13727, high load)	1 min
<b>Virucidal activities (EN 14476)</b>	
Limited virucide PLUS	60 min
Vaccinia Ankara (incl. HIV, HBV, HCV, high load)	1 min
Noroviruses (high load)	30 min
Adenoviruses (high load)	60 min

REF	Article	Wipe size	Pack contents	Shipping unit	PZN
00969+	Wipes in a container	26 x 29 cm	70 wipes	1 container	16015179
00970+	Refill pack	26 x 29 cm	70 wipes	6 rolls	16015162



## Composition


Active substances in 100 g:  
0.1 g 2-phenoxyethanol, 0.35 g alkyl (C12-16)dimethyl benzyl ammonium chloride, 0.25 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine

## Regulatory requirements

Class IIa medical device

# Medizid® Alpha+ Flowpack

## Surface disinfection



### Product type

Medizid® Alpha+ Flowpack are ready-to-use, aldehyde- and alcohol-free hygiene wipes with a broad-based mechanism of action for the disinfection of surfaces not resistant to alcohol, medical devices and objects in healthcare facilities. Ideal when contact times need to be kept short for practical reasons.

### Applications

Suitable for the rapid disinfection of all surfaces and medical devices not resistant to alcohol in the patient's immediate environment in hospitals, general practice, outpatient units and nursing homes, as well as the food industry, manufacturing and public organisations. Test material compatibility on an unobtrusive portion if necessary. Not for use in the final disinfection of semi-critical and invasive medical devices

### Usage

Remove visible soiling before performing disinfection. Then wipe down the surfaces carefully with the impregnated cloth. Ensure wetting is complete and keep moist during the entire contact time.

### Product properties

- comprehensive bactericide and yeasticide
- "Limited virucide PLUS" rating
- material compatibility for surfaces sensitive to alcohol
- wipes made from 100% biodegradable cellulose fibres
- listed by the German Association for Applied Hygiene (VAH)



### Shelf life

3 years (see label). Shelf life after breaking seal: 28 days. Keep container tightly closed. Use moist wipes only. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).

### Composition

Active substances in 100 g: 0.10 g 2-phenoxyethanol 0.10 g, 0.35 g alkyl (C12-16)dimethyl benzyl ammonium chloride, 0.25 g N-(3-aminopropyl)-N-dodecylpropane-1,3-diamine

### Regulatory requirements

Class IIa medical device

### Contact time

As per VAH guidelines	
Bactericide, yeasticide (EN 16615; high load)	1 min
Bactericide, yeasticide (EN 13624, EN 13727) (high load)	1 min
Virucidal activities (EN 14476)	
Limited virucide PLUS	60 min
Vaccinia Ankara (incl. HIV, HBV, HCV, high load)	1 min
Noroviruses (high load)	30 min
Adenoviruses (high load)	60 min

REF	Wipe size	Pack contents	Shipping unit	PZN
00964+FP	22 x 20 cm (approx.)	80 pcs	12 packs	17198644



# Instrusol® AF +

## Instrument disinfection for manual reprocessing

### Product type

Instrusol® AF + is a highly effective instrument disinfectant for the reprocessing of industry-standard medical and surgical instruments, as well as thermally sensitive anaesthetic accessories and rigid endoscopes, either by an immersion method or in an ultrasonic bath.

### Applications

For all applications where medical instruments are reprocessed by the use of immersion procedures and in ultrasonic cleaning systems.

### Usage

Always follow the usage instructions and reprocessing guidelines offered by instrument manufacturers; note in particular information about the preparatory steps to take, cleaning and disinfection, rinsing and drying, and storage.

### Product properties

- high-performance concentrate
- aldehyde- and phenol-free
- bactericide, yeasticide, tuberculocide
- "Limited virucide" rating
- listed by the German Association for Applied Hygiene (VAH)

### Shelf life

3 years (see label). Shelf life after breaking seal: 12 days. Keep container tightly closed. Warehousing and storage under standard storage conditions (protect from moisture, direct sunlight, heat and frost).

### Composition

Active substances in 100 g: 5 g alkyl dimethyl ethyl ammonium ethyl sulphate, 1 g polyhexamethylene biguanide hydrochloride, 1.8 g cocospropylene guanidinium diacetate

### Regulatory requirements

Class IIb medical device



### Contact time

Activities	5 min	15 min	30 min	60 min
Bactericide, yeasticide (as per VAH) (high load)		4.0 %	3.0 %	2.0 %
Tuberculocide (as per DGHM 2001) (high load)		4.0 %	3.0 %	2.0 %
<b>Virucidal activity:</b>				
Rotavirus (DVV/RKI) (high load)	2.0 %			
SV 40 (DVV/RKI)			3.0 %	
BVDV/Vaccinia (DVV/RKI) (no load)		1.0 %		

### Reprocessing instructions

#### Immersion method

Immerse instruments after use. All surfaces and cavities must be covered and/or filled. Remove instruments, rinse thoroughly with water and dry.

#### Ultrasonic bath

The prescribed frequency range and the prescribed disinfection duration must be maintained. Water temperatures above 45 °C are to be avoided.

#### Product changeover

Instrusol® AF + is not compatible with disinfectant products containing aldehyde. On first use and when changing products, instruments and accessories, baths, reprocessing equipment and any application aids used must be thoroughly cleaned and rinsed off.

Instrusol® AF + ready-to-use solution is protein-resistant and can be used more than once if required. The ready-to-use solution is usable for up to 7 days if instruments immersed are properly pre-cleaned. If clouding or visible contamination occurs, we recommend renewing the ready-to-use solution earlier.

#### Notes

As a ready-to-use solution, Instrusol® AF + is not a hazardous formulation according to the German Ordinance on Hazardous Substances. The ready-to-use solution is not caustic. However, we do recommend taking the normal safety precautions for handling chemicals and wearing the usual kinds of protective clothing.

REF	Pack contents	Shipping unit	PZN
00979D	500 ml	12 bottles	–
00973	2 l	6 bottles	11169699
00974	5 l	2 canisters	11169707

#### Accessories (see page 147)

Canister pump for dosing

Canister tap

Canister key

# Gentle Med<sup>®</sup> wash lotion

## pH skin neutral and gentle washing lotion

### Features

Dermatologically tested and rated a "Excellent" washing lotion for daily, gentle all-over body use, with a fresh scent. Prevents the skin from drying out. Paraben-free. Alkali-free, free of silicone and paraffin oil. Can be used sparingly thanks to the optional metering pump.

### Usage

Hand washing, body washing, showering and bathing.



REF	Pack contents	Shipping unit	PZN
00980D	500 ml	12 bottles	11170314
00981D	1 l	10 bottles	11170320
00982D	5 l	2 canisters	11170337



## Gentle Med<sup>®</sup> skin care lotion (O/W)

The skin care lotion to protect and care for stressed skin, with almond and jojoba oil

### Features

Skin care lotion for sensitive and dry skin, pH skin neutral, enriched with urea. Creates a non-greasy protective film and enriches stressed skin with valuable natural oils (almond and jojoba oil) plus vitamin E. The ingredients urea and hyaluronic acid care for the skin and provide moisture.

Dermatologically tested as "Excellent". PEG-free. Silicone- and paraffin oil-free.

REF	Pack contents	Shipping unit	PZN
01011	100 ml	30 tubes	11170308
00990D	500 ml	12 bottles	11170231



## Gentle Med<sup>®</sup> cream lotion (W/O)

Skin protection balsam with panthenol and vitamin E

### Features

Cream lotion with an advanced skincare formula to protect and promote the regeneration of extremely stressed skin. The Gentle Med<sup>®</sup> cream lotion in this product has been dermatologically tested and rated "Excellent": suitable for all-over body care, it is pH skin neutral and enriched with panthenol and vitamin E. Free from silicone oil, PEG and parabens, with marigold extract.

REF	Pack contents	Shipping unit	PZN
01013	100 ml	30 tubes	11169794
00988D	500 ml	12 bottles	11170053



## Gentle Med<sup>®</sup> active gel

Vitalising and cooling active gel with essential oils

### Features

Active rubbing gel to stimulate the circulation and enliven the skin, while being skin-friendly. Has a deodorising and cooling effect. Well-suited for massage to relieve tension. Dermatologically tested as "Excellent". Alkali-free. Paraben-free. PEG-free. Silicone- and paraffin oil-free. With essential oils from medicinal plants: menthol, camphor and rosemary.

REF	Pack contents	Shipping unit	PZN
00991D	500 ml	12 bottles	11169788



## Gentle Med<sup>®</sup> oil bath

Oil bath with natural oils

### Features

Medical oil bath enriched with high-quality soya oil. Promotes the physiological regeneration of the skin. Offers relief for sensitive and dry skin. Suitable for use in showers and baths.

REF	Pack contents	Shipping unit	PZN
00992D	500 ml	12 bottles	11170260





# Gentle Med®

## skin cleansing and care foam

Gentle cleaning and care

### Features

Fast and gentle cleaning and care of the skin, also for incontinence, patients with a stoma, and for bed-ridden patients. This skin-cleansing care foam contains moisturising and odour absorbing substances. Paraben-free. PEG-free. Silicone- and paraffin oil-free. Contains no allergenic perfumes.

### Usage

To lather, hold the can (valve) down at an angle, thus ensuring that the container is completely emptied. Leave to act briefly and then wipe off with a disposable cloth. A protective film continues to care for the skin and provide a deodorising effect. Shake the can vigorously before use.



REF	Pack contents	Shipping unit	PZN
00994D	500 ml	12 bottles	11170248

# Gentle Med® moist care wipes

Our care wipes for normal and sensitive skin

### Features

Moist care wipes with a witch hazel extract texture, based on mild, coconut oil surfactants. Gentle Med® moist care wipes are PEG-, paraben-, paraffin- and silicone-free.

### Usage

Ready-to-use, high-quality care wipes for the effective, rapid and especially mild care and cleansing of soiled areas of the skin. Disposable to prevent germs spreading any further. For use without water. Please re-seal the pack carefully after removing wipes.



REF	Size	Pack contents	Shipping unit	PZN
2356	18 x 20 cm (approx.)	80 pcs	12 packs	12357658

# Application aids

Our application aids for the metering and dosing of disinfectants and care products

### Metering pumps/attachments

Enable products to be pumped out of the original containers comfortably and hygienically.

REF	Article		PZN
00929	Metering pump	500 ml	–
00930	Metering pump	1 l	–
00931	Canister metering pump	2 l, 5 l and 10 l	–

### (Foam) spray heads

Ensure even wetting, even in difficult-to-access areas, such as for extensive spraying with surface disinfectants. Also suitable for overhead use.

REF	Article	PZN
00934	Foam spray head 1 l	–

### Canister tap

Ensure products can be removed from cans cleanly and without dripping.

REF	Article	PZN
00996	Canister tap 2 l, 5 l and 10 l	–

### Canister key

For effortlessly opening cans.

REF	Article	PZN
00932	Canister key 5 l and 10 l	–

### Dosing bottle/graduated flask

For removing the product in a user-friendly way when decanting or for making ready-to-use solutions or concentrates.

REF	Article		PZN
00943	Dosing bottle	1 l	–
00945	Graduated flask	1,000 ml	–
00927	Graduated flask	2,000 ml	–



# Disinfectant dispenser

For the dosed dispensing of hand disinfectants, liquid soaps and lotions

## Aluminium disinfectant dispenser

Soap/disinfectant dispenser for disposable bottles and refill containers, aluminium housing, short operating lever, 500 ml or 1,000 ml. Suitable for wall mounting. Models with reusable or disposable pump, can be replaced from the front.

REF	Product name	Model
1417549	Aluminium dispenser 500 ml	with reusable pump
4402102	Aluminium dispenser 500 ml	with DHP disposable pump
1417589	Aluminium dispenser 1,000 ml	with reusable pump
4402107	Aluminium dispenser 1,000 ml	with DHP disposable pump



## Automatic disinfectant dispenser

Automatic soap/disinfectant dispenser for an even greater level of hygiene. With drip tray, made from rugged plastic, 500 ml, with rear wall, lockable. Suitable for wall mounting. Models with reusable or disposable pump, can be replaced from the front.

REF	Product name	Model
4402237	Automatic plastic dispenser 500 ml	with reusable pump
4402283	Automatic plastic dispenser 500 ml	with DHP disposable pump



## Plastic disinfectant dispenser

Manual dispenser for disinfectants, liquid soaps and low-viscosity lotions, suitable for wall mounting, white plastic housing with white drip tray, stainless steel operating lever. Model with DHP disposable pump or reusable pump, can be replaced from the front.

REF	Product name	Model
4401841	Plastic dispenser 500 ml	with reusable pump
4401974	Plastic dispenser 500 ml	with DHP disposable pump



## DHP disposable pump

For even greater hygiene and safety!  
Recyclable plastic disposable pump for alcohol-based hand disinfectants, liquid soaps and lotions. Suitable for use with all our dispensers.

### Features:

- packed individually for good hygiene, quick and easy use
- simple and hygienic installation: packaging is used as a glove
- no need to clean the pump or the dispenser
- cost optimisation: time-consuming dispenser reprocessing is not necessary
- no chance of accidental reuse: pump is destroyed when container is changed
- good for the environment: completely recyclable and produced using 100% renewable energy

REF	Product name
4401857	manual dispenser 500 ml
4401858	manual dispenser 1,000 ml
4401860	touchless dispenser 500 ml







## Everyday products – and therefore essential!

Work processes in the healthcare sector should primarily be patient-centred. For this reason, we have a comprehensive range of ward and patient supplies in our portfolio to properly meet this everyday need. From daily physical hygiene for bed-ridden patients with our care products and washing gloves, to essential items such as cotton buds, wooden tongue depressors or protective towels.



# BagBath® haircare

## Our dry shampoo

### Features

A water-free hair washing system with no rinsing. Can be used either heated up or cooled down. Contains no perfumes. Can also be used in combination with the BagBath body wash system.

### Usage

Hair washing system for bed-ridden patients and people requiring care.



REF	Pack contents	Shipping unit	PZN
APO021012	100 ml	30 bottles	-

# BagBath®

## Disposable body wash system without water or soap

### Features

The wipes in the BagBath® have a pleasant and subtle scent. Moisturising ingredients and clinically tested formulas for these ingredients work to protect the skin from irritation.



### Usage

The packaging can be opened easily and then resealed, even with gloves. Wash bowls, soap, water, flannels, towels and care lotion are no longer required.


REF	Size	Pack contents	Shipping unit	PZN
2359	25 x 20 cm	8 wipes	30 packs	06149252

BagBath® enables hospitals and nursing care centres to improve the quality of skin care for their patients, and also while saving time and money.

BagBath® makes it possible to wash the body without soap, washcloths and towels in situations where there is either no bath available or it is not possible to get the patient to the bath.

### All body parts are covered:

BagBath® ensures full and complete body hygiene. A pack contains 8 wipes – one for each area of the body.



### Benefits for the patient:

- Avoidance of contact infections caused by traditional washing utensils such as flannels and towels.
- Prevention of cross-contamination by dispensing with reusable washing materials.
- Washing with BagBath® has fewer risks for patients who have undergone surgery, been intubated, suffered trauma or who are immobile, since there is no need to move the patient as is the case when the bed clothing becomes wet.

### Special pack that can be warmed up:

#### ➔ Microwave + BagBath®

- Most patients prefer the feeling of a warm bath.
- Warmth prevents the patient losing body heat.



### Special pack that can be chilled:

#### ➔ Refrigerator + BagBath®

- A refreshing and cooling effect on patients



### Simple to dispose of:

#### ➔ rubbish bin

As an adisposable product, BagBath® can be disposed of with normal household waste. When washing patients who have been quarantined, internal disposal regulations must be respected.



## BeeSana® Wash glove

Hygienic, disposable wash gloves

### Features

- disposable wash glove in various designs



### Usage

For skin cleansing and personal care. Can be used in all areas of nursing, with a particular focus on the care of the elderly and sick.

REF	Product name	Model	Pack contents	Shipping unit	PZN
3709	Film	Ultrasonically welded	50 pcs	20 packs	05461042
3723	Molton/soft	Ultrasonically welded	50 pcs	20 packs	11703466
3724	Molton/soft	Sewn	25 pcs	40 packs	05529001
3725	Molton 80 g	Sewn	50 pcs	20 packs	09336931

## Protective towel

### Features

- protective towels made from a layer of cellulose and PE film
- PE inner coating lessens wetness leaking through onto the clothing
- with punched straps for fixing around the neck, plus a movable catch pouch lower down
- thanks to the soft cellulose outer layer, can also be used as a serviette after the meal
- latex-, PVC- and DEHP-free
- colour: Apricot

### Usage

Disposable.

REF	Size	Pack contents	Shipping unit	PZN
2997	38 x 68 cm	100 pcs	6 packs	11727314



## Protective bib

### Features

- protective bib made from medical gauze
- 6-ply
- white

REF	Size	Pack contents	Shipping unit	PZN
4239	20 x 30 cm	10 pcs	10 packs	12724401



## BeeSana® Disposable briefs

Briefs to give a feeling of secure privacy

### Features

- disposable briefs made from non-woven material
- with rubber drawstring
- white with coloured waistband for easy sizing

### Usage

To protect the privacy of patients during examinations and after surgery.

REF	Size	Waist colour	Waist size	Pack contents	Shipping unit	PZN
2251	M	Blue	51–55 cm	20 pcs	30 packs	11857318
2252	L	Brown	54–58 cm	20 pcs	30 packs	11857324
2253	XL	Green	58–62 cm	20 pcs	30 packs	11857330
2254	XXL	Red	62–64 cm	20 pcs	30 packs	11857301



# Anti-slip sock universal

Our sock to prevent slips and falls

## Features

- disposable ABS sock (laid flat)
- with grips as an anti-slip coating
- made from cotton
- black
- universal size
- washable at 30 °C



## Usage

Disposable, usable for many different applications in hospitals and nursing homes.

REF	Size	Pack contents	Shipping unit	PZN
9605	40 x 9 cm	20 pcs	8 packs	12671811

# Surgical sock

Our cotton sock

## Features

- disposable sock (laid flat), without heel
- top with rubber drawstring
- white
- washable at 30 °C



## Usage

Suitable for many applications. Disposable, for use in the operating theatre, lab or clean rooms.

REF	Size	Pack contents	Shipping unit	PZN
9837	39 x 8 cm	40 pcs	20 packs	12671805
9836	44 x 8 cm	40 pcs	20 packs	12671797

# PE shoe covers

Shoe covers made from CPE with rubber drawstring

## Features

- with rubber drawstring
- disposable
- colour: blue

Dimensions (approx.): Length: 41 cm  
Height: 15 cm  
Thickness: 0.03 mm



REF	Pack contents	Shipping unit	PZN
70-151	1000 pcs	1 box	11363007

# Cellulose sponge sheets

## Features

- highly bleached (German pharmacopoeia (DAB) standard)
- absorbent

## Usage

Cleaning and hygiene cloth for medical applications.

REF	Size	Pack contents	Shipping unit	PZN
12431	20 x 20 cm	300 g	26 packs	13655576



# Layered cellulose

Crepe cuts in loose layers for multiple uses

## Features

- layered cellulose, unbleached, made from 100% recycled cellulose
- many different uses

## Usage

Can be used as an underlay for non-sterile activities, as an absorbent material or for cleaning and wiping.

REF	Size	Pack contents	Shipping unit	PZN
102200D	20 x 20 cm	5 kg	3 bags	06803349
102300D	20 x 30 cm	5 kg	3 bags	06803355
104600D	40 x 60 cm	5 kg	3 bags	06803361





# Rollicel®/medical crepe

Exam table covers on a roll

## Features

- exam table covering made from cellulose
- absorbent and micro-perforated
- several available variants
- colour: white

## Usage

For treatment tables, and partially-clothed examinations and treatments.



REF	Size	Layers	Coating	Shipping unit	PZN
3739-50	39 cm x 50 m	2-ply	No	9 rolls	11350476
3739-50B	39 cm x 50 m	2-ply	Yes	6 rolls	11350571
3739-100	39 cm x 100 m	2-ply	No	6 rolls	11350536
3745-100B	45 cm x 100 m	2-ply	Yes	6 rolls	11350884
3750-50	50 cm x 50 m	2-ply	No	9 rolls	11350482
3750-50B	50 cm x 50 m	2-ply	Yes	6 rolls	11350855
3750-100	50 cm x 100 m	2-ply	No	6 rolls	11350542
3755-50	55 cm x 50 m	2-ply	No	9 rolls	11350499
3755-50B	55 cm x 50 m	2-ply	Yes	6 rolls	11350861
3755-100	55 cm x 100 m	2-ply	No	6 rolls	11350559
3759-50	59 cm x 50 m	2-ply	No	9 rolls	11350513
3759-50B	59 cm x 50 m	2-ply	Yes	6 rolls	11350878
3759-100	59 cm x 100 m	2-ply	No	6 rolls	11350565
3759-100B	59 cm x 100 m	2-ply	Yes	6 rolls	11350890

# Protective sheets

Disposable sheet as a protective layer for sitting or lying

## Features

- 2-ply cellulose with transparent PE film that prevents fluids from seeping through
- latex-, DHEP- and PVC-free.
- reinforced lengthwise with rugged polyester stitching and side barriers.

## Usage

Disposable protective sheets to cover stretchers, examination tables and gurneys, work tables and mattresses.

REF	Size (cm)	Shipping unit	PZN
1988	80 x 175	100 pcs	04353787
1992	80 x 200	100 pcs	04353793
1995	80 x 210	100 pcs	04353801



# Protective covering

Disposable cover as a protective layer for sitting or lying, the cost-effective alternative to protective sheets

## Features

- 1-ply cellulose with transparent PE film that prevents fluids from seeping through
- latex-, DHEP- and PVC-free.

## Usage

Disposable protective sheets to cover stretchers, examination tables and gurneys, work tables and mattresses.

REF	Size (mm) length x width	Side barriers	Shipping unit	PZN
1993	390 x 375	No	250 pcs	04353681
1990	900 x 800	Yes	200 pcs	04353735
2982	1400 x 800	Yes	100 pcs	01082431
2984	2100 x 800	Yes	100 pcs	01082454



# Multi-purpose cloths

## Features

- multi-purpose cloth made from thick felted non-woven
- tear-resistant, absorbent and lint-free

## Usage

For all cleaning work in the kitchen, laboratory, hospitals and hotel rooms.

REF	Cloth format	Pack contents	Shipping unit	PZN
3539	32 x 38 cm	50 pcs	20 packs	–



# Protective bed coverings

Transparent HDPE

## Features

- protective HDPE (high-density polyethylene) cover

## Usage

For hospital and nursing beds, to protect the fresh, newly-made bed.

REF	Cloth format	Pack contents	Shipping unit	PZN
1985	320 x 95 cm	200 pcs	1 roll	–



# Universal cloth

## Features

- disposable product for cleaning and hygienic care of intact skin
- high-viscose version also optionally available

## Usage

Disposable product for use in all areas of care.



### Universal cloth

REF	Size	Pack contents	Shipping unit	PZN
6051	30 x 30 cm (approx., unfolded) 15 x 15 cm (approx., folded)	50 pcs	10 packs	04355415

### Universal soft wipe (thicker quality)

REF	Size	Pack contents	Shipping unit	PZN
6052	30 x 30 cm (approx., unfolded) 15 x 15 cm (approx., folded)	50 pcs	10 packs	11722759

# Wipe cloths

Recyclable wipes made from cellulose

## Features

- multi-layered, semi-bleached wipes made from cellulose
- latex-, PVC- and DEHP-free

## Usage

Disposable for use as a cleaning material.



REF	Size	Number of layers	Pack contents	Shipping unit	PZN
6020	39 x 30 cm	3-ply	50 pcs	20 packs	11857270
6021	39 x 30 cm	4-ply	50 pcs	20 packs	11857264

# Cosmetic wipes

## Features

- made from 2-ply cellulose
- colour: brilliant white

REF	Size	Pack contents	Shipping unit	PZN
6701	20.5 x 21 cm	100 pcs	35 packs	-



# V-fold paper hand towels

## Features

- dimensions: 25 x 23 cm

### 1-ply, natural

REF	Pack contents	Shipping unit	PZN
5722	20 x 250 pcs	1 pack	-

### 2-ply, white

REF	Pack contents	Shipping unit	PZN
5721	20 x 160 pcs	1 pack	-



# BaSick Bag

## Sick bag

### Features

- for use with the onset of nausea
- partially transparent, tear-proof sick bag made from polyethylene
- with graduation and seal ring

### Usage

Disposable, for use with nausea.



REF	Size, height	Capacity	Pack contents	Shipping unit	PZN
9401	16 x 40 cm	1,500 ml	25 pcs	20 bags	13818825
SP9401	47 cm	3 packs	1 pc.	1 pc.	–

# Kidney dish

Non-sterile, paperboard

REF	Size	Colour	Shipping unit	PZN
9400	255 x 150 x 48 mm	Grey	300 pcs	01082164



Double sterile, plastic, pack of 1 pc.

REF	Size	Pack contents	Shipping unit	PZN
4249	700 ml	12 x 1 pc.	2 packs	–



# Disposable comb

Plastic comb, individually wrapped in cellophane

REF	Colour	Pack contents	Shipping unit	PZN
9679	Brown	100 x 1 pc.	20 packs	11857672



# Disposable razor

### Features

- stainless steel, non-sterile disposable razor

REF	Razor blades	Pack contents	Shipping unit	PZN
9610	Single-edged	100 pcs	10 packs	02501189



# Urine bottle

Urine bottle for men 1,000 ml

### Features

- opaque, graduated urine bottle with clip fastener

REF	Shipping unit	PZN
3006	25 pcs	–



# Urine bottle holder

Urine bottle holder with chain and lid

### Features

- urine bottle holder for securely picking up the urine bottle

### Usage

The urine bottle holder should be attached (hung on) within the reach of the patient to the bed frame of the patient's bed. Clean as necessary. **IMPORTANT:** When using disinfectants, always respect the manufacturer's instructions for use.

REF	Unit	Shipping unit	PZN
3001	1 pc.	20 pcs	–





# Nurse's stethoscope

Nurse's stethoscope for auscultation

### Features

- acoustic flat-head stethoscope
- mechanical
- ear hooks with earpieces, tube and diaphragm
- tube length: approx. 60 cm
- latex-free

### Usage

The stethoscope must only be used by authorised medical personnel. The stethoscope is placed on the part of the body relevant for the exam and used by exerting only a little pressure.



REF	Pack contents	Shipping unit	PZN
9805	1 pc.	100 pcs	11857459

# Tourniquet

High-quality tourniquet with innovative method of use

### Features

- usable with one hand, no "reaching around"
- with safety fastener
- contains no latex or metal



REF	Article	Colour	PZN
5020	Meditrade tourniquet	Colourful designer band	06177892
N20200	Jetpull	Blue designer band	–

# Dosicair®

Our infusion regulator with injection adapter

### Features

- luer taper connection system
- can be connected to infusion systems

### Usage

For infusion treatments.



REF	Pack contents	Shipping unit	PZN
DF100	50 pcs	8 packs	–

# Urine beaker

REF	Size	Pack contents	Shipping unit	PZN
3009	125 ml	50 pcs	10 packs	04354427
3010	100 ml	500 pcs	1 pack	–



# Ultrasound scan gel

### Features

- ultrasound gel for scans using ultrasound probes

### Usage

For all types of ultrasound scans.

REF	Pack contents	Shipping unit	PZN
6831	250 g	30 bottles	01904778



# Hot-Cold Pack

## Pad for cold and warm applications

### Features

- permanently gel-filled pad
- can be easily shaped to fit
- suitable for the freezer (stays flexible to  $-30\text{ }^{\circ}\text{C}$ ), water bath and microwave
- follow information on label

### Usage

For strains, pulled muscles and swelling.  
Also for fever, headaches and bruising.



	REF	Size	Pack contents	Shipping unit	PZN
Gel colour: transparent	9668	13 x 14 cm	50 pcs	1 pack	04352500
	9870	12 x 29 cm	25 pcs	1 pack	04352463
	9666	16 x 26 cm	25 pcs	1 pack	04352517
	9661	30 x 40 cm	6 pcs	1 pack	04352546
	9663	7.5 x 35 cm	40 pcs	1 pack	04352552
Gel colour: blue	9465	16 x 25 cm	20 pcs	1 pack	–

# Hot-Cold Pack Sleeve

## Non-woven sleeve for hot-cold packs

### Features

- non-woven covers for hot-cold packs
- used to protect the skin, since direct contact between the pack and skin must be avoided
- for disposable use only

REF	Size	Pack contents	Shipping unit	PZN
9682	13 x 14 cm	500 pcs	1 pack	–
9683	12 x 29 cm	500 pcs	1 pack	–
9684	16 x 26 cm	500 pcs	1 pack	–
9687	21 x 38 cm	500 pcs	1 pack	–
9693	30 x 40 cm	500 pcs	1 pack	–

# Rescue blanket

### Features

- colour: gold/silver
- reflects bodily heat or solar radiation
- reduces the risk of hypothermia or works to shield the body from the heat of the sun's rays
- packed individually with instructions for use.

### Usage

Cold and heat protection for injured parties.  
Usable for accidents in windy, rainy or hot conditions.



REF	Size	Pack contents	Shipping unit	PZN
4070	160 x 210 cm	50 pcs	5 packs	13510895

# Medicine Pod

## Medicine Pod, 30 ml, graduated



### Features

- made from transparent polypropylene
- contents visible
- with graduation at side



REF	Colour	Pack contents	Shipping unit	PZN
<b>Medicine Pod</b>				
9729	Transparent	75 pcs	50 packs	–
9727	Blue	75 pcs	50 packs	–
9725	Red	75 pcs	50 packs	–
9726	Yellow	75 pcs	50 packs	–
9728	Green	75 pcs	50 packs	–
<b>Medicine Pod lid</b>				
9700	White	750 pcs	5 packs	–

# Medication dispenser

## Features

- tablet dispenser made from PP/PVC with transparent sliding lid
- helps patients to take their medication regularly and on time in the morning, at lunchtime, in the evening and at night



REF	Colour	Pack contents	PZN
9707	White	200 pcs	–
9708	Blue	200 pcs	–

# Sharps Bin

Container for securely disposing of needles and syringes

## Features

- round disposal container
- with lid and round opening
- flat closure
- pull-off and screw-off aids



## Usage

For disposing of needles and syringes.

REF	Size	Pack contents	PZN
9694	0.5 l	1 pc.	–
9702	2.5 l	1 pc.	–

# Glove dispenser

Wall mounted glove dispenser for boxes of all sizes

## Features

- acrylic or polystyrene
- easy to fit



REF	Article	Product name	Dimensions (in mm)	Shipping unit	PZN
4040	Dracula	PS wall holder, without overprint (white)	Universal	1 pc.	–
4050	Ghost	Acrylic wall holder with overprint	250 x 145 x 100	1 pc.	–
4059	Ghost	Acrylic wall holder without overprint	250 x 140 x 70	1 pc.	–
4060	Ghost	Acrylic wall holder without overprint	250 x 145 x 100	1 pc.	–



## Rubbish bags HD

### PE rubbish bags

#### Features

- rubbish bags made from HDPE (high-density polyethylene)
- for collecting general waste
- disposable



REF	Size	Colour	Pack contents	Shipping unit	PZN
9747	49 x 60 cm (approx. 24 l)	Transparent	50 pcs	40 rolls	–
9748	63 x 74 cm (approx. 60 l)	Transparent	50 pcs	40 rolls	–

## Rubbish bags LD

### PE rubbish bags

#### Features

- rubbish bags made from LDPE (low-density polyethylene)
- for collecting general waste
- disposable



REF	Size	Colour	Pack contents	Shipping unit	PZN
9651	70 x 110 cm (approx. 120 l)	Blue	25 pcs	10 rolls	–
9746	50 x 60 cm (approx. 30 l)	Grey	50 pcs	20 rolls	–
9749	63 x 72 cm (approx. 60 l)	Transparent	50 pcs	20 rolls	–
9750	63 x 85 cm (approx. 65 l)	White	40 pcs	9 rolls	–
9760	57.5 x 100 cm (approx. 70 l)	Blue	25 pcs	10 rolls	–
9761	70 x 110 cm (approx. 120 l)	Blue	25 pcs	10 rolls	–
9762	70 x 110 cm (approx. 120 l)	Transparent	25 pcs	10 rolls	–
9764	60 x 80 cm (approx. 60 l)	Grey	25 pcs	10 rolls	–
9765	70 x 110 cm (approx. 120 l)	Red	25 pcs	10 rolls	–
9766	57.5 x 100 cm (approx. 70 l)	Yellow	25 pcs	10 rolls	–
9768	70 x 110 cm (approx. 120 l)	Blue	15 pcs	10 rolls	–

## Kitchen towel

#### Features

- brilliant white
- made from cellulose
- 50 sheets per roll
- roll length: 11.2 m
- sheet width: 25.7 cm
- sheet length: 22.4 cm



### 3-ply

REF	Pack contents	Shipping unit	PZN
57301	4 rolls	8 packs	–

## Toilet paper

#### Features

- brilliant white
- made from cellulose
- 250 sheets per roll
- roll length: 28 m
- sheet width: 9.5 cm
- sheet length: 11.2 cm



### 2-ply

REF	Pack contents	Shipping unit	PZN
9200	8 rolls	8 packs	–

### 3-ply

REF	Pack contents	Shipping unit	PZN
9201	8 rolls	9 packs	–



# Index

## A

ABE® Absorbent adhesive dressing 79  
 ABE® body bandage 97  
 ABE® cannula retention dressing 81  
 ABE® derm (plus) 80  
 ABE® drainal 69  
 ABE® extra 92  
 ABE® extra cohesive 93  
 ABE® fix color 87  
 ABE® gauze bandage 84  
 ABE® Gilchrist's bandage 96  
 ABE® heart catheter pressure bandage 98  
 ABE® Ideal 100 89  
 ABE® injection plaster 78  
 ABE® Last 85  
 ABE® Last cohesive 88  
 ABE® Last crepe 86  
 ABE® lastic 91  
 ABE® light 78  
 ABE® limb holder 98  
 ABE® non-woven material 81  
 ABE® nose bandage 72  
 ABE® padded bandage 92  
 ABE® pressure bandage roll 99  
 ABE® ribbon gauze 74  
 ABE® silk 82  
 ABE® tracheal 69  
 ABE® tubular bandage 94  
 ABE® tubular net bandage 95  
 ABE® universal 90  
 ABE® Ven 93  
 ABE® X-Ray ribbon gauze 105  
 Absorbent cushions 121  
 Absorbent dressing pads 70  
 Absorbent pads 122  
 Active gel 145  
 Alcohol swabs 75  
 Anti-slip sock universal 156  
 Application aids 148-149  
 Arm bandages 96  
 Astronaut-shaped surgical cap 53

## B

BagBath® (Haircare) 152  
 Bandage change sets 112  
 Bandages 84-97  
 BaSick Bag 162  
 BeeSana® wash gloves 154  
 BeeSana® absorbent cushion 121  
 BeeSana® absorbent dressing pad 70  
 BeeSana® absorbent pad 122  
 BeeSana® bandage change set 114  
 BeeSana® cellulose swab 75  
 BeeSana® Comflock Super 125  
 BeeSana® disposable brief 155  
 BeeSana® eye pad 71  
 BeeSana® female incontinence pad 123  
 BeeSana® gauze balls 73  
 BeeSana® medical gauze 62  
 BeeSana® gauze cutting 62  
 BeeSana® gauze swab 63 - 65  
 BeeSana® male incontinence pad 123  
 BeeSana® medical gauze 62  
 BeeSana® Medi-Pants 117  
 BeeSana® non-woven swab 67  
 BeeSana® PEG bandage set 113  
 BeeSana® pointed swab 74  
 BeeSana® PP gown t 54  
 BeeSana® PP scrubs 54  
 BeeSana® PP/PE gown 55  
 BeeSana® pre-cut gauze 62  
 BeeSana® rectangular pad 122  
 BeeSana® SafeCut 83  
 BeeSana® Slit gauze swab 66  
 BeeSana® Slit non-woven swab 68  
 BeeSana® SMS gown 55 / 56  
 BeeSana® suture removal set 114  
 BeeSana® ward catheter set 112  
 BeeSana® ward protection set 115  
 BeeSana® Wash glove 154  
 BeeSana® X-Ray gauze balls 104 / 107  
 BeeSana® X-Ray gauze swab 102-103  
 BeeSana® X-Ray lap sponge 100-101  
 BeeSana® X-Ray pointed swabs 106  
 BeeSana® X-Ray preparation swabs 105

Body bandages 97  
 Body wash systems 152  
 Box holder 169

## C

Canister keys 147  
 Canister metering pump 147  
 Canister/outlet taps 147  
 Cannula fixations/plasters 81  
 Cellulose sponge sheets 157  
 Cellulose swabs 75  
 Chemical protective suit 56  
 Combs 162  
 Comflock Super 125  
 Compression bandages 89-93  
 Copolymed® 28  
 Copolymer gloves 28  
 Cosmetic wipes 161  
 Cotton buds 76  
 Cotton gloves 24  
 Cream lotion 144

## D

Dental masks 46-52  
 Disinfectant dispenser 148-149  
 Disinfectant wipes 133-134, 137-138  
 Disinfection (products) 130-141  
 Disposable aprons 57  
 Disposable body wash system 152  
 Disposable briefs 155  
 Disposable comb 162  
 Disposable gloves 9-45  
 Disposable razors 163  
 Disposal container 168  
 Dosicair® 165  
 Dosing aids 147  
 Dosing bottle 147  
 Dressing material 62-107  
 Dressing pads 70

## E

Ethasept® 130  
 Ethasept® Gel 131  
 Exam table covers 158  
 Examination gloves, (non-)sterile 9-23, 28-30  
 Eye pads 71

## F

Face masks 47-52  
 Female incontinence pad 123  
 Female pads 123  
 Film wound bandages 80  
 Finessis® Aegis® 43  
 Finessis® Corium® 45  
 Finessis® Zero® 44  
 Finger cots 26  
 Fixation bandages 83-88  
 Fixation briefs 118 / 119  
 Fixation dressings 95  
 Fixation plasters 81-82  
 Flexylon gloves 42-45  
 Flock (patient) undersheets 124 / 125  
 Flock incontinence briefs 118  
 Foam spray heads 147  
 Froggy™ 25

## G

Gauze balls 73, 104, 107  
 Gauze bandage 84  
 Gauze swabs 62-72, 102-103  
 Gentle Med® active gel 145  
 Gentle Med® care foam 146  
 Gentle Med® cream lotion (W/O) 144  
 Gentle Med® moist care wipes 146  
 Gentle Med® oil bath 145  
 Gentle Med® skin care lotion (O/W) 143  
 Gentle Med® skin cleansing and care foam 146  
 Gentle Med® wash lotion 143  
 Gentle Skin® Aloecare™ 11  
 Gentle Skin® black 12

Gentle Skin® compact + 10  
 Gentle Skin® grip 10  
 Gentle Skin® Isopretex® OP 40  
 Gentle Skin® Micro OP® 37  
 Gentle Skin® Premium OP 34  
 Gentle Skin® Securitex OP 39  
 Gentle Skin® sensitive 9  
 Gentle Skin® Smooth OP 36  
 Gentle Skin® sterile 29  
 Gentle Skin® Superior OP 35  
 Gilchrist's bandage 96  
 Glove box holder 169  
 Glove dispenser 169  
 Glove liners/undergloves 24, 28, 37, 45  
 Gloves 8-45  
 Gowns 54-56  
 Graduated flasks 147

## H

Hair wash system 152  
 Hand disinfection 130 / 131  
 Hand spray head 147  
 Heart catheter pressure bandage 98  
 Hot-Cold Pack 166  
 Hot-Cold Pack Sleeve 166  
 Household gloves 27  
 Hygiene gloves 24

## I

Ideal bandages 89  
 Incontinence briefs 116 / 117  
 Incontinence briefs (with leg) 118 / 119  
 Incontinence lines 120  
 Incontinence pads 120-123  
 Incontinence products 115-125  
 Incontinence undersheets 124-125  
 Infusion regulators 187  
 Injection plasters 78  
 Instrument disinfection 140 / 141  
 Instrusol® AF + 140 / 141

## K

Kidney dishes 162  
 Kitchen towel 171

## L

Lap sponge with X-ray contrast 100-101  
 Latex gloves 9-12, 26-27, 29, 33-39  
 Layered cellulose 157  
 Layered undersheets 118  
 Leg cuffs 98  
 Liberty Light 120  
 Light fixation briefs 118  
 Limb holder 98  
 Long-stretch bandages 93

## M

Male incontinence pads 123  
 Masks 48-52  
 Maxitex® UG 38  
 Medication dispenser 172  
 Medi-Pants 117  
 Medicine Pod 167  
 Medical crepe 158  
 Medical gauze 62  
 Medication dispenser 168  
 Medium-stretch bandages 90  
 Medizid® AF+ 135  
 Medizid® Alpha+ 136  
 Medizid® Alpha+ Flowpack 138  
 Medizid® Alpha+ XXL (disinfectant wipes) 137  
 Medizid® Rapid QF Flowpack 134  
 Medizid® Rapid+/XXL (disinfectant wipes) 135  
 MEDlzid® sanitary cleaner 155  
 Metering pumps 147  
 Moisturising wipes 146  
 Mouth tongue depressor 77  
 Multi-purpose cloths 159  
 MyClean® easy touch 27



# Index

## N

Nasal slings 72  
 Needle collection boxes 168  
 Neoprene gloves 40  
 Neopretex® OP 41  
 Nitril® 3000 13  
 Nitril® 3000 X-Long Blue 15  
 Nitril® black 16  
 Nitril® LIOX 20 / 21  
 Nitril® Magenta 18  
 Nitril® NextGen® 13  
 Nitril® Sensory® 14  
 Nitril® sterile 30  
 Nitril® Viola 17  
 Nitrile gloves 13-21, 30  
 Non-slip socks 160  
 Non-woven (slit) swab 67-68  
 Non-woven material 81  
 Nose bandages 72  
 Nurse's caps 53  
 Nurse's stethoscopes 164  
 Nursing pads 126

## O

Oil bath 145

## P

Padded bandages 92  
 Pads for post-natal women 122-123  
 Panty fixation briefs 119  
 Paper hand towels 157  
 PE aprons 57  
 PE disposable aprons 57  
 PE gloves 25, 28  
 PE shoe covers 157  
 PEG bandage set 111  
 Plasters 78-82  
 Pointed swabs 74, 106  
 Polychloroprene gloves 41  
 Polyisoprene gloves 40  
 PP/(PE) gowns/scrubs/aprons 52-54  
 Pressure bandages 99  
 Protective bed coverings 159  
 Protective bibs 155  
 Protective coverings 159  
 Protective Respirator Mask 52  
 Protective scrubs/gowns/aprons 54-59  
 Protective sheets 158  
 Protective towel 154

## R

Razors 163  
 Ready-to-use finger bandages 83  
 Rectangular pads 122  
 Reference™ 9  
 Reference™ OP 33  
 Rescue blankets 167  
 Respima FFP2 52  
 Ribbon gauze with X-ray 105  
 Roll plasters 82  
 Rollicel®/medical crepe 158  
 Rubbish bags 170

## S

SafeCut 83  
 Scrotum roll 99  
 Set solutions 110-113  
 Sharps Bin 168  
 Shoe covers 157  
 Short-stretch bandages 89, 91-93  
 Shoulder bandage 98  
 Sick bags (BaSick Bag) 162  
 Skin care lotion 144  
 Skin cleansing foam 146  
 Skin disinfection 134  
 Slit gauze swab 66  
 Slit (hole) dressing pads 69  
 SMMS surgical scrubs 58 / 59  
 SMS gowns 55 / 56  
 Soap dispensers 148-149  
 Socks 156  
 Stethoscopes 164  
 Suavel® Antifluid 51  
 Suavel® Antifog 51  
 Suavel® Apollo 53  
 Suavel® Astrid 53  
 Suavel® Comfort 50  
 Suavel® Comfort Plus 50  
 Suavel® OP PRO Expert 59  
 Suavel® OP PRO Standard 58  
 Suavel® Protec 49  
 Suavel® Protec Black 49  
 Suavel® Protec Plus 48  
 Suavel® Sensima 52  
 Suavel® Yeti Coverall 56  
 Superform® 120  
 Superform® Air 121  
 Superslip® 116  
 Superslip® Air Comfort 117  
 Support bandages 95-99

## V

V-fold paper hand towels 161  
 Vinyl 2000 PF 22  
 Vinyl Blue Gen 23  
 Vinyl gloves 22 / 23  
 Visitor gowns 54-56

## W

Wall holders 169  
 Ward catheter sets 110  
 Ward gowns/aprons 54-56  
 Ward protection set 113  
 Wash lotions 143  
 Wet (moisturising) wipes 146  
 Wipe cloths 160  
 Wipes 133-134, 137-138  
 Wooden tongue depressor 77  
 Wound bandages 80-82  
 Wound plasters 78-80

## X

X-ray wound care 100-107

Surface disinfection 135-141  
 Surgical caps 53  
 Surgical caps 53  
 Surgical dental masks 47-52  
 Surgical face masks 47-52  
 Surgical gloves 33-45  
 Surgical gowns/scrubs/aprons 54-59  
 Surgical masks 47-52  
 Surgical sets 108-113  
 Surgical socks 156  
 Suture removal sets 112  
 Swabs 73-75, 104-107

## T

Tablet dispensers 168  
 Toilet paper 171  
 Tourniquet 164  
 Tube bandages 96 / 97  
 Tubular bandage 94  
 Tubular net bandage 95

## U

Ultrasound scan gel 165  
 Undersheets 124 / 125  
 Universal bandages 90 / 91  
 Universal cloth 160  
 Universal wipes 160  
 Urine beaker 165  
 Urine bottle (holder) 163



**Senator Hans-Hermann Rösner-Mautby**

(\* 27 June 1917 in Beuthen, † 24 June 1999 in Oberaudorf)

For his services towards the establishment of the hospital in Oberaudorf, he was made honorary senator of Ludwig-Maximilian University, Munich. He was also awarded the German Federal Cross of Merit and the Bavarian Order of Merit for his entrepreneurial achievements in the healthcare sector.

## SENATOR RÖSNER FOUNDATION



### Touch. Help. Create opportunities.

With its wide-ranging activities and carefully targeted public relations work, the Senator Rösner Foundation collects funds for people who have fallen into financial difficulty because of cancer. Founders Peter Rösner and his wife Eveline Rösner use their foundation to support cancer research and the training of medical personnel.

As a successor to the "Hilfswerks der Klinik Bad Trissl" (Friends of Bad Trissl Hospital), the Senator Rösner Foundation can now look back on over 25 years of experience as a charitable organisation.

The following projects are currently receiving long-term support.

#### HEAT study

The Senator Rösner Foundation is supporting the clinical research study led by Prof. Dr. med. Dipl.-Biochem. Rolf D. Issels at Munich University Hospital with its innovative approaches to the field of pancreatic cancer; funding volume here is several hundred thousand euros.

#### Rapid help where it's needed most

Cancer can leave some people facing unexpected financial and social difficulties. Here, the Senator Rösner Foundation "Leben mit Krebs" (Living with Cancer) provides quick and straightforward help.

#### Promoting skills and professional development

The Senator Rösner Foundation takes advantage of its proximity to all the major oncological centres in Munich and Bavaria. The Foundation maintains regular contact with scientists and researchers while also providing funds for research and development.

#### We run to help

The charity run organised by the Senator Rösner Foundation is held under this banner, plus the patronage of FISWorld Champion Felix Neureuther, and is now well-known throughout Germany. In 2013 alone, almost 800 runners together covered the amazing distance of 6,756 km, and managed to amass a record-breaking sponsorship figure of € 105,000.

#### WSG-Swarovski-Wattens

#### SOS Kinderdörfer

A pioneering aid project that goes beyond national boundaries. Thanks to the bonds of friendship between the Rösner family and the charity president, Diana Langes-Swarovski, children and young people are helped via SOS Kinderdörfer to get through difficult life crises, including by other means than sporting development.

[www.wsg-wattens-fussball.at](http://www.wsg-wattens-fussball.at)

Contact: [senator-roesner-stiftung.de](http://senator-roesner-stiftung.de)

# Company development

1982	Formation of Hermes Trans Atlantik Trading as a subsidiary of our Canadian companies	1994–1995	Breakeven point reached
1988	Formation of Meditrade GmbH, with its present-day business focus and slogan: "360° protection" – hand protection, face masks, caps, operating scrubs, disinfection, incontinence products, bandages, dressing pads, swabs, abdominal dressings, etc.	1996	Shareholding acquired in Eurocare GmbH and Mediclean Bayern GmbH
1989	50% participation in Meditrade Austria (divested 1992), and shareholding in Meditrade Switzerland (divested 1992)	1996	Formation of Meditrade International Thailand Ltd. with subsequent purchase of shareholding in joint venture Siamount Plastic Ltd. Bangkok with Mitsubishi Group (vinyl factory)
1990	Turnover transfer from Cosmotec and Sauter Trading	Up to 2002	Range continuously extended with exclusivity agreements and minority shareholdings in Italy
1990–1992	Acquisition of first glove factory in Thailand, Hycare Ltd. Hatyai	2003	Construction of another new factory in Thailand Hatyai, Rösner-Mautby Thailand Ltd. (DVD)
1992	Formation of Meditrade Dresden (now Praximed Meditrade DD)	Up to 2006	Sale of all factories in Thailand
1994	Acquisition of AiA Leinberger GmbH Großbränschen	2008	Share investments in Advent Health Konzern AG, Malaysia (global leader in surgical glove production)
1992	Relocation to new company site at Tierseeestr. 196–198, Kiefersfelden	2009–2011	Acquisition of Beese GmbH & Co. KG in Hamburg with production subsidiary in China
1993	Acquisition of Hese GmbH in Feucht (Nuremberg), cellulose and undersheet maker, plus VP-Schickedanz/Temca factory	2012	Move to new company premises at Medipark 1, Kiefersfelden
1992–1994	Formation of new companies Hycare Ltd. UK/Barrier Healthcare Ltd. UK in Cambridge	2016	Majority shareholding in MaiMed Group
1995	Asset deal with MCM Landhut	2018	Expansion of logistics centre in Neuenkirchen, Meditrade GmbH celebrates 30 years in business
		2020	Formation of Meditrade UK, a Meditrade subsidiary based in London: Meditrade UK Ltd., 176 – 178 Brompton Road, SW3 1HQ, London, England <a href="http://www.meditradeuk.com">www.meditradeuk.com</a>







### **Meditrade UK Ltd.**

176 – 178 Brompton Road  
SW3 1HQ, London  
England  
Phone: +44 (0) 20 4513 2390  
Email: [info@meditradeuk.com](mailto:info@meditradeuk.com)

### **Meditrade GmbH**

Medipark 1  
83088 Kiefersfelden  
Germany  
Phone: +49 (0) 80 33/97 60-0  
Fax: +49 (0) 80 33/97 60-70

**Please contact our Customer Service  
if you have any questions about our  
products or require further information:**

Phone: +44 (0) 20 4513 2390  
Email: [info@meditradeuk.com](mailto:info@meditradeuk.com)

**You can also visit our website at  
[www.meditradeuk.com](http://www.meditradeuk.com)**

The website provides full details of our product  
portfolio as well as other information.

